

**Yes, you can:**
**Matrx MX2 Back  
June 2025**

Tel: 01656 776222  
Email: [ordersuk@invacare.com](mailto:ordersuk@invacare.com)  
[www.invacare.co.uk](http://www.invacare.co.uk)

<b>Account Name :</b>	<b>Account No :</b>
<b>Address :</b>	<b>Order No :</b>
<b>Tel No:</b>	<b>Order Qty :</b>
<b>Delivery address :</b>	<b>Order Date :</b>

## Matrx MX2

### Retail order form

Max.user weight: 113kg / 17 stone





		WIDTH				
HEIGHT	mm	355	380	405	430	455
	inch	14	15	16	17	18
	150	6				
	230	9				
	310	12				
	410	16				



#### Standard Features:

- Lightweight carbon-fibre back shell
- Anatomical contouring with 76 mm depth for enhanced lateral stability
- Angle and depth adjustable
- Cutouts for breathability to help manage microclimate
- Slotted design to allow height adjustments even when in the chair
- Breathable Spacer cover
- Backrests 230mm or higher come with leatherette side panels
- Oversized foam for enhanced comfort and protection

Code	X	Description	Price	
MX2 Back shell				
1629002		MX2 Back (35.5w x 15h)	14x6"	£519.00
1629003		MX2 Back (35.5w x 23h)	14x9"	£519.00
1629004		MX2 Back (35.5w x 31h)	14x12"	£519.00
1629005		MX2 Back (35.5w x 41h)	14x16"	£519.00
1629006		MX2 Back (38w x 15h)	15x6"	£519.00
1629007		MX2 Back (38w x 23h)	15x9"	£519.00
1629008		MX2 Back (38w x 31h)	15x12"	£519.00
1629009		MX2 Back (38w x 41h)	15x16"	£519.00
1629010		MX2 Back (40.5w x 15h)	16x6"	£519.00
1629011		MX2 Back (40.5w x 23h)	16x9"	£519.00
1629012		MX2 Back (40.5w x 31h)	16x12"	£519.00
1629013		MX2 Back (40.5w x 41h)	16x16"	£519.00
1629014		MX2 Back (43w x 15h)	17x6"	£519.00
1629015		MX2 Back (43w x 23h)	17x9"	£519.00
1629016		MX2 Back (43w x 31h)	17x12"	£519.00
1629017		MX2 Back (43w x 41h)	17x16"	£519.00

1629018		MX2 Back (45.5w x 15h)	18x6"	£519.00
1629019		MX2 Back (45.5w x 23h)	18x9"	£519.00
1629020		MX2 Back (45.5w x 31h)	18x12"	£519.00
1629021		MX2 Back (45.5w x 41h)	18x16"	£519.00
<b>Hardware</b>				
		Hardware is width-adjustable allowing the backrest to be fitted to wheelchairs that are up to 25mm narrower or wider than the actual backrest width for a precise fit and positioning in rigid or foldable wheelchairs with backrest tube diameter of 19 mm, 22 mm and 25 mm.	 <div>Fixed hardware</div>	 <div>Mini set hardware</div>
1629027		Fixed carbon fibre hardware		£153.00
1629028		Mini set (quick release) hardware		£206.00
<b>Privacy flap</b>				
1551799		Extra Small	Fits 35.5cm width back	£49.00
1551800		Small	Fits 40.5cm width back	£49.00
1551903		Medium	Fits 43 and 45.5cm width back	£49.00