



Invacare  Modulte™

Invacare® | **Fox™**

compact, versatile and  
highly adjustable

**LINX®**

**INVACARE**  
Yes, you can.™

The *Invacare Fox* is packed with innovations to make life more manageable for everyday powerchair use. **Fox** has excellent indoor and outdoor capabilities and now, with the integrated **Modulite**<sup>™</sup> seating system, this remarkable powerchair features even more adjustment possibilities to accommodate a wide range of different needs.

- Compact yet efficient
- Quick to dismantle and fold
- Easy set-up and fast delivery
- **Modulite**<sup>™</sup> Seating and **LiNX**<sup>®</sup> Controls
- Fantastic range



### Compact and efficient

The *Invacare Fox* was designed to maximise independence for individuals through its compact build. The **Fox**'s efficient motors make it ideal for those who want to travel those extra few miles in a comfortable, safe powerchair. In addition the **Fox** can be easily dismantled and quickly folded for greater transportability.

### Quick to set up and to store

The **Fox** can be changed from a short to an extended wheelbase, counteracting any stability limitations a client may have. Settings can be stored and retrieved, making it easy to set-up and adapt to each person. The entire unit can be manually tipped when folded and manoeuvred, making it easy to store and convenient to deliver\*.

\* not with Modulite seating



### Advanced Modulite<sup>™</sup> seating

Modulite seating has now been added to the **Fox** to create even more adjustment possibilities, bringing more support and ensuring a greater level of comfort. Benefits like a powered recline, lateral supports and a range of comfortable padded backs have all been added.

## Simply Smart LiNX control system

*Invacare LiNX* is our new insight inspired control system that learns over time to bring a revolutionary driving experience to individuals. Thanks to its unique modular design, this futureproof system can be seamlessly expanded and tailored to complement changing needs. This Simply Smart technology has programming and maintenance updates that can be made wirelessly through bluetooth and in real time, making it easy to service. With easy access to statistical and diagnostic information, **LiNX** allows complete visibility for professionals. Features include:



### Enhanced control system

Ergonomically designed **LiNX** user remotes bring a new level of simplicity and control. Also available is the **LiNX** attendant control allowing professionals and carers to benefit from the steady and predictable drive experience.



### Real-time programming and diagnostics

Take advantage of quicker and easier programming with **LiNX** wireless functionality. Professionals can connect to the controller on the **Fox** via Bluetooth allowing programming and maintenance updates to be made wirelessly, as well as accessing live chair diagnostics, in real time.

## Features and options



### Folding mechanism

More links than a traditional crossbar for better stress distribution and reduced wear.



### 4-point tie downs

The base on the **Fox** includes built-in, crash tested tie-downs.

### Modulite telescopic inner frame

Allows for easy seat depth, width and center of gravity adjustments.



### Aluminium rear wheels

A choice of stylish aluminium rear wheels are available with matching silver rim castors.



### Powered recline\*

The 90° - 120° recline provides a high pivot point to minimise shearing forces.



### Short and extended wheelbase settings

The modular design caters for flexible handling and installation. The short wheelbase can be set to an extended version without any additional parts.



### Battery pack

A unique, easy to remove and stylish battery pack accepts different battery sizes.



### Kerb climber

An attachment for the kerb climber is built-in, so it can be retrofitted with very little effort.



### Castor locks

The castor lock enables easy folding and is useful for those who have difficulties reversing down ramps.



### Headrest\*

Rea® neck and headrests are available as optional extras.



### Lateral support\*

Laterals and a padded back offer increased comfort and support.

\*These options are only available with Modulite.

## Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please visit [www.invacare.co.uk](http://www.invacare.co.uk)

	Seat width	Seat depth	Seat height	Backrest height	Armrest height	Length of footrests	Seat angle
<b>Modulite™</b>	330-530 + 50 mm	410 - 510 mm	480 / 510 <sup>3</sup> mm	480 - 540 mm 560 - 660 mm	245 - 310 mm	350 - 450 mm	0, 2, 4, 6, 8°
<b>Fox Seat Unit</b>	350 <sup>1</sup> , 400, 450, 500 <sup>2</sup> mm	400, 430, 460, 490 mm	480 / 510 <sup>3</sup> mm	490 / 570 mm	230 - 280 mm 205 - 255 <sup>1</sup> mm	350 - 450 mm	0, 2, 4, 6, 8°

  

	Backrest angle	Driving unit width <sup>4</sup>	Total height	Total length inc. footrests <sup>4</sup>	Total width, folded	Total weight	Max. user weight
<b>Modulite™</b>	90 - 120°	590 <sup>4</sup> mm	960 - 1060 mm	1015 / 1045 <sup>5</sup> mm	540 x 590 x 540 mm <sup>4</sup>	ca. 83,2 kg	127 kg
<b>Fox Seat Unit</b>	85, 95, 105, 115°	590 <sup>4</sup> mm	1010-1120 <sup>5</sup> mm	1015 / 1045 <sup>5</sup> mm	540 x 590 x 540 mm <sup>4</sup>	ca. 79 kg <sup>5</sup>	127 kg

  

	Battery capacity	Turning diameter	Max. climbable obstacle height	Max. safe slope	Motor capacity	Speed	Driving range
<b>Modulite™</b>	ca. 34, 40 Ah	1600 / 1630 <sup>5</sup> mm	60, 80 mm <sup>6</sup>	10,5% / 6°	2 x 200 W	3, 6 km/h	up to 40 km <sup>7</sup>
<b>Fox Seat Unit</b>	ca. 34, 40 Ah	1600 / 1630 <sup>5</sup> mm	60, 100 mm <sup>6</sup>	10,5% / 6°	2 x 200 W	3, 6 km/h	up to 40 km <sup>7</sup>



1. With inrigged armrest
2. By armrest adjustment
3. Without seat cushion at 4° tilt (450 / 480 mm at 0° tilt)
4. Wheelbase
5. Depending on configuration
6. With kerb climber
7. According to ISO 7176-4:40 km - 45 Ah (3 / 6 km/h)  
32 km - 40 Ah (3 / 6 km/h)  
27 km - 34 Ah (3 / 6 km/h)

The LINX word, mark and logos are registered trademarks owned by Dynamic Controls and any use of such marks by Invacare is under license.

Fulfils the requirements of ISO 7176-19

## Shroud colours

Silky Matt Black    Pearl Night Blue    Metallic Silver    Transparent Cherry Red    White



FOX - UK - 06/2017 - 1576887



**Like us on Facebook:**  
Invacare UK



**Follow us on Twitter:**  
@Think\_Mobility



**Follow us on YouTube:**  
Invacare Ltd - UK

Invacare Limited  
Pencoed Technology Park - Pencoed  
Bridgend CF35 5AQ - United Kingdom  
Tel: +44 1 656 776 222  
Fax: +44 1656 776 220  
E-mail: [uk@invacare.com](mailto:uk@invacare.com)  
Sales Order E-mail: [ordersuk@invacare.com](mailto:ordersuk@invacare.com)  
[www.invacare.co.uk](http://www.invacare.co.uk)

Invacare Ireland Ltd  
Unit 5 - Seatown Business Campus  
Seatown Road - SWORDS - County Dublin - Ireland  
Tel. +353 1 8107084  
Fax +353 1 8107085  
Email: [Ireland@invacare.com](mailto:Ireland@invacare.com)  
[www.invacare.ie](http://www.invacare.ie)

