


## **Iconic design that takes you anywhere**

TDX® SP2's iconic design looks every bit as good as its performance, effortlessly manoeuvring you through daily activities. Featuring a true centre wheel drive, it easily turns on the spot without deviation to manoeuvre you into and out of those confined spaces.














# Technology with a human touch

Our range of LiNX remotes have harnessed innovative, touchscreen technology to offer a fresh approach to powerchair controls. LiNX is all about connectivity, using Bluetooth simply connect to PCs, laptops and Macs - keeping you in touch with your world.

Whether you're looking for basic drive only functionality or the latest touchscreen remotes, our LiNX Technology offers you smart controls, in tune with your life.

















### **Meet the Range – Up close and personal**

#### TDX® SP2 Ultra

Our elite powerchair features the Ultra Low Maxx seating for the more complex seating requirements and the revolutionary LiNX smart technology. The perfect mix of functionality, customisation and style.



- Suitable for a range of specialist controls including head controls
- Ultra Low Maxx seating system
- LiNX technology
- Range of accessories
- Max user weight 136kg
- Battery capacity 73.5Ah
- Speed 3, 6, 8, 10 & 12km/h
- Drive range up to 31km\*

#### TDX® SP2 Ultra NB

This narrow base version of TDX SP2 Ultra features a slimline power base ideal for confined environments whilst offering the same great features as the standard version.



- Suitable for a range of specialist controls including head controls
- Ultra Low Maxx seating system
- LiNX technology
- Range of accessories
- Max user weight 136kg
- Battery capacity 50Ah
- Speed 3, 6 & 10km/h
- Drive range up to 26km\*

#### TDX® SP2 & TDX® SP2 Max

Our standard version of the TDX SP2 range offers exceptional durability and performance. Versatile for use indoors and outdoors with Modulite seating and a wide range of options and accessories.



- Suitable for a range of specialist controls
- Modulite Seating system (TDX® SP2 only)
- LiNX technology
- Wide range of options and accessories
- Max user weight 180kg
- Battery capacity 73.5Ah
- Speed 3, 6, 8 & 12km/h\*
- Drive range up to 31km\*\*
- TDX SP2 Max version with Max seating available

<sup>\*</sup> TDX SP2 Max version; speeds 3, 6 & 8 km/h.

<sup>\*\*</sup> Dependent on chair configuration, terrain and individual use.

#### TDX® SP2 NB

Our narrow base version has exceptional drive performance with LiNX technology combined with safety and stability making this chair the popular choice for use both indoors and out. A highly functional and practical mobility solution with great styling.



- Suitable for a range of specialist controls
- Modulite Seating system
- LiNX technology
- Wide range of options and accessories
- Max user weight 150kg
- Battery capacity 50Ah
- Speed 3, 6 & 10km/h
- Drive range up to 26km\*

#### TDX® SP2 Low Rider

Featuring a drop down seat providing very low seat to ground heights of either 403 mm or 420 mm. This allows for practical driving positions and can easily access Wheelchair Accessible Vehicles (WAVs), fit under tables and align with a variety of transfer surfaces.



Ideal for young teens and petite adults offering smaller seat sizes, it easily adapts as a person grows to provide a truly customised fit. TDX SP2 Teen has the same great features on our stylish TDX power-base.



- 30cm riser unit with centre of gravity shift
- Modulite Seating system
- LiNX Technology
- Wide range of options and accessories
- Max user weight 160kg
- Battery capacity 60Ah or 73.5Ah
- Speed 3, 6, 8 & 12km/h
- Drive range up to 39km\*\*



- Matrx Mini-series seating or Modulite Seating system
- Bodypoint positioning belts
- LiNX Technology
- Wide range of options and accessories
- Max user weight 150kg
- Battery capacity 50Ah
- Speed 3, 6 & 10km/h
- Drive range up to 26km\*

<sup>\*</sup> Dependent on chair configuration, terrain and individual use.

<sup>\*\* 8</sup>km/h version.



Invacare LiNX is our insight inspired control system with advanced technology to provide a consistent driving experience. It's quick and simple to configure to suit your requirements.

#### **Built-in Bluetooth**

Both the REM400 and REM500 controls have built-in bluetooth which allows you to be connected to your computer, laptop and Mac or up to three devices simultaneously.





#### **REM500**

#### **Key Features**

- Display only version of REM400
- 3.5" interactive LCD touchscreen
- Allows alternative controls to be used to operate the powerchair.

#### **Touchscreen Remote**

The REM400 brings next generation technology to your Invacare powerchair through its impressive 3.5" colour touchscreen. This innovative remote is able to work alongside your physical and cognitive ability.

#### **REM400**

#### **Key features**

- Swipe or tap-only operation
- Direct access or step by step menu operation
- Adjustable brightness
- Configurable for left or right hand use
- · Glove mode setting

- Screen lock out ability
- Profiles and functions accessible via the touchscreen, buttons or joystick
- Mouse mover with on touchscreen mouse clicks
- Programmable Toggle Switches available





# 



This innovative app has been designed for use with LiNX controls to provide you with a quick and easy way to see how your powerchair is performing. No more second guessing when you plan your day!

MyLiNX connects your smart device to your LiNX powerchair control using Bluetooth. Once connected it sends data direct from the powerchair controls to your smartphone to show clear and simple graphics for battery charge, number of charging cycles, total drive time and much more.

In addition, the app highlights any faults on the powerchair with suggestions to fix them. If any issue can not be resolved, you simply contact your Service Provider directly through the app.

Now you can drive with confidence with just a quick tap of the app!


**Free Download** 











### **Colour your world**

Make it personal! Mix and match shroud colours and wheel rims giving you a choice of over 100 colour combinations.

Coloured wheel rims easily clip into the centre of the wheel and can be simply replaced or changed to suit your personal style.

What colour will you choose?

#### Rockstar Red

**Island Blue** 

**Lights Out Matt Black** 

**Invacare Blue** 

SpaceStation Silver

**Grape Jellybean** 

Wet Black



## **Specialist Controls & Accessories**

We offer a wide range of specialist controls to meet your specific requirements.

The LiNX system can be added to as technology develops or as and when your needs change. LiNX controls use Bluetooth which makes set up and programming quick, easy and in real-time to ensure drive profiles exactly match your requirements.

Talk to your local Invacare Retailer for further details.

www.invacarelinx.com











Invacare Limited
Unit 4 Pencoed Technology Park
Bridgend CF35 5AQ
United Kingdom
T: 01656 776222
F: 01656 776220
E: uk@invacare.com

Invacare Ireland Limited
Unit 5 Seatown Business Campus
Seatown Road
SWORDS
County Dublin - Ireland
T: +353 1 8107084
F: +353 1 8107085
E: Ireland@invacare.com

www.invacare.ie

The LINX word, mark and logos are registered trademarks owned by Dynamic Controls, and any use of LiNX such marks by Invacare is under license.

Apple App store logo is a trademark of Apple Inc., registered in the US and other countries.

Google Play and the Google Play logo are trademarks of Google LLC.

LiNX - UK 07/22 - Part number: 1630261













www.invacare.co.uk