



**Yes, you can:**

# Manual Retail Prescription Form

Tel: 01656 776222  
Email: [ordersuk@invacare.com](mailto:ordersuk@invacare.com)  
Online Spares available at: [www.invacare.co.uk](http://www.invacare.co.uk)

DECEMBER 2025  
Customer Ref:

<b>Customer Ref*</b>	<b>Account No*</b>
<b>Contact No*</b>	<b>Order No*</b>
<b>Delivery Address*</b>	<b>Order Qty*</b>
	<b>Order date*</b>

## Rea Clematis Pro Stock

### Standard Features

- Max user weight: 135kg (21.5 stone)
- Crash tested ISO 7176-Part 19
- 12" Transit wheels black solid tyres, 8" x 1 1/4" Black solid castors
- Choice of 3 widths 390mm, 440mm and 490mm with +25mm spacer on armrests
- Seat depth, stepless adjustable 420mm-480mm
- Seat height 450mm
- Flosshape seat cushion with multi stretch PU coated fabric TR26
- Plastic rigid backrest plate with Swingaway trunk supports
- Laguna back cushion with multi stretch PU coated fabric TR26
- Height adjustable push handles
- Assistant operated gas strut 30 Deg backrest recline and 25 Deg seat tilt with tilt lock
- Angle adjustable legrests with depth and angle adjustable footplates and heelstraps
- Concave armrests with armrest spacer adds additional 25mm of seat space
- Headrest, adjustable with handwheel
- Pearl Beige frame colour
- Quick release castor attachment with 200mm x 45mm medium hard black castors
- Foot operated brakes, Pelvic belt, Flip up Anti tip devices
- Transportation kit

**Standard Base price: £1,650**



Please make your selection with a cross in the box    ☐ Free of charge    ☒ standard feature    POA special request

Customer weight =                      KG

Code	Description	
<b>Rea Clematis Pro*</b>		
DV1673889	W 390mm (Small chassis) x D 430-500mm (Depth adjustable)	£1,650
DV1673890	W 440mm (Medium chassis) x D 430-500mm (Depth adjustable)	£1,650
DV1673891	W 490mm (Large chassis) x D 430-500mm (Depth adjustable)	£1,650

	TOTAL PRICE	
--	-------------	--

**Non Conforming Options:**



For a more complex customer requirement  
please refer to our Invacare Unique special order form.  
email [specialsuk@invacare.com](mailto:specialsuk@invacare.com)

Special Instructions  
Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):