



Yes, you can:

Tel: 01656 776222 Fax: 01656 776220
 Email: ordersuk@invacare.com
 Online Spares available at: www.invacare.co.uk

Manual Retail Prescription Form

NOVEMBER 2023
 Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

Rea Clematis Pro Stock

Standard Features

- Max user weight: 135kg (21.5 stone)
- Crash tested ISO 7176-Part 19
- 12" Transit wheels black solid tyres, 8" x 1 1/4" Black solid castors
- Choice of 3 widths 390mm, 440mm and 490mm with +25mm spacer on armrests
- Seat depth, stepless adjustable 420mm-480mm
- Seat height 450mm
- Flosshape seat cushion with multi stretch PU coated fabric TR26
- Plastic rigid backrest plate with Swingaway trunk supports
- Laguna back cushion with multi stretch PU coated fabric TR26
- Height adjustable push handles
- Assistant operated gas strut 30 Deg backrest recline and 25 Deg seat tilt with tilt lock
- Angle adjustable legrests with depth and angle adjustable footplates and heelstraps
- Concave armrests with armrest spacer adds additional 25mm of seat space
- Headrest, adjustable with handwheel
- Pearl Beige frame colour
- Quick release castor attachment with 200mm x 45mm medium hard black castors
- Foot operated brakes, Pelvic belt, Flip up Anti tip devices
- Transportation kit



Standard Base price: £1,650

Please make your selection with a cross in the box Free of charge standard feature POA special request

Customer weight = KG

Code	Description	
Rea Clematis Pro*		
DV1673889	W 390mm (Small chassis) x D 430-500mm (Depth adjustable)	£1,650
DV1673890	W 440mm (Medium chassis) x D 430-500mm (Depth adjustable)	£1,650
DV1673891	W 490mm (Large chassis) x D 430-500mm (Depth adjustable)	£1,650

TOTAL PRICE	
--------------------	--



For a more complex customer requirement please refer to our Invacare Unique special order form. email specialsuk@invacare.com or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: _____ Email: _____

Phone: _____ Fax: _____

Description of special request (please include relevant option code or part number):