

Cetus[®]

PREMIUM PERFORMANCE | ADVANCED COMFORT | COMPLETE INDEPENDENCE



The *Invacare* **Cetus** scooter has been specifically designed to offer a durable, premium product for those who desire complete independence to travel further afield. With high performance to support those with an active, independent lifestyle, the **Cetus** can travel up to 27 Miles with its high performance 100Ah batteries. Advanced comfort, adjustability and front and rear suspension, the **Cetus** will ensure a smooth ride across almost all terrains.

Complete independence

Longer distances

The **Cetus** is the ideal scooter for those who live life to the full. A maximum range of up to 27 miles, the **Cetus** is able to support individuals on those longer journeys. With speeds up to 8 mph and a high-low speed switch on the control panel to effortlessly switch between pavement and road speeds, the **Cetus** is ideal for all environments.

Suitable for all

A 22" wide Captains swivel seat, the **Cetus** is built to support all shapes and sizes. The **Cetus** can be configured easily with a built-in sliding seat mechanism, back angle adjustment and height adjustment. Maximum user weight of 35.5st, the **Cetus** is truly a diverse product for a range of individuals.

All terrain mobility scooter

Perfect for those rougher terrains, the **Cetus** is equipped with 15" tyres and front and rear suspension to provide a smoother ride for those looking for a more adventurous journey.











Available accessories

Front basket

USB charger adaptor

Oxygen cylinder holder

- Mobile phone holder
- Rollator holder
 - Protective storage cover
- Cane holder
- Front and rear bumpers
- Electric seat raiser

Technical data 🗸

For more information about this product, including the product's user manual, please visit www.invacare.co.uk

			L'S			<u>L</u>
	Backrest width ▼	Seat depth 🔻	Seat to floor height* ▼	Backrest height** ▼	Driving unit width ▼	Total length 🔻
Cetus	559 mm	510 mm	670 mm	500 mm	720 mm	1160 mm

	A		*			<u>L</u>
	Empty weight*** ▼	Max. user weight ▼	Battery capacity ▼	Turning diameter 🔻	Max. climbable obstacle height ▼	Ground clearance ▼
Cetus	178 kg	35.5 st (226kg)	2 x 12 V/100Ah	1860 mm	90 mm	110 mm

					*Seat post adjusted at lowest hole	
	5	КМ		W	** Without headrest height ***Including	
	Max. safe slope ▼	Max. driving distance ▼	Max. driving speed ▼	Motor capacity 🔻	batteries	
Cetus	9°	27miles (44km)	8 mph (12.8km/h)	800W	•	

Colour of shrouds -

Upholstery -





Invacare Limited Pencoed Technology Park - Pencoed Bridgend CF35 5AQ - United Kingdom Tel: +44 1656 776 222 Fax: +44 1656 776 220 E-mail: uk@invacare.com Sales Order E-mail: ordersuk@invacare.com www.invacare.co.uk Invacare Ireland Ltd Unit 5 - Seatown Business Campus Seatown Road - SWORDS County Dublin - Ireland Tel. +353 1 8107084 Fax +353 1 8107085 Email: Ireland@invacare.com www.invacare.ie





© 2020 Invacare International GmbH. All rights reserved. All information quoted is believed to be correct at time of print. Invacare reserves the right to alter product specifications without prior consultation. Cetus - UK - 05/2020 - 1632194