



Yes, you can:

Tel: 01656 776222
Email: ordersuk@invacare.com
www.invacare.co.uk

RETAIL PRESCRIPTION FORM

Matrx Elite E2 back
June 2025

Account Name :	Account No :
Address :	Order No :
Tel No:	Order Qty :
Delivery address :	Order Date :

Matrx Elite E2 (New)

Retail Order form
Maximum User Weight: 136kg

Standard Features:

3" Contour Depth
40° angle adjustment
4" depth adjustment
20° rotation

136kg (21 stone) user weight

All E2S backs 6 - 16" high include mini set, QR mounting hardware & hardware interface. All E2S backs 18" & 20" high include Easyset QR mounting hardware & easy set mounting interface

Important - Standard hardware is included



Height	Width							
	Inch	14	15	16	17	18	19	20
	cm	36	38	41	43	46	48	51
	FITS W/C	13-15"	14-16"	15-17"	16-18"	17-19"	18-19"	19-21"
	Inch	cm						
6	15	M	M	M	M	M	M	M
8	20	M	M	M	M	M	M	M
10	26	M	M	M	M	M	M	M
12	31	M	M	M	M	M	M	M
14	36	M	M	M	M	M	M	M
16	41	M	M	M	M	M	M	M
18	46	E	E	E	E	E	E	E
20	51	E	E	E	E	E	E	E


MiniSet Hardware

EasySet Hardware

Code	X	Description	Price	
Matrx Elite E2 Back		Dimension in cm (width/height) Hardware included	Fits wheelchair widths: 1" wider or 1" narrower	
		Description in brackets refers to previous Matrx Elite code	Inches	cm
1632513		ELITE E2 BACK 3615 (new size) includes std hardware	14	36
1632514		ELITE E2 BACK 3620 (new size) includes std hardware	14	36
1632515		ELITE E2 BACK 3626 (PBE1410NH) includes std hardware	14	36
1632516		ELITE E2 BACK 3631 (PBE1412NH) includes std hardware	14	36
1632517		ELITE E2 BACK 3636 (PBE1414NH) includes std hardware	14	36
1632518		ELITE E2 BACK 3641 (PBE1416NH) includes std hardware	14	36
1632519		ELITE E2 BACK 3646 (PBE1418NH) includes std hardware	14	36
1632520		ELITE E2 BACK 3651 (PBE1420NH) includes std hardware	14	36
1632521		ELITE E2 BACK 3815 (new size) includes std hardware	15	38
1632522		ELITE E2 BACK 3820 (new size) includes std hardware	15	38
1632523		ELITE E2 BACK 3826 (PBE1510NH) includes std hardware	15	38
1632524		ELITE E2 BACK 3831 (PBE1512NH) includes std hardware	15	38
1632525		ELITE E2 BACK 3836 (PBE1514NH) includes std hardware	15	38
1632526		ELITE E2 BACK 3841 (PBE1516NH) includes std hardware	15	38
1632527		ELITE E2 BACK 3846 (PBE1518NH) includes std hardware	15	38
1632528		ELITE E2 BACK 3851 (PBE1520NH) includes std hardware	15	38
1632529		ELITE E2 BACK 4115 (new size) includes std hardware	16	41
1632530		ELITE E2 BACK 4120 (new size) includes std hardware	16	41
1632531		ELITE E2 BACK 4126 (PBE1610NH) includes std hardware	16	41
1632532		ELITE E2 BACK 4131 (PBE1612NH) includes std hardware	16	41
1632533		ELITE E2 BACK 4136 (PBE1614NH) includes std hardware	16	41
1632534		ELITE E2 BACK 4141 (PBE1616NH) includes std hardware	16	41
1632535		ELITE E2 BACK 4146 (PBE1618NH) includes std hardware	16	41
1632536		ELITE E2 BACK 4151 (PBE1620NH) includes std hardware	16	41
1632537		ELITE E2 BACK 4315 (new size) includes std hardware	17	43
1632538		ELITE E2 BACK 4320 (new size) includes std hardware	17	43
1632539		ELITE E2 BACK 4326 (PBE1710NH) includes std hardware	17	43
1632540		ELITE E2 BACK 4331 (PBE1712NH) includes std hardware	17	43
1632541		ELITE E2 BACK 4336 (PBE1714NH) includes std hardware	17	43
1632542		ELITE E2 BACK 4341 (PBE1716NH) includes std hardware	17	43

1632543	ELITE E2 BACK 4346 (PBE1718NH) includes std hardware	17	43	£494.00
1632544	ELITE E2 BACK 4351 (PBE1720NH) includes std hardware	17	43	£494.00
1632545	ELITE E2 BACK 4615 (new size) includes std hardware	18	46	£494.00
1632546	ELITE E2 BACK 4620 (new size) includes std hardware	18	46	£494.00
1632547	ELITE E2 BACK 4626 (PBE1810NH) includes std hardware	18	46	£494.00
1632548	ELITE E2 BACK 4631 (PBE1812NH) includes std hardware	18	46	£494.00
1632549	ELITE E2 BACK 4636 (PBE1814NH) includes std hardware	18	46	£494.00
1632550	ELITE E2 BACK 4641 (PBE1816NH) includes std hardware	18	46	£494.00
1632551	ELITE E2 BACK 4646 (PBE1818NH) includes std hardware	18	46	£494.00
1632552	ELITE E2 BACK 4651 (PBE1820NH) includes std hardware	18	46	£494.00
1632553	ELITE E2 BACK 4815 (new size) includes std hardware	19	48	£494.00
1632554	ELITE E2 BACK 4820 (new size) includes std hardware	19	48	£494.00
1632555	ELITE E2 BACK 4826 (PBE1910NH) includes std hardware	19	48	£494.00
1632556	ELITE E2 BACK 4831 (PBE1912NH) includes std hardware	19	48	£494.00
1632557	ELITE E2 BACK 4836 (PBE1914NH) includes std hardware	19	48	£494.00
1632558	ELITE E2 BACK 4841 (PBE1916NH) includes std hardware	19	48	£494.00
1632559	ELITE E2 BACK 4846 (PBE1918NH) includes std hardware	19	48	£494.00
1632560	ELITE E2 BACK 4851 (PBE1920NH) includes std hardware	19	48	£494.00
1632561	ELITE E2 BACK 5115 (new size) includes std hardware	20	51	£494.00
1632562	ELITE E2 BACK 5120 (new size) includes std hardware	20	51	£494.00
1632563	ELITE E2 BACK 5126 (PBE2010NH) includes std hardware	20	51	£494.00
1632564	ELITE E2 BACK 5131 (PBE2012NH) includes std hardware	20	51	£494.00
1632565	ELITE E2 BACK 5136 (PBE2014NH) includes std hardware	20	51	£494.00
1632566	ELITE E2 BACK 5141 (PBE2016NH) includes std hardware	20	51	£494.00
1632567	ELITE E2 BACK 5146 (PBE2018NH) includes std hardware	20	51	£494.00
1632568	ELITE E2 BACK 5151 (PBE2020NH) includes std hardware	20	51	£494.00

Laterals



Fixed laterals

1632755	Easy fix laterals for Elite E2 backs 6.5" H, 5" Deep (Pair) can be fixed on inside or outside of back shell	£159.00
1545167	Standard Left Lateral (191mm long x 127mm high) PBLS-L	£131.00
1545168	Standard Right Lateral (191mm long x 127mm high) PBLS-R	£131.00
1545172	Small Left Lateral (140mm long x 76mm high) PBSS-L	£131.00
1545225	Small Left Lateral (140mm long x 76mm high) PBSS-R	£131.00
1545169	Offset Left Lateral (For 76mm height adjustment) PBOL-L	£131.00
1545170	Offset Right Lateral (For 76mm height adjustment) PBOL-R	£141.00

Swing away laterals

1551913	Swing Away Hardware for Left or Right Side SAHEB (Swings completely behind back. Does not interfere with flip-back arms)	£141.00
1544798	Offset Left Swing Away Hardware (Each) OSAH-L (Allows 76mm higher or lower pad position relative to attachment slot)	£141.00
1544799	Offset Right Swing Away Hardware (Each) OSAH-R (Allows 76mm higher or lower pad position relative to attachment slot)	£141.00
1544292	Extra Small Lateral Pad 82.5mm x 108mm (Each) EPAD-XS	£86.00
1544291	Small Lateral Pad 89mm x 133mm (Each) EPAD-S	£86.00
1544239	Medium Lateral Pad 95mm x 146mm (Each) EPAD-M	£86.00
1544238	Large Lateral Pad 108mm x 165mm (Each) EPAD-L	£86.00

PSIS support

1632759	E2 Backbone PSIS support	£50.00
---------	--------------------------	--------

PostureFlex Abdominal Support

1632760	Postureflex Abdominal Support H:6.5",Length 16" To 23"	£168.00
1632761	Postureflex Abdominal Support H:7",Length 20" To 32"	£168.00
1632796	Quad Release Tab for Removal of PostureFlex	£21.00

Harness interface

1632757	E2 Chest Harness Interface For E2 Backs 14/15"Wide	£101.00
1632758	E2 Chest Harness Interface For E2 Backs 16-20"Wide	£101.00

Please note: Shoulder harness back guides are also available on the Body point price list

Privacy flap attachment to fit all backs from 30.5cm to 50.5cm wide

1551799	Extra Small	26cm x 41cm	£49.00
1551800	Small	33cm x 41cm	£49.00
1551903	Medium	38cm x 41cm	£49.00
1551904	Large	43cm x 41cm	£49.00

Headrest



1634173	MatrX Loxx soft headrest pad kit	£269.00
---------	----------------------------------	---------

Headrest Adaptor

1632756	E2 Headrest adaptor - When required with back heights 12" - 13" - must be selected as an option	£21.00
---------	---	--------

TOTAL PRICE: