

Tel: 01656 776222

Email: ordersuk@invacare.com

Online spares available at www.invacare.co.uk

Account Name :	Account No :
Address :	Order No :
Tel No:	Order Qty :
Delivery address :	Order Date :

Matrx Posture Back

Retail order form
Maximum User Weight: 140kg

Features
5" Contour Depth
40° angle adjustment
4" depth adjustment
20° rotation
136kg (21 stone) user weight
Hardware not included



Standard MatRX Posture Back Backrest (no hardware)

Chest Width		Backrest Width	31cm	41cm	51cm	Fits W/ch Widths		With Extended Mounting Pins EMP*	
			(12")	(16")	(20")				
≤26cm	(10")	36cm (14")	MPB1412NH	MPB1416NH		33cm - 38cm	(13" - 15")	41cm	(16")
≤31cm	(12")	41cm (16")	MPB1612NH	MPB1616NH	MPB1620NH	38cm - 43cm	(15" - 17")	46cm	(18")
≤33cm	(13")	43cm (17")	MPB1712NH	MPB1716NH	MPB1720NH	41cm - 46cm	(16" - 18")	48cm	(19")
≤36cm	(14")	46cm (18")	MPB1812NH	MPB1816NH	MPB1820NH	43cm - 48cm	(17" - 19")	51cm	(20")
≤41cm	(16")	51cm (20")	MPB2021NH	MPB2016NH	MPB2020NH	48cm - 53cm	(19" - 21")	56cm	(22")

*not suitable for HD users (max. 136kg/300lb weight capacity)

Code	X	Description	Fits wheelchairs		Price
Matrx PB			Width	Height	
1544783		Matrx Posture Back MPB1412NH	330mm-380mm	310mm	£353.00
1544784		Matrx Posture Back MPB1416NH	330mm-380mm	410mm	£353.00
1544785		Matrx Posture Back MPB1612NH	380mm-430mm	310mm	£353.00
1515727		Matrx Posture Back MPB1616NH	380mm-430mm	410mm	£353.00
1515728		Matrx Posture Back MPB1620NH	380mm-430mm	510mm	£353.00
1544786		Matrx Posture Back MPB1712NH	410mm-460mm	310mm	£353.00
1515729		Matrx Posture Back MPB1716NH	410mm-460mm	410mm	£353.00
1515730		Matrx Posture Back MPB1720NH	410mm-460mm	510mm	£353.00
1544787		Matrx Posture Back MPB1812NH	430mm-480mm	310mm	£353.00
1515731		Matrx Posture Back MPB1816NH	430mm-480mm	410mm	£353.00
1515732		Matrx Posture Back MPB1820NH	430mm-480mm	510mm	£353.00
1544788		Matrx Posture Back MPB2012NH	480mm-510mm	310mm	£353.00
1518131		Matrx Posture Back MPB2016NH	480mm-510mm	410mm	£353.00
1518132		Matrx Posture Back MPB2020NH	480mm-510mm	510mm	£353.00

Hardware



Easy Set hardware



Fixed hardware



Split hardware



Heavy Duty hardware

1545275	Easy Set mounting hardware (pair) PBMH	Fits 19mm, 22mm and 25mm back canes	£182.00
1545276	Easy Set mounting hardware (pair) PBMHQ	Fits 29mm back canes	£182.00
1545277	Fixed Space Saver Mounting Hardware FSSH		£200.00
1545278	Split Mounting hardware SMH		£182.00
1544776	Heavy Duty Easy Set Mounting Hardware HDMH	Fits 19mm, 22mm and 25mm back canes	£260.00

Harness interface

CHIPB	Harness Interface to allow attachment of various chest supports without drilling of the back shell	Width, angle and height adjustability for precise orientation of shoulder straps	£29.00
-------	--	--	--------

Privacy flap attachment to fit all backs from 30.5cm to 50.5cm wide

1551799	Extra Small	26cm x 41cm	£49.00
1551800	Small	33cm x 41cm	£49.00
1551903	Medium	38cm x 41cm	£49.00
1551904	Large	43cm x 41cm	£49.00

Spares

EMP	Extended Mounting Pins	Allows Matrx backs to be installed on wider wheelchairs	£50.00
-----	------------------------	---	--------

Headrest

1634173	Matrx Loxx soft headrest pad kit		£269.00
---------	----------------------------------	--	---------

Headrest Receiver to suit backrest width:

1544189	Small Adapter for backs 12" wide EHRA-S		£45.00
1544188	Medium Adapter for backs 13-17" wide EHRA-M		£45.00
1544187	Large Adapter for backs 18-20" wide EHRA-L		£45.00

TOTAL PRICE:



Special Instructions - Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code):