



Yes, you can:

Tel: 01656 776222

Email: ordersuk@invacare.com

Online spares available at www.invacare.co.uk

ROBIN PRESCRIPTION FORM

April 2025

Account Name :	Account No :
Address :	Order No :
Tel No:	Order Date :
Delivery address :	

Robin Retail Order Form

Please indicate details of your installation on the bottom page

Step 1. Choose a Robin, Trolley and Charger

Product code	Qty	Motor Unit	RRP	Total
1588287		Invacare Robin (manual version) N.B Remember to select a trolley and charger	1411	£ -

	Qty	Trolley - The below trolleys are approved for use with the tracks as indicated below	RRP	Total
1444008-9010		Trolley for Robin motor EC Track in white	78	£ -
1485183-9010		Trolley for Robin motor room-to room in white	76	£ -
1451871-9010		Trolley for Guldmann tracks for min / maxi / jumbo	76	£ -
1446450-9010		Trolley for Liko track	101	£ -
1505320-9010		Trolley for Chiltern / Prism / Handicare rails	76	£ -
1530370-9010		Trolley for Oxford Easy Track Gantry / Arjo / BHM rails	141	£ -
1530293-9010		Trolley for Ergolet Travers Rail 135 mm rails	141	£ -
1525009-9010		Trolley for Ergolet Loftpro Rail 75 / 103 mm rails	141	£ -
1530434-9010		Trolley for ETAC / LIC Rails	89	£ -
60132135-9010		Trolley for Handicare MilkyWay rails	141	£ -
60132136-9010		Trolley for Etac / Molift rails	141	£ -
	Qty	Charger Set	RRP	Total
1493471-7035		Charger set (1.2/3m) IP67 with holder and bracket	203	£ -

Step 2. Now choose your track profile and accessories

White	Qty	EC Track - S-profile	Total Qty of Track	Length of track 1	Length of track 2	Length of track 3	Length of track 4	RRP	Total
1450426-9010		2m in Length						160	£ -
1450428-9010		3m in Length						238	£ -
1450430-9010		4m in Length						315	£ -
1450432-9010-SPECIAL		5m in Length						364	£ -
White	Qty	EC Track - M-profile	Total Qty of Track	Length of track 1	Length of track 2	Length of track 3	Length of track 4	RRP	Total
1446406-9010		3m in Length						353	£ -
1446408-9010		4m in Length						466	£ -
1446409-9010-SPECIAL		4.8m in Length						558	£ -
White	Qty	EC Track - L-profile	Total Qty of Track	Length track 1	Length track 2	Length of track 3	Length track 4	RRP	Total
1446413-9010-SPECIAL		5m in Length						705	£ -
1446415-9010-SPECIAL		6m in Length						790	£ -
1446417-9010-SPECIAL		7m in Length						900	£ -
1446418-9010-SPECIAL		7.8m in Length						1225	£ -
White (Standard)	Qty	Curves radius 400 for small track						RRP	Total
1450465-9010		30 degree, right						510	£ -
1450466-9010		30 degree, left						510	£ -
1450496-9010		45 degree, right						510	£ -
1450497-9010		45 degree, left						510	£ -
1450498-9010		60 degree, right						510	£ -
1450499-9010		60 degree, left						510	£ -
1450500-9010		90 degree, right						542	£ -
1450501-9010		90 degree, left						542	£ -
Fixtures - Standard	Qty	White						RRP	Total
1444045-9010		Curve 30 mm						15	£ -
1444046-9010		Standard Straight Fixture 70mm						17	£ -
1450509-9010		Standard Connection 120mm m/lock bolt						33	£ -
1450508-9010		Standard Straight 70mm m/lock bolt						28	£ -
Fixture - Quick	Qty	White						RRP	Total
1446295-9010		Straight 70 mm						20	£ -
1446296-9010		Curve 30 mm						23	£ -
1446297-9010		Connection 120 mm						23	£ -
White	Qty	Pendels						RRP	Total
1449402-9010		20-30 cm, adjustable, white						60	£ -
1449403-9010		28-46 cm, adjustable, white						56	£ -
1449404-9010		44-74 cm, adjustable, white						58	£ -

	Qty	Pendels (Special Order Item)	Length Required	RRP	Total
White/Alu					
1449401-9010-SPECIAL		Fixed, state length		71	£ -
1450291-9010-SPECIAL		Pendel with Tube, Fixed, White - State Length		121	£ -
White	Qty	Wall Fixture Kit		RRP	Total
1449669-9010		Standard, M- profile		84	£ -
1449670-9010		Standard, L-profile		113	£ -
1449671-9010		Diagonal, M-profile		84	£ -
1449672-9010		Diagonal, L- profile		141	£ -
1451647-9010		Side Hung Wall Fixture		45	£ -
	Qty	Floor and wall support kit		RRP	Total
1450557-9010		2,6 m for M profile		206	£ -
1450559-9010		3,0 m for M profile		218	£ -
1450558-9010		2,6 m for L profile		206	£ -
1450560-9010		3,0 m for L profile		218	£ -
1449334-9010		2,6 m for M profile, diagonal		275	£ -
1449336-9010		3,0 m for M profile, diagonal		240	£ -
1449335-9010		2,6 m for L profile, diagonal		309	£ -
1449337-9010		3,0 m for L profile, diagonal		281	£ -
1441488-9010		Feet to wall support, std		6	£ -
1441416-9010		Feet for wall support, skirting board		9	£ -
	Qty	Insert Bar for joining Track		RRP	Total
1446357-9010		Insert Bar for joining M Profile track		34	£ -
1446358-9010-SPECIAL		Insert Bar for joining L Profile track		74	£ -
	Qty	Traverse trolley - Pair (standard)		RRP	Total
1446422-9010		White Robin Traverse trolley set - partly between hung/under hung		235	£ -
1449667-9010		White Robin Traverse trolley set - between hung		140	£ -
Transit Coupling	Qty	Motor Box Track		RRP	Total
1491160-9010-SPECIAL		White		1700	£ -

N.B Remember to order two pendels per coupling

Placement of Motor	The motor will be placed in		Room 1	Room 2		
	Please tick box					
White	Qty	Fixtures for transitchoupling (Special order)	Qty Room 1	Qty Room 2	RRP	Total
1478763-9010-SPECIAL		Transit Coupling Lock Fixture Traverse Track S profile			290	£ -
1483864-9010-SPECIAL		Transit Coupling Lock Fixture M/L profile			296	£ -
1483862-9010-SPECIAL		Transit Coupling adjusting fixture S profile (parallel track profile)			59	£ -
1485118-9010-SPECIAL		Transit Coupling adjusting fixture M/L profile (parallel track profile)			52	£ -
		Standard Traverse - Description				
		Under hung				
		Between hung				
		Special Order Traverse - Description				
		Partly between hung				
	Qty	EC-Track Switch			RRP	Total
1451885-9010-SPECIAL		White			1213	£ -
White	Qty	Standard Shims			RRP	Total
1446327-9010		Shim, 1 mm for 30 mm for quick ceiling bracket			8	£ -
1446392-9010		Shim, 1 mm for 70 mm for quick ceiling bracket			8	£ -
1446330-9010		Shim, 1 mm for 120 mm for quick ceiling bracket			8	£ -
1446328-9010		Shim, 3 mm for 30 mm for quick ceiling bracket			8	£ -
1446393-9010		Shim, 3 mm for 70 mm for quick ceiling bracket			8	£ -
1446331-9010		Shim, 3 mm for 120 mm for quick ceiling bracket			8	£ -
1446329-9010		Shim, 5 mm for 30 mm for quick ceiling bracket			8	£ -
1446394-9010		Shim, 5 mm for 70 mm for quick ceiling bracket			8	£ -
1446332-9010		Shim, 5 mm for 120 mm for quick ceiling bracket			8	£ -
1446333-9010		Shim, 1 mm for 30 mm for standard ceiling bracket			8	£ -
1446395-9010		Shim, 1 mm for 70 mm for standard ceiling bracket			8	£ -
1449118-9010		Shim, 1 mm for 120 mm for standard ceiling bracket			8	£ -
1446334-9010		Shim, 3 mm for 30 mm for standard ceiling bracket			8	£ -
1446396-9010		Shim, 3 mm for 70 mm for standard ceiling bracket			8	£ -
1449119-9010		Shim, 3 mm for 120 mm for standard ceiling bracket			13	£ -
1449117-9010		Shim, 5 mm for 30 mm for standard ceiling bracket			13	£ -
1446397-9010		Shim, 5 mm for 70 mm for standard ceiling bracket			9	£ -
1449120-9010		Shim, 5 mm for 120 mm for standard ceiling bracket			13	£ -
White	Qty	End Cover (Each)			RRP	Total
1441420-9010		S-profile			5	£ -
1441421-9010		M-profile			5	£ -
1441422-9010-SPECIAL		L-profile			5	£ -
	Qty	End Stop (Pair)			RRP	Total
1449563-9010		White			35	£ -
	Qty	Accessories			RRP	Total
1446455-7012		Strap, pair, incl. hooks and locks			174	£ -
1493630-7035		Charger for External Battery			114	£ -
1446269-7012		Flexi Scope, Bariatric			289	£ -
1492967-7012		2 Point Spreader Bar 350 mm			199	£ -
1492968-7012		2 Point Spreader Bar 450 mm			248	£ -
1451870-7012		4 Point Spreader Bar 450 mm			400	£ -
1493132-SPECIAL		Spreader Bar with Weighing Scales			2837	£ -
1493227-7012		Transport Handle			114	£ -
					TOTAL	
					£	-
					TOTAL	

If you have selected items from the light grey shaded boxes, please use this page to indicate details of your installation
This will allow Invacare to ensure that track is cut according to the layout of the installation

State length between parallel tracks, measured from centre of track (T)
State traverse track's profile

Profile T=

cm

Mounting	
Ceiling	p
Wall	p
Wall/floor	p
Pendle	p
Portal	p
Other	p

Make these marks on your system;

P = Parking/charging station
T = Distance between parallel tracks (traverse system only)

	A	B	C	D
Wall material:				
Length of track:				
Track profile:				

Material	Length
Concrete	
Wood	
Light concrete	
Gypsum	
False ceiling	

A diagram of a square with vertices labeled A, B, C, and D. The square is oriented with its sides parallel to the axes. Vertex A is at the top-left corner, B is at the top-right corner, C is at the bottom-right corner, and D is at the bottom-left corner. Each vertex is marked with a small circle containing its respective letter.

Name of technical contact:

Telephone no.of technical contact:

Special Requirements:

--