



Yes, you can:

Tel: 01656 776222 Fax: 01656 776220
 Email: ordersuk@invacare.com
 Online Spares available at: www.invacare.co.uk

January 2022
 Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

Colibri 3w

Standard Features

- Colibri 3 wheeled scooter, 8km/h (6mph)
- "Lite-Lock" one handed, take-apart for easy transportation
- 5 colours (Sapphire Blue, Jasper Red, Marble White, Diamond Silver, Jade Green)
- Stylish swivel seat with width adjustable arms
- Adjustable seat height
- Delta handlebar with tiller adjustment in 3 steps
- 200mm x 50mm solid Black tyres
- Front basket in Black
- 2 x 18Ah batteries
- 2 year warranty (12 months on batteries)
- Max User Weight 136Kg (21.4 stone)



Please make your selection with a cross in the box ○ Free of charge ◆ standard feature POA special request

Customer weight = KG

Code	Description	
BASE PRODUCT*		
000000	Colibri 3 Wheeled	○
SPECIFICATION*		
1605406	210mm x 65mm Wheels, 18Ah	£1,571
COLOURED SHROUDS*		
SP1604351	Sapphire Blue	○
SP1604313	Jasper Red	○
SP1604312	Marble White	○
SP1604352	Diamond Silver	○
SP1604314	Jade Green	○
SEAT SUSPENSION		
SP1577212	Low, seat height 425mm	£49
SP1577213	Mid, seat height 450mm	£49
SP1577214	High, seat height 475mm	£49

SPARE / ADDITIONAL SHROUD KITS		
SP1604351		Sapphire Blue £37
SP1604313		Jasper Red £37
SP1604312		Marble White £37
SP1604352		Diamond Silver £37
SP1604314		Jade Green £37
ACCESSORIES		
SP1543153		Cylinder with 2 Personal Keys £24
SP1509434		Storage Cover £86
SP1507388		Accessory Bracket £15
SP1528274		Crutch & Cane Holder £45
SP1572971		Oxygen Bottle Holder £56
SP1505727		Rear Basket £49

TOTAL PRICE	
--------------------	--

Non Conforming Options:



For a more complex customer requirement please refer to our Invacare Unique special order form. email specialsuk@invacare.com or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: _____ Email: _____

Phone: _____ Fax: _____

Description of special request (please include relevant option code or part number):