



Yes, you can:

Tel: 01656 776222
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

Power Retail Prescription Form

01 DECEMBER 2025
Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

Esprit Action

Standard Features

- Max User Weight 125kg (19.7 Stone)
- Crash tested ISO7176 - 19*
- Double crossbar
- Seat widths 380mm to 605mm (15" - 24")
- Seat depths 400mm, 450mm or 500mm (16", 18" or 20")
- Seat heights 460mm or 485mm (18" or 19")
- Nylon seat upholstery
- Fixed 7° 430mm backrest
- Lumbar Support (+1 STRAP) BI Material Padded backrest
- 80° swing in / swing out hangers
- Flip up footplates including heel strap
- Removable full length armrests
- Choice of 4 standard colours (with 9 special colour options)
- Anti tippers
- Buckle belt
- *
- * Battery MUST be removed when used as a seat in a vehicle



Standard Base price: £3,706



Please make your selection with a cross in the box ☐ Free of charge ☒ standard feature POA special request

Customer weight = KG

Code	Description	
BASE PRODUCT*		
SPE0060	ESPRIT ACTION ADULT (DOUBLE CROSSBAR)	£3,706
CHAIR TYPE*		
SPESPEC	SPECIAL	POA
SPE0100	BASE MODEL ESPRIT ACTION	◆
CHASSIS*		
SPE0210	DOUBLE CROSSBAR (SW 380 MM TO 605 MM / SD FROM 400 MM TO 500 MM) USER MAX WEIGHT = 125Kg	◆

BRACKET TYPE*			
SPE0300		STANDARD SWING AWAY BRACKET	◆
SPE0310		ALBER E35 SWIVEL ARM	£360
USER REMOTE*			
SPE0400		USER REMOTE CONTROL RH	○
SPE0410		USER REMOTE CONTROL LH	○
DUAL REMOTE			
SPE0440		INTUITIVE DUAL REMOTE CONTROL RH	£853
SPE0450		INTUITIVE DUAL REMOTE CONTROL LH	£853
SPE0460		INTUITIVE SINGLE DUAL REMOTE CONTROL (ONE HANDLE) (ADJUSTABLE RH TO LH SIDE)	£608
SPEED*			
SPE0500		6 KM/H	◆
BATTERY*			
SPE0600		BATTERY 7,65AH/36V (NON SPILLABLE LEAD ACID BATTERY)	◆
POWER OPTIONS			
SPE0610		ADDITIONAL BATTERY 7,65AH/36V	£271
SPE0630		REMOTE PROTECTOR	£62
SEAT WIDTH*			
SPE1040		380 MM SEAT WIDTH (NC 1250)	○
SPE1050		405 MM SEAT WIDTH (NC 1250)	○
SPE1060		430 MM SEAT WIDTH	○
SPE1070		455 MM SEAT WIDTH	○
SPE1080		480 MM SEAT WIDTH	○
SPE1090		505 MM SEAT WIDTH	○
SPE1100		555 MM SEAT WIDTH (REINFORCED SEAT UPHOLSTERY) (NC 2000, 2010, 2030, 2200, 2210, 2220, 2230, 5020, 5040, 5080, 8200)	£103
SPE1110		605 MM SEAT WIDTH (REINFORCED SEAT UPHOLSTERY) (NC 2000, 2010, 2030, 2200, 2210, 2220, 2230, 5020, 5040, 5060, 5070, 5080, 8200)	£103
SEAT DEPTH*			
SPE1210		400 MM SEAT DEPTH (NC 5090)	○
SPE1220		400 TO 450 MM ADJUSTABLE SEAT DEPTH SET AT 430 MM	○
SPE1230		450 MM SEAT DEPTH (NC 5080, 5090)	○
SPE1240		450 TO 500 MM ADJUSTABLE SEAT DEPTH (NC 5080)	£30
SPE1250		500 MM SEAT DEPTH (NC 5080)	○
SEAT HEIGHT*			
SPE1300		460 MM SEAT HEIGHT	○
SPE1310		485 MM SEAT HEIGHT	○
SEAT UPHOLSTERY TYPE*			
SPE1410		NYLON WITH LOOP BANDS STRIPS SEAT UPHOLSTERY (DOUBLE CROSSBAR ONLY)	◆

BACKREST TYPE*			
SPE2000		FIXED 7° BACKREST 2 HEIGHTS 400/430 MM (DELIVERED 430 MM)	◆
SPE2010		FIXED 7° BACKREST 510 MM HEIGHT	○
SPE2020		FOLDING BACKREST 430 MM HEIGHT	£32
SPE2030		FOLDING BACKREST 510 MM HEIGHT	£32
SPE2040		ANGLE ADJUSTABLE AND FOLD DOWN (-15° TO + 15°) FIXED 7° HEIGHT 425-450 MM / 17"-18" SET AT 425 MM / 17" (2 POSITIONS; 25 MM / 1" STEPS THROUGH UPHOLSTERY)	£84
SPE2050		RECLINING/FOLDING (GAS SPRING) 0/30° BACKREST 3 HEIGHTS 390/410/430 MM (DELIVERED 430 MM)	(NC 2210, 2220, 2230) £147
SPE2060		RECLINING (RATCHET) 0/30° BACKREST 3 HEIGHTS 390/410/430 MM (DELIVERED 430 MM)	£124
SPE2070		RECLINING (RATCHET) 0/30° BACKREST 3 HEIGHTS 470/490/510 MM (DELIVERED 510 MM)	£124
SPE2080		RECLINING (GAS SPRING) 0/30° BACKREST 3 HEIGHTS 390/410/430 MM (DELIVERED 430 MM)	£51
SPE2090		RECLINING (GAS SPRING) 0/30° BACKREST 3 HEIGHTS 470/490/510 MM (DELIVERED 510 MM)	£64
BACKREST UPHOLSTERY TYPE*			
SPE2200		MATRX ELITE E2 BACK WITH HARDWARE 405 MM BACK HEIGHT ADDITIONAL DELIVERY TIME + 5 DAYS	(NC 2310, 2320, 5140) £557
SPE2210		MATRX ELITE E2 BACK WITH HARDWARE 455 MM BACK HEIGHT ADDITIONAL DELIVERY TIME + 5 DAYS	(NC 2310, 2320, 5140) £557
SPE2220		MATRX ELITE DEEP E2 BACK WITH HARDWARE 405 MM BACK HEIGHT ADDITIONAL DELIVERY TIME + 5 DAYS	(NC 2310, 2320, 5140) £647
SPE2230		MATRX ELITE DEEP E2 BACK WITH HARDWARE 455 MM BACK HEIGHT ADDITIONAL DELIVERY TIME + 5 DAYS	(NC 2310, 2320, 5140) £647
SPE2300		LUMBAR SUPPORT (+1 STRAP) BI MATERIAL PADDED COVER	◆
SPE2310		CONTOUR ADJUSTABLE (4+1 STRAPS) BI MATERIAL COVER	○
SPE2320		CONTOUR ADJUSTABLE (FULL STRAPS) BI MATERIAL COVER	○
HANGERS RH*			
SPE3010		80° LONG, UL FROM 340 TO 470MM	○
SPE3020		ANGLE ADJUSTABLE (0° TO 75°) RH SWING IN/SWING OUT HANGER, UL FROM 340 TO 500MM (INCLUDING CALFPAD)	(NC 5150) £88
SPE3030		AMPUTEE LEGREST RH	(NC 3210, 3220) £129
HANGERS LH*			
SPE3110		80° LONG, UL FROM 340 TO 470MM	○
SPE3120		ANGLE ADJUSTABLE (0° TO 75°) LH SWING IN/SWING OUT HANGER, UL FROM 340 TO 500MM (INCLUDING CALFPAD)	(NC 5150) £88
SPE3130		AMPUTEE LEGREST LH	(NC 3300, 3310) £129
FOOTPLATES RH*			
SPE3200		FLIP UP RH FOOTPLATE INCLUDING HEEL STRAP	○
SPE3210		FLIP UP ANGLE AND DEPTH ADJUSTABLE RH FOOTPLATE, INCLUDING HEEL STRAP	£42
SPE3220		FLIP UP ONE PIECE FOOTPLATE, DEPTH & ANGLE ADJUSTABLE	(NC 3310) £76
SPZ7000		AMPUTEE LEGREST OR FOOTBOARD SELECTED	(NC 3010, 3020) ○

FOOTPLATES LH*			
SPE3300		FLIP UP LH FOOTPLATE, INCLUDING HEEL STRAP	○
SPE3310		FLIP UP ANGLE AND DEPTH ADJUSTABLE LH FOOTPLATE, INCLUDING HEEL STRAP	£42
SPZ7001		AMPUTEE LEGREST OR FOOTBOARD SELECTED	○
ARMRESTS RH*			
SPE4000		REMOVABLE RH DESK, ADJ IN HEIGHT (235-280 for FSTF 460 & 210-280 for FSTF 485MM) DEPTH (+35 MM FRONTWARD) & WIDTH (+25MM OUTWARDS)	○
SPE4010		REMOVABLE RH FULL LENGTH, ADJ IN HEIGHT (210 TO 280 MM) DEPTH (+35 MM BACKWARD) & WIDTH (+25MM OUTWARDS)	◆
SPE4020		REMOVABLE RH ARMPAD HEMI, BLACK PU TR26, MULTIPOSITION AUTOLOCK SAFETY DEVICE INCLUDED, ADJ IN HEIGHT (210 TO 345 MM) DEPTH (35 MM FRONTWARD) & WIDTH (+25MM OUTWARDS)	£32
ARMRESTS LH*			
SPE4100		REMOVABLE LH DESK, ADJ IN HEIGHT (235-280 for FSTF 460 & 210-280 for FSTF 485MM) DEPTH (+35 MM FRONTWARD) & WIDTH (+25MM OUTWARDS)	○
SPE4110		REMOVABLE LH FULL LENGTH, ADJ IN HEIGHT (210 TO 280 MM) DEPTH (35 MM BACKWARD) & WIDTH (+25MM OUTWARDS)	◆
SPE4120		REMOVABLE LH ARMPAD HEMI, BLACK PU TR26, MULTIPOSITION AUTOLOCK SAFETY DEVICE INCLUDED, ADJ IN HEIGHT (210 TO 345 MM) DEPTH (35 MM FRONTWARD) & WIDTH (+25MM OUTWARDS)	£32
SEAT CUSHION*			
SPE5010		BASIC 50 MM FOAM CUSHION - 50 MM HEIGHT, NYLON & ANTI-SLIP COVER	○
SPE5020		MATRX FLO-TECH LITE CUSHION - 78 MM HEIGHT	£117
SPE5040		MATRX CONTOUR VISCO NG CUSHION - 100 MM HEIGHT MAXIMUM USER WEIGHT 100 KG	£164
SPE5060		MATRX VI CUSHION - 82 MM HEIGHT, ADDITIONAL DELIVERY TIME + 5 DAYS	£359
SPE5070		MATRX LIBRA CUSHION - 70 MM HEIGHT, ADDITIONAL DELIVERY TIME + 5 DAYS	£387
SEAT CUSHION DEPTH*			
SPE5080		FOR 400 MM SEAT DEPTH	○
SPE5085		FOR 450 MM SEAT DEPTH	○
SPE5090		FOR 500 MM SEAT DEPTH	○
SEATING OPTIONS			
SPE5100		NECKREST BLACK PU TR26	£127
SPE5110		HEADREST BLACK PU TR26	£153
SPE5120		HEADREST WITH CHEEK SUPPORT BLACK PU TR26 (INCLUDING BACKBRACE EXCEPT FOR LAGUNA)	£177
SPE5140		LATERAL POSITIONING PADS (PAIR)	£84
SPE5150		ADJUSTABLE CALF STRAP	○
BACK BRACE			
SPE5170		BACK BRACE FIXED	○

COLOURS*			
SPE6000	☐	ANTHRACITE BLACK	○
SPE6015	☐	DARK PURPLE	○
SPE6035	☐	NAVY BLUE	○
SPE6040	☐	SLATE GREY	○
SWO6925	☐	WATER BLUE - 3 week lead time	○
SWO0806	☐	Electric Green - 3 week lead time	£164
SWO6033	☐	Red Matt - 3 week lead time	£164
SWO0607	☐	Dark Red - 3 week lead time	£164
SWO0515	☐	Pearl Grey - 3 week lead time	£164
SWO0647	☐	Intense Yellow - 3 week lead time	£164
SWO0648	☐	Vibrant Orange - 3 week lead time	£164
SWO0651	☐	Polar Blue - 3 week lead time	£164
SWO0652	☐	Candy Pink - 3 week lead time	£164
SPE6050	☐	SPECIAL COLOUR (ON DEMAND ONLY)	POA
ACCENT COLOUR*			
SPE6100	☐	BLACK	○
SPE6120	☐	LIME	○
FRONT WHEELS*			
SPE7000	☐	200 X 50 MM (8" X 2") BLACK SOLID TYRE	£47
SPE7010	☐	200 X 50 MM (8" X 2") BLACK PNEUMATIC TYRE	◆
REAR WHEELS*			
SPE7100	☐	355 X 70 MM (14" X 27") BLACK SOLID WITH PUNCTURE PROOF TYRE	£82
SPE7110	☐	355 X 70 MM (14" X 27") BLACK PNEUMATIC TYRE	◆
ANTI TIPPER*			
SPE8000	☐	REAR ANTI-TIPPERS, FIXED, PAIR	◆
SPE8010	☐	REAR ANTI-TIPPERS REMOVABLE WITH JACK UP, PAIR	£276
LAP BELT*			
SPE8020	☐	POSTURE BELT	◆
TRANSPORTATION BRACKET*			
SPE8030	☐	TRANSPORTATION BRACKET	◆
OTHER OPTIONS			
SPE8100	☐	CANE HOLDER	£20
SPE8200	☐	ACRYLIC TRAY FLIP UP	£308
TOTAL PRICE			

Non Conforming Options:	
SPE1250	500 MM SEAT DEPTH
SPE2000	FIXED 7° BACKREST 2 HEIGHTS 400/430 MM (DELIVERED 430 MM)
SPE2010	FIXED 7° BACKREST 510 MM HEIGHT
SPE2030	FOLDING BACKREST 510 MM HEIGHT
SPE2200	MATRX ELITE E2 BACK WITH HARDWARE 405 MM BACK HEIGHT ADDITIONAL DELIVERY TIME + 5 DAYS
SPE2210	MATRX ELITE E2 BACK WITH HARDWARE 455 MM BACK HEIGHT ADDITIONAL DELIVERY TIME + 5 DAYS
SPE2220	MATRX ELITE DEEP E2 BACK WITH HARDWARE 405 MM BACK HEIGHT ADDITIONAL DELIVERY TIME + 5 DAYS
SPE2230	MATRX ELITE DEEP E2 BACK WITH HARDWARE 455 MM BACK HEIGHT ADDITIONAL DELIVERY TIME + 5 DAYS
SPE2310	CONTOUR ADJUSTABLE (4+1 STRAPS) BI MATERIAL COVER
SPE2320	CONTOUR ADJUSTABLE (FULL STRAPS) BI MATERIAL COVER
SPE3010	80° LONG, UL FROM 340 TO 470MM
SPE3020	ANGLE ADJUSTABLE (0° TO 75°) RH SWING IN/SWING OUT HANGER, UL FROM 340 TO 500MM (INCLUDING CALFPAD)
SPE3210	FLIP UP ANGLE AND DEPTH ADJUSTABLE RH FOOTPLATE, INCLUDING HEEL STRAP
SPE3220	FLIP UP ONE PIECE FOOTPLATE, DEPTH & ANGLE ADJUSTABLE
SPE3300	FLIP UP LH FOOTPLATE, INCLUDING HEEL STRAP
SPE3310	FLIP UP ANGLE AND DEPTH ADJUSTABLE LH FOOTPLATE, INCLUDING HEEL STRAP
SPE5020	MATRX FLO-TECH LITE CUSHION - 78 MM HEIGHT
SPE5040	MATRX CONTOUR VISCO NG CUSHION - 100 MM HEIGHT MAXIMUM USER WEIGHT 100 KG
SPE5060	MATRX VI CUSHION - 82 MM HEIGHT, ADDITIONAL DELIVERY TIME + 5 DAYS
SPE5070	MATRX LIBRA CUSHION - 70 MM HEIGHT, ADDITIONAL DELIVERY TIME + 5 DAYS
SPE5080	FOR 400 MM SEAT DEPTH
SPE5090	FOR 500 MM SEAT DEPTH
SPE5140	LATERAL POSITIONING PADS (PAIR)
SPE5150	ADJUSTABLE CALF STRAP
SPE8200	ACRYLIC TRAY FLIP UP



For a more complex customer requirement,
please email: orders@invacare.com.au

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):