



**Yes, you can:**

# Powerchairs Retail Prescription Form

Tel: 01656 776222 Fax: 01656 776220  
 Email: ordersuk@invacare.com  
 Online Spares available at: www.invacare.co.uk

APRIL 2024  
 Customer Ref:

<b>Customer Ref*</b>	<b>Account No*</b>
<b>Contact No*</b>	<b>Order No*</b>
<b>Delivery Address*</b>	<b>Order Qty*</b>
	<b>Order date*</b>

## Esprit Action Junior

### Standard Features

- Max User Weight : 80kg (12.6 Stone)
- Single crossbar (Width Adjustable)
- Seat widths 305mm to 380mm 12" - 15"
- Seat depth 325mm (13") to 375mm (15") adjustable (set at 350mm)
- Seat height 460mm (18")
- Nylon seat upholstery
- FIXED 7° BACKREST 430mm backrest
- Standard nylon sling type backrest upholstery
- 80° swing in / swing out legrest hangers
- Flip up angle and depth adjustable footplates
- Height adjustable, removable arm armrests
- Choice of 3 standard colours (with 10 Special colour options)

**Standard Base price: £3,706**



Please make your selection with a cross in the box    ○ Free of charge    ◆ standard feature    POA special request

Customer weight =                      KG

Code	Description	
<b>BASE PRODUCT*</b>		
SPD0050	ESPRIT ACTION JUNIOR (SINGLE CROSSBAR) (NC 2300)	£3,706
<b>CHAIR TYPE*</b>		
SPDSPEC	SPECIAL	POA
SPD0100	BASE MODEL ESPRIT ACTION	◆
<b>CHASSIS*</b>		
SPD0200	SINGLE CROSSBAR (SW 305 MM TO 380 MM / SD 325 MM TO 375 MM) USER MAX WEIGHT = 80Kg	◆

<b>BRACKET TYPE*</b>			
SPD0300	<input type="checkbox"/>	STANDARD SWING AWAY BRACKET	○
SPD0310	<input type="checkbox"/>	ALBER E35 SWIVEL ARM	£360
<b>USER REMOTE*</b>			
SPD0400	<input type="checkbox"/>	USER REMOTE CONTROL RH	○
SPD0410	<input type="checkbox"/>	USER REMOTE CONTROL LH	○
<b>DUAL REMOTE</b>			
SPD0440	<input type="checkbox"/>	INTUITIVE DUAL REMOTE CONTROL RH	£853
SPD0450	<input type="checkbox"/>	INTUITIVE DUAL REMOTE CONTROL LH	£853
SPD0460	<input type="checkbox"/>	INTUITIVE SINGLE DUAL REMOTE CONTROL (ONE HANDLE) (ADJUSTABLE RH TO LH SIDE)	£608
<b>SPEED*</b>			
SPD0500	<input type="checkbox"/>	6 KM/H	◆
<b>BATTERY*</b>			
SPD0600	<input type="checkbox"/>	BATTERY 7,65AH/36V (NON SPILLABLE LEAD ACID BATTERY)	◆
<b>POWER OPTIONS</b>			
SPD0610	<input type="checkbox"/>	ADDITIONAL BATTERY 7,65AH/36V	£271
SPD0630	<input type="checkbox"/>	REMOTE PROTECTOR	£62
<b>SEAT WIDTH*</b>			
SPD1000	<input type="checkbox"/>	305 MM SEAT WIDTH	○
SPD1010	<input type="checkbox"/>	330 MM SEAT WIDTH	○
SPD1020	<input type="checkbox"/>	355 MM SEAT WIDTH	○
SPD1030	<input type="checkbox"/>	380 MM SEAT WIDTH	○
<b>SEAT DEPTH*</b>			
SPD1200	<input type="checkbox"/>	325 TO 375 MM ADJUSTABLE SEAT DEPTH SET AT 350 MM	○
<b>SEAT HEIGHT*</b>			
SPD1300	<input type="checkbox"/>	460 MM SEAT HEIGHT	○
SPD1310	<input type="checkbox"/>	485 MM SEAT HEIGHT	○
<b>SEAT UPHOLSTERY TYPE*</b>			
SPD1400	<input type="checkbox"/>	NYLON WITH LOOP BANDS STRIPS SEAT UPHOLSTERY, WIDTH & DEPTH ADJUSTABLE (SINGLE CROSSBAR ONLY)	◆
<b>BACKREST TYPE*</b>			
SPD2000	<input type="checkbox"/>	FIXED 7° BACKREST 2 HEIGHTS 400/430 MM (DELIVERED 430 MM)	◆
SPD2020	<input type="checkbox"/>	FOLDING BACKREST 430 MM HEIGHT	£32
SPD2040	<input type="checkbox"/>	ANGLE ADJUSTABLE AND FOLD DOWN (-15° TO + 15°) FIXED 7° HEIGHT 425-450 MM / 17"-18" SET AT 425 MM / 17" (2 POSITIONS; 25 MM / 1" STEPS THROUGH UPHOLSTERY)	£97
SPD2050	<input type="checkbox"/>	RECLINING/FOLDING (GAS SPRING) 0/30° BACKREST 3 HEIGHTS 390/410/430 MM (DELIVERED 430 MM)	£64
SPD2060	<input type="checkbox"/>	RECLINING (RATCHET) 0/30° BACKREST 3 HEIGHTS 390/410/430 MM (DELIVERED 430 MM)	£124
SPD2080	<input type="checkbox"/>	RECLINING (GAS SPRING) 0/30° BACKREST 3 HEIGHTS 390/410/430 MM (DELIVERED 430 MM)	£51
<b>BACKREST UPHOLSTERY TYPE*</b>			
SPD2300	<input type="checkbox"/>	LUMBAR SUPPORT (+1 STRAP) BI MATERIAL PADDED COVER	◆
SPD2320	<input type="checkbox"/>	CONTOUR ADJUSTABLE (FULL STRAPS) BI MATERIAL COVER	○

<b>HANGERS RH*</b>		
SPD3000	<input type="checkbox"/>	80° SHORT, UL FROM 200 TO 330MM <span style="float: right;">◆</span>
<b>HANGERS LH*</b>		
SPD3100	<input type="checkbox"/>	80° SHORT, UL FROM 200 TO 330MM <span style="float: right;">◆</span>
<b>FOOTPLATES RH*</b>		
SPD3200	<input type="checkbox"/>	FLIP UP RH FOOTPLATE INCLUDING HEEL STRAP <span style="float: right;">○</span>
SPD3210	<input type="checkbox"/>	FLIP UP ANGLE AND DEPTH ADJUSTABLE RH FOOTPLATE, INCLUDING HEEL STRAP <span style="float: right;">£42</span>
SPD3220	<input type="checkbox"/>	FLIP UP ONE PIECE FOOTPLATE, DEPTH & ANGLE ADJUSTABLE (NC 3310) <span style="float: right;">£76</span>
<b>FOOTPLATES LH*</b>		
SPD3300	<input type="checkbox"/>	FLIP UP LH FOOTPLATE, INCLUDING HEEL STRAP <span style="float: right;">○</span>
SPD3310	<input type="checkbox"/>	FLIP UP ANGLE AND DEPTH ADJUSTABLE LH FOOTPLATE, INCLUDING HEEL STRAP <span style="float: right;">£42</span>
SPD9925	<input type="checkbox"/>	No LH footplate required due to FLIP UP ONE PIECE BEING SELECTED (NC 3210) <span style="float: right;">○</span>
<b>ARMRESTS RH*</b>		
SPD4000	<input type="checkbox"/>	REMOVABLE RH DESK, ADJ IN HEIGHT (235-280 for FSTF 460 & 210-280 for FSTF 485MM) DEPTH (+35 MM FRONTWARD) & WIDTH (+25MM OUTWARDS) <span style="float: right;">○</span>
SPD4010	<input type="checkbox"/>	REMOVABLE RH FULL LENGTH, ADJ IN HEIGHT (210 TO 280 MM) DEPTH (+35 MM BACKWARD) & WIDTH (+25MM OUTWARDS) <span style="float: right;">◆</span>
<b>ARMRESTS LH*</b>		
SPD4100	<input type="checkbox"/>	REMOVABLE LH DESK, ADJ IN HEIGHT (235-280 for FSTF 460 & 210-280 for FSTF 485MM) DEPTH (+35 MM FRONTWARD) & WIDTH (+25MM OUTWARDS) <span style="float: right;">○</span>
SPD4110	<input type="checkbox"/>	REMOVABLE LH FULL LENGTH, ADJ IN HEIGHT (210 TO 280 MM) DEPTH (35 MM BACKWARD) & WIDTH (+25MM OUTWARDS) <span style="float: right;">◆</span>
<b>SEAT CUSHION*</b>		
SPD5000	<input type="checkbox"/>	ESSENTIAL CUSHION 375 MM DEPTH <span style="float: right;">○</span>
<b>SEATING OPTIONS</b>		
SPD5110	<input type="checkbox"/>	HEADREST BLACK PU TR26 <span style="float: right;">£153</span>
SPD5120	<input type="checkbox"/>	HEADREST WITH CHEEK SUPPORT BLACK PU TR26 (INCLUDING BACKBRACE EXCEPT FOR LAGUNA) <span style="float: right;">£177</span>
SPD5140	<input type="checkbox"/>	LATERAL POSITIONING PADS (PAIR) <span style="float: right;">£59</span>
SPD5150	<input type="checkbox"/>	ADJUSTABLE CALF STRAP <span style="float: right;">○</span>
<b>BACK BRACE</b>		
SPD5160	<input type="checkbox"/>	BACK BRACE WIDTH ADJUSTABLE <span style="float: right;">○</span>

<b>COLOURS*</b>			
SPD6000	<input type="checkbox"/>	ANTHRACITE BLACK	○
SPD6015	<input type="checkbox"/>	DARK PURPLE	○
SPD6020	<input type="checkbox"/>	WATER BLUE	○
SWO6529	<input type="checkbox"/>	POP BLUE - 3 week lead time	£164
SWO6259	<input type="checkbox"/>	SLATE GREY - 3 week lead time	£164
SWO0806	<input type="checkbox"/>	Electric Green - 3 week lead time	£164
SWO6033	<input type="checkbox"/>	Red Matt - 3 week lead time	£164
SWO0607	<input type="checkbox"/>	Dark Red - 3 week lead time	£164
SWO0515	<input type="checkbox"/>	Pearl Grey - 3 week lead time	£164
SWO0647	<input type="checkbox"/>	Intense Yellow - 3 week lead time	£164
SWO0648	<input type="checkbox"/>	Vibrant Orange - 3 week lead time	£164
SWO0651	<input type="checkbox"/>	Polar Blue - 3 week lead time	£164
SWO0652	<input type="checkbox"/>	Candy Pink - 3 week lead time	£164
SPD6050	<input type="checkbox"/>	SPECIAL COLOUR (ON DEMAND ONLY)	POA
<b>ACCENT COLOUR*</b>			
SPD6100	<input type="checkbox"/>	BLACK	○
SPD6120	<input type="checkbox"/>	LIME	○
<b>FRONT WHEELS*</b>			
SPD7000	<input type="checkbox"/>	200 X 50 MM (8" X 2") BLACK SOLID TYRE	£47
SPD7010	<input type="checkbox"/>	200 X 50 MM (8" X 2") BLACK PNEUMATIC TYRE	◆
<b>REAR WHEELS*</b>			
SPD7100	<input type="checkbox"/>	355 X 70 MM (14" X 27") BLACK SOLID WITH PUNCTURE PROOF TYRE	£82
SPD7110	<input type="checkbox"/>	355 X 70 MM (14" X 27") BLACK PNEUMATIC TYRE	◆
<b>ANTI TIPPER*</b>			
SPD8000	<input type="checkbox"/>	REAR ANTI-TIPPERS, FIXED, PAIR	◆
SPD8010	<input type="checkbox"/>	REAR ANTI-TIPPERS REMOVABLE WITH CRUTCHES (E25), PAIR	£276
<b>LAP BELT*</b>			
SPD8020	<input type="checkbox"/>	POSTURE BELT	◆
<b>TRANSPORTATION BRACKET*</b>			
SPD8030	<input type="checkbox"/>	TRANSPORTATION BRACKET	◆
<b>OTHER OPTIONS</b>			
SPD8100	<input type="checkbox"/>	CANE HOLDER	£20
SPD8200	<input type="checkbox"/>	ACRYLIC TRAY FLIP UP	£308
<b>TOTAL PRICE</b>			

**Non Conforming Options:**

SPD2300	LUMBAR SUPPORT (+1 STRAP) BI MATERIAL PADDED COVER
SPD3210	FLIP UP ANGLE AND DEPTH ADJUSTABLE RH FOOTPLATE, INCLUDING HEEL STRAP
SPD3310	FLIP UP ANGLE AND DEPTH ADJUSTABLE LH FOOTPLATE, INCLUDING HEEL STRAP



For a more complex customer requirement please refer to our Invacare Unique special order form. email [specialsuk@invacare.com](mailto:specialsuk@invacare.com) or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: \_\_\_\_\_ Email: \_\_\_\_\_

Phone: \_\_\_\_\_ Fax: \_\_\_\_\_

Description of special request (please include relevant option code or part number):