



Yes, you can.®

Tel: 01656 776222  
Email: ordersuk@invacare.com  
www.invacare.co.uk

|                           |                     |
|---------------------------|---------------------|
| <b>Account Name :</b>     | <b>Account No :</b> |
| <b>Address :</b>          | <b>Order No :</b>   |
| <b>Tel No:</b>            | <b>Order Qty :</b>  |
| <b>Delivery address :</b> | <b>Order Date :</b> |




**Matrx Mini back**  
Retail order form  
For paediatric applications

**Standard Features:**

Lightweight, contoured and durable backrest  
Adjustable mounting pins allow it to fit chairs with the same or 25mm wider Mini set hardware for easy installation and quick release with 100mm of depth adjustment, recline up to 30° and rotate  
Hardware fits 25, 22 & 19mm diameter back canes with included shims  
Choice of additional positioning options  
Water resistant and breathable covers with Kidabra print



| Code                   | X   | Description       | Price   |
|------------------------|---|-------------------|---------|
| <b>Matrx Mini Back</b> |   |                   |         |
| MINI1007               | Mini Back with Kidabra print inc. Mini Set hardware | 26Wx18H cm 10x07" | £519.00 |
| MINI1010               | Mini Back with Kidabra print inc. Mini Set hardware | 26Wx26H cm 10x10" | £519.00 |
| MINI1012               | Mini Back with Kidabra print inc. Mini Set hardware | 26Wx31H cm 10x12" | £519.00 |
| MINI1107               | Mini Back with Kidabra print inc. Mini Set hardware | 28Wx18H cm 11x07" | £519.00 |
| MINI1110               | Mini Back with Kidabra print inc. Mini Set hardware | 28Wx26H cm 11x10" | £519.00 |
| MINI1112               | Mini Back with Kidabra print inc. Mini Set hardware | 28Wx31H cm 11x12" | £519.00 |
| MINI1207               | Mini Back with Kidabra print inc. Mini Set hardware | 31Wx18H cm 12x07" | £519.00 |
| MINI1210               | Mini Back with Kidabra print inc. Mini Set hardware | 31Wx26H cm 12x10" | £519.00 |
| MINI1212               | Mini Back with Kidabra print inc. Mini Set hardware | 31Wx31H cm 12x12" | £519.00 |
| MINI1214               | Mini Back with Kidabra print inc. Mini Set hardware | 31Wx36H cm 12x14" | £519.00 |
| MINI1307               | Mini Back with Kidabra print inc. Mini Set hardware | 33Wx18H cm 13x07" | £519.00 |
| MINI1310               | Mini Back with Kidabra print inc. Mini Set hardware | 33Wx26H cm 13x10" | £519.00 |
| MINI1312               | Mini Back with Kidabra print inc. Mini Set hardware | 33Wx31H cm 13x12" | £519.00 |
| MINI1314               | Mini Back with Kidabra print inc. Mini Set hardware | 33Wx36H cm 13x14" | £519.00 |
| MINI1407               | Mini Back with Kidabra print inc. Mini Set hardware | 36Wx18H cm 14x07" | £519.00 |
| MINI1410               | Mini Back with Kidabra print inc. Mini Set hardware | 36Wx26H cm 14x10" | £519.00 |
| MINI1412               | Mini Back with Kidabra print inc. Mini Set hardware | 36Wx31H cm 14x12" | £519.00 |
| MINI1414               | Mini Back with Kidabra print inc. Mini Set hardware | 36Wx36H cm 14x14" | £519.00 |

| <b>Additional options</b>  |                             |   |                             |
|--|-----------------------------|---|-----------------------------|
|   | Swing away lateral          |  | Chest harness interface kit |
|  | Mini-Flex support band      |   |                             |
| 1627466  | Swing-away lateral (each)   | Hardware only   | £140.00                     |
| EPAD-XS  | Extra small lateral pad     | Hardware not included   | £86.00                      |
| 1627465  | Chest harness interface kit |   | £59.00                      |
| 1631855  | Mini-Flex support band      | 4.5in high  | £158.00                     |
| 1631856  | Mini-Flex support band      | 6in high  | £158.00                     |

**Mini Elan headrest**



|         |                                  |                |         |
|---------|----------------------------------|----------------|---------|
| 1603251 | Mini Elan headrest (28cm)        | Hardware only  | £184.00 |
| 1603255 | Standard Pad - Small (16Wx9H cm) | Standard cover | £124.00 |
| 1603294 | Standard Pad - Small (16Wx9H cm) | Wipedown cover | £124.00 |

**Replacement options**

|         |                                    |            |        |         |
|---------|------------------------------------|------------|--------|---------|
| 1626831 | Mini Back cover with Kidabra print | 26Wx18H cm | 10x07" | £116.00 |
| 1626832 | Mini Back cover with Kidabra print | 26Wx26H cm | 10x10" | £116.00 |
| 1626833 | Mini Back cover with Kidabra print | 26Wx31H cm | 10x12" | £116.00 |
| 1626834 | Mini Back cover with Kidabra print | 28Wx18H cm | 11x07" | £116.00 |
| 1626835 | Mini Back cover with Kidabra print | 28Wx26H cm | 11x10" | £116.00 |
| 1626836 | Mini Back cover with Kidabra print | 28Wx31H cm | 11x12" | £116.00 |
| 1626837 | Mini Back cover with Kidabra print | 31Wx18H cm | 12x07" | £116.00 |
| 1626838 | Mini Back cover with Kidabra print | 31Wx26H cm | 12x10" | £116.00 |
| 1626839 | Mini Back cover with Kidabra print | 31Wx31H cm | 12x12" | £116.00 |
| 1626840 | Mini Back cover with Kidabra print | 31Wx36H cm | 21x14" | £116.00 |
| 1626841 | Mini Back cover with Kidabra print | 33Wx18H cm | 13x07" | £116.00 |
| 1626842 | Mini Back cover with Kidabra print | 33Wx26H cm | 13x10" | £116.00 |
| 1626843 | Mini Back cover with Kidabra print | 33Wx31H cm | 13x12" | £116.00 |
| 1626844 | Mini Back cover with Kidabra print | 33Wx36H cm | 13x14" | £116.00 |
| 1626845 | Mini Back cover with Kidabra print | 36Wx18H cm | 14x07" | £116.00 |
| 1626846 | Mini Back cover with Kidabra print | 36Wx26H cm | 14x10" | £116.00 |
| 1626847 | Mini Back cover with Kidabra print | 36Wx31H cm | 14x12" | £116.00 |
| 1626848 | Mini Back cover with Kidabra print | 36Wx36H cm | 14x14" | £116.00 |
| 1626822 | Matrx Mini-Set Hardware            |            |        | £122.00 |

**TOTAL PRICE:**

**Special Instructions - Please supply the following information:**



**Name:**

**Email:**

**Phone:**

**Description of special request (please include relevant option code):**