



Yes, you can:

Tel: 01656 776222
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

Manual Retail Prescription Form

DECEMBER 2025
Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

Rea Azalea Transit Stock

Standard Features

- Max user weight: 135kg (21.5 stone)
- Crash tested ISO 7176-Part 19
- Transit Propulsion 16" Black solid tyre, attendant hub brakes with +50mm pre-width adjustment
- Small, medium or large chassis
- Seat height 450mm/18"
- Depth adjustable seat 430mm-500mm/17"-20"
- Flex 3 Backrest Plate
- Multi-stretch Black PU coated Laguna backrest cushion
- Multi-stretch Black PU coated Flo-shape seat cushion
- Connected push handles
- Assistant operated gas piston seat and backrest recline with tilt lock
- Removable, angle adjustable legrests with calf pads and angle adjustable footplates
- Height adjustable wide pad armrest (including autolock)
- Sideways adjustable headrest
- Trunk supports, multi-adjustable extended height
- 200 X 45mm Black medium hard front castors
- User manoeuvred tyre brakes
- Pelvic belt, seat belt attachment bracket
- Front transportation brackets
- Anti-tippers
- Silver grey frame



Please make your selection with a cross in the box ☐ Free of charge ☒ standard feature POA special request

Customer weight = KG

Code	Description	
Rea Azalea Transit with +5cm width adjustment		
1587805	<input type="checkbox"/> W 39cm x D 43-50cm	£2,393
1587806	<input type="checkbox"/> W 44cm x D 43-50cm	£2,393
1587807	<input type="checkbox"/> W 49cm x D 43-50cm	£2,393

TOTAL PRICE

Non Conforming Options:



For a more complex customer requirement
please refer to our Invacare Unique special order form.
email specialsuk@invacare.com

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):