

**AQUATEC®**

# Toilet seat raisers

**a range of adaptable toileting solutions**



**Providing dependable toileting support**

Adjustable, reliable, comfortable



**Yes, you can.®**



## Introducing the Aquatec range of raised toilet seats

An adaptable range offering supportive toileting solutions.

Providing a reliable and supportive toileting solution for individuals is very important to aid independence.

The range offers a variety of possible seat heights to accommodate individual sizes and to ensure optimal comfort.

The dependable Aquatec range offers a secure and stable environment with fixations that support the placement of the raised toilet seats in situ.

Featuring models that offer supportive armrests for stability and standing transfers, including the swing-away options for ease of functionality.



### Aquatec 900

- > Individual adjustment at any time, no tools required:
  - three different seat heights (6 cm, 10 cm, 15 cm),
  - three different angles, lower at the front.
- > Swing-away, ergonomic armrests for additional support when sitting down and standing up.
- > Available with or without armrests.
- > Seat can be folded up in all setting positions - ease of cleaning.



### Aquatec 90000

- > Available in three different seat heights: 2 cm, 6 cm, 10 cm.
- > Can be ordered with or without armrests.
- > Swing-away, ergonomic armrests for additional support when sitting down and standing up.
- > Innovative clip system: The lid and toilet seat can be removed from the base; no tools are required.



# Technical data

For more comprehensive pre-sales information about this product range, including the product's user manuals, please see your local Invacare website.



Seat width ▼



Seat depth ▼



Seat height ▼



Total width ▼

	Seat width ▼	Seat depth ▼	Seat height ▼	Total width ▼
<b>Aquatec 900</b>	360 mm	425 mm	60 / 100 / 150 mm	590 mm (with armrests) 430 mm (without armrests)
<b>Aquatec 90000</b>	365 mm	390 mm	20 / 60 / 100 mm	540 mm (with armrests) 380 mm (without armrests)



Total depth ▼



Total height ▼



Width between armrests ▼



Width of opening ▼

	Total depth ▼	Total height ▼	Width between armrests ▼	Width of opening ▼
<b>Aquatec 900</b>	490 mm	380 mm (with armrests) 180 mm (without armrests) 200 mm (seat to armrests)	510 mm	225 mm
<b>Aquatec 90000</b>	475 mm	230 / 270 / 310 mm (with armrests) 30 / 70 / 110 mm (without armrests) 200 mm (seat to armrests)	450 mm	205 mm



Depth of opening ▼



Weight ▼



Load capacity ▼

	Depth of opening ▼	Weight ▼	Load capacity ▼
<b>Aquatec 900</b>	285 mm	5.0 kg (with armrests) 2.7 kg (without armrests)	120 kg
<b>Aquatec 90000</b>	275 mm	4.4 kg (SH:100 mm) 4 kg (SH:60 mm) 3.6 kg (SH:20 mm) (with armrests)	150 kg



Frame colour

Invacare Limited  
 Pencoed Technology Park - Pencoed  
 Bridgend CF35 5AQ - United Kingdom  
 Tel: +44 1 656 776 222  
 Fax: +44 1656 776 220  
 E-mail: [uk@invacare.com](mailto:uk@invacare.com)  
 Sales Order E-mail:  
[ordersuk@invacare.com](mailto:ordersuk@invacare.com)  
[www.invacare.co.uk](http://www.invacare.co.uk)

Invacare Ireland Ltd  
 Unit 5 - Seatown Business Campus  
 Seatown Road - SWORDS County  
 Dublin - Ireland Tel. +353 1 8107084  
 Fax +353 1 8107085  
 Email: [Ireland@invacare.com](mailto:Ireland@invacare.com)  
[www.invacare.ie](http://www.invacare.ie)

