



Yes, you can:

Tel: 01656 776222 Fax: 01656 776220
 Email: ordersuk@invacare.com
 Online Spares available at: www.invacare.co.uk

APRIL 2024
 Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

Action 3NG Transit One Solution

Standard Features

- Max user weight 125kg (19.7 stone)
- Crash tested ISO 7176-Part 19
- Attendant Propulsion
- Choice of 6 seat widths 380mm to 505mm (15" to 20")
- Choice of 2 seat depths 400mm or 450mm (16" or 18)
- Choice of 3 seat heights 460mm to 510mm (18" to 20")
- Black nylon upholstery
- Fixed 7° 430mm backrest
- Lumbar support padded backrest with shaped cover
- 80° Swing in / swing out legrest hangers
- Flip up footplates
- Swing away and detachable desk armrests
- Choice of 8 colours (with 7 special colours)
- Push to lock brakes
- Transportation kit and lap belt (adjustable both sides)
- Adjustable heel straps
- User tool kit

Standard Base price: £641



Please make your selection with a cross in the box ○ Free of charge ◆ standard feature POA special request

Customer weight = KG

Code	Description	
BASE PRODUCT*		
CEF0004	ACTION 3NG TRANSIT (NC 1074)	○
CEFSPEC	SPECIAL	POA
PROPULSION*		
CEF0110	TRANSIT PROPULSION WITH RH TIPPING AID	£641

SEAT WIDTH*			
CEF44444		13" / 330 MM SEAT WIDTH - SELECT SEAT WIDTH OPTION CEE0180 AND IN-RIGGED ARMREST OPTION CEE0952 (lapbelt adjustment 420mm - 1520mm)	£41
CEF0180		15" / 380 MM SEAT WIDTH (NC 0485, 0486, 0487, 0488, 0852)	○
CEF0190		16" / 405 MM SEAT WIDTH (NC 0482, 0486, 0487, 0488)	○
CEF0200		17" / 430 MM SEAT WIDTH (NC 0482, 0483, 0487, 0488)	◆
CEF0210		18" / 455 MM SEAT WIDTH (NC 0482, 0483, 0484, 0488)	○
CEF0220		19" / 480 MM SEAT WIDTH (NC 0482, 0483, 0484, 0485)	○
CEF0230		20" / 505 MM SEAT WIDTH (NC 0482, 0483, 0484, 0485, 0486)	○
SEAT DEPTH*			
CEF0240		16" / 400 MM SEAT DEPTH (NC 0879)	○
CEF0250		18" / 450 MM SEAT DEPTH (NC 0878)	◆
CEF0252		16" TO 18" / 400 TO 450 MM ADJUSTABLE SEAT DEPTH (NC 0330, 0350)	£21
SEAT HEIGHT*			
CEF0260		20" / 510 MM SEAT HEIGHT	○
CEF0270		19" / 485 MM SEAT HEIGHT	○
CEF0280		18" / 460 MM SEAT HEIGHT	○
SEAT UPHOLSTERY*			
CEF0330		NYLON W/O HOOK AND LOOP STRIPS	◆
CEF0340		NYLON WITH HOOK AND LOOP STRIPS	○
CEF0345		TENSION ADJUSTABLE SEAT UPHOLSTERY & DEPTH ADJUSTABLE -50MM	○
CEF0350		PADDED NYLON	£30
BACKREST - Anti-tippers should be selected with reclining backrests 0435, 0440, 0445 and 0450*			
CEF0380		FIXED 0° BACKREST 2 HEIGHTS 16" / 400 MM or 17" / 430 MM	○
CEF0390		FIXED 0° BACKREST 20" / 510 MM HEIGHT	○
CEF0410		FIXED 7° BACKREST 2 HEIGHTS 16" / 400 or 17" / 430 MM FIX PUSH HANDLES	◆
CEF0420		FIXED 7° BACKREST 20" / 510 MM HEIGHT FIX PUSH HANDLES	○
CEF0432		ANGLE ADJUSTABLE AND FOLD DOWN (15° to + 15°) 17"	£86
CEF0450		RECLINING (GAS SPRING) 0/30° 3 HEIGHTS 18" 19" 20"	£233

MATRX BACKREST (LOSS OF APPROX 20 MM SEAT DEPTH WITH RECLINING GAS SPRING BACKREST)*			
CEF0451		MATRX ELITE E2 BACK WITH HARDWARE - 160 MM BACK HEIGHT	£494
CEF0452		MATRX ELITE E2 BACK WITH HARDWARE - 210 MM BACK HEIGHT	£494
CEF0453		MATRX ELITE E2 BACK WITH HARDWARE - 260 MM BACK HEIGHT	£494
CEF0454		MATRX ELITE E2 BACK WITH HARDWARE - 310 MM BACK HEIGHT	£494
CEF0456		MATRX ELITE E2 BACK WITH HARDWARE - 360 MM BACK HEIGHT	£494
CEF0461		MATRX ELITE E2 BACK WITH HARDWARE - 410 MM BACK HEIGHT	£494
CEF0462		MATRX ELITE E2 BACK WITH HARDWARE - 460 MM BACK HEIGHT	£494
CEF0457		MATRX ELITE E2 BACK WITH HARDWARE - 510 MM BACK HEIGHT	£494
CEF0458		MATRX ELITE DEEP E2 BACK WITH HARDWARE - 160 MM BACK HEIGHT	£546
CEF0459		MATRX ELITE DEEP E2 BACK WITH HARDWARE - 210 MM BACK HEIGHT	£546
CEF0466		MATRX ELITE DEEP E2 BACK WITH HARDWARE - 260 MM BACK HEIGHT	£546
CEF0467		MATRX ELITE DEEP E2 BACK WITH HARDWARE - 310 MM BACK HEIGHT	£546
CEF0468		MATRX ELITE DEEP E2 BACK WITH HARDWARE - 360 MM BACK HEIGHT	£546
CEF0463		MATRX ELITE DEEP E2 BACK WITH HARDWARE - 410 MM BACK HEIGHT	£546
CEF0464		MATRX ELITE DEEP E2 BACK WITH HARDWARE - 460 MM BACK HEIGHT	£546
CEF0469		MATRX ELITE DEEP E2 BACK WITH HARDWARE - 510 MM BACK HEIGHT	£546
MATRX BACKREST WIDTH CHOICE*			
CEF0482		360mm BACKREST WIDTH, TO FIT CHEST WIDTH 265mm	○
CEF0483		380mm BACKREST WIDTH, TO FIT CHEST WIDTH 290mm	○
CEF0484		410mm BACKREST WIDTH, TO FIT CHEST WIDTH 315mm	○
CEF0485		430mm BACKREST WIDTH, TO FIT CHEST WIDTH 340mm	○
CEF0486		460mm BACKREST WIDTH, TO FIT CHEST WIDTH 370mm	○
CEF0487		480mm BACKREST WIDTH, TO FIT CHEST WIDTH 390mm	○
CEF0488		510mm BACKREST WIDTH, TO FIT CHEST WIDTH 420mm	○
LEGREST HANGER RH*			
CEF0510		80° RH SWING IN/SWING OUT (NC 1111, 66666)	◆
CEF0520		ELEVATING LEGREST RH SWING IN/SWING OUT (NC 0611, 0930, 0935, 0936, 66666)	£116
CEF0545		AMPUTEE LEGREST RH (NC 0600, 0610, 0611, 0930, 0936, 1111)	£133
LEGREST HANGER LH*			
CEF0560		80° LH SWING IN/SWING OUT (NC 2222)	◆
CEF0570		ELEVATING LEGREST LH SWING IN/SWING OUT (NC 0611, 0930, 0935, 0936, 55555)	£116
CEF0595		AMPUTEE LEGREST LH (NC 0611, 0620, 0630, 0930, 0936, 2222)	£133

FOOTREST RH*			
CEF0600		FLIP UP RH FOOTPLATE (NC 0937, 0938)	◆
CEF0610		FLIP UP DEPTH AND ANGLE AJUSTABLE RH FOOTPLATE - NB: IF SELECTED WITH A FIXED FOOTREST - FOOTPLATE WILL BE ANGLE AND DEPTH ADJ AND COMPATIBLE WITH ANKLE HUGGERS. IF SELECTED WITH AN ELEVATING LEGREST - FOOTPLATE WILL BE ANGLE, DEPTH AND WIDTH ADJ - NOT COMPATIBLE WITH ANKLE HUGGERS.	£43
SWO1111		FLIP UP DEPTH AND ANGLE ADJUSTABLE RH FOOTPLATE - NB: COMPATIBLE WITH ANKLE HUGGERS (FOR ELEVATING LEGRESTS ONLY)	£43
CEF0611		FLIP UP ONE PIECE (NC 0620, 0630, 0935, 0936, 0937, 0938)	£78
CEF66666		NO FOOTREST	○
FOOTREST LH*			
CEF0620		FLIP UP LH FOOTPLATE (NC 0937, 0938)	◆
CEF0630		FLIP UP DEPTH AND ANGLE AJUSTABLE LH FOOTPLATE - NB: IF SELECTED WITH A FIXED FOOTREST - FOOTPLATE WILL BE ANGLE AND DEPTH ADJ AND COMPATIBLE WITH ANKLE HUGGERS. IF SELECTED WITH AN ELEVATING LEGREST - FOOTPLATE WILL BE ANGLE, DEPTH AND WIDTH ADJ - NOT COMPATIBLE WITH ANKLE HUGGERS.	£43
SWO2222		FLIP UP DEPTH AND ANGLE ADJUSTABLE RH FOOTPLATE - NB: COMPATIBLE WITH ANKLE HUGGERS (FOR ELEVATING LEGRESTS ONLY)	£43
CEF55555		NO FOOTREST	○
ARMREST RH*			
CEF0640		REMOVABLE RH (N°1) DESK ADJ. 200 TO 280 MM HEIGHT AUTOLOCK SAFETY DEVICE INCLUDED (NC 0760, 0770, 0780, 0790, 0795, 0796, 1480, 1483)	○
CEF0650		REMOVABLE RH (N°1) FULL LENGTH ADJ. 200 TO 280 MM HEIGHT AUTOLOCK SAFETY DEVICE INCLUDED (NC 0760, 0770, 0780, 0790, 0795, 0796)	○
CEF0660		SWING AWAY DETACHABLE RH (N°2) DESK 240 MM HEIGHT (NC 0740, 0750, 0780, 0790, 0795, 0796, 0955, 1480, 1483)	◆
CEF0670		SWING AWAY DETACHABLE RH (N°2) FULL LENGTH 240 MM HEIGHT (NC 0740, 0750, 0780, 0790, 0795, 0796, 0955)	○
CEF0680		SWING AWAY DETACHABLE RH (N°3) DESK ADJ. IN HEIGHT (170-295 MM) (NC 0740, 0750, 0760, 0770, 0795, 0796, 0955, 1480, 1483)	£56
CEF0690		SWING AWAY DETACHABLE RH (N°3) FULL LENGTH ADJ. IN HEIGHT (170-295 MM) AND DEPTH (NC 0740, 0750, 0760, 0770, 0795, 0796, 0955, 1483)	£56
CEF0695		REMOVABLE RH (N°4) DESK ADJ. IN HEIGHT (210-345 MM) + SKIRTGUARD HEIGHT ADJ. AUTOLOCK SAFETY DEVICE INCLUDED (NC 0740, 0750, 0760, 0770, 0780, 0790, 0955, 1480, 1483)	£60
CEF0696		REMOVABLE RH (N°4) FULL LENGTH ADJ. IN HEIGHT (210-345 MM) AND DEPTH + SKIRTGUARD HEIGHT ADJ. AUTOLOCK SAFETY DEVICE INCLUDED (NC 0740, 0750, 0760, 0770, 0780, 0790, 0955, 1483)	£60
CEF0732		REMOVABLE RH (N°4) ARMPAD HEMI BLACK MULTI STRETCH PU COATED FABRIC MULTI POSITION AUTOLOCK SAFETY DEVICE INCLUDED (NC 1480, 1483)	£59

ARMREST LH*			
CEF0740		REMOVABLE LH (N°1) DESK ADJ. 200 TO 280 MM HEIGHT AUTOLOCK SAFETY DEVICE INCLUDED (NC 1480, 1484)	○
CEF0750		REMOVABLE LH (N°1) FULL LENGTH ADJ. 200 TO 280 MM HEIGHT AUTOLOCK SAFETY DEVICE INCLUDED	○
CEF0760		SWING AWAY DETACHABLE LH (N°2) DESK 240 MM HEIGHT (NC 0955, 1480, 1484)	◆
CEF0770		SWING AWAY DETACHABLE LH (N°2) FULL LENGTH 240 MM HEIGHT (NC 0955)	○
CEF0780		SWING AWAY DETACHABLE LH (N°3) DESK ADJ. IN HEIGHT (170-295 MM) (NC 0955, 1480, 1484)	£56
CEF0790		SWING AWAY DETACHABLE LH (N°3) FULL LENGTH ADJ. IN HEIGHT (170-295 MM) AND DEPTH (NC 0955, 1484)	£56
CEF0795		REMOVABLE LH (N°4) DESK ADJ. IN HEIGHT (210-345 MM) + SKIRTGUARD HEIGHT ADJ. AUTOLOCK SAFETY DEVICE INCLUDED (NC 0955, 1480, 1484)	£60
CEF0796		REMOVABLE LH (N°4) FULL LENGTH ADJ. IN HEIGHT (210-345 MM) AND DEPTH + SKIRTGUARD HEIGHT ADJ. AUTOLOCK SAFETY DEVICE INCLUDED (NC 0955, 1484)	£60
CEF0832		REMOVABLE LH N°4 ARMPAD HEMI BLACK MULTI STRETCH PU COATED FABRIC MULTI POSITION ARMPAD HEMI MULTI POSITION (NC 1480, 1484)	£59
SEAT CUSHION			
CEF0849		BASIC 50 MM FOAM CUSHION - NYLON AND ANTI-SLIP COVER	£24
CEF0869		MATRX FLO-TECH LITE CUSHION - Low Risk	£121
CEF0871		MATRX FLO-TECH CONTOUR VISCO CUSHION - High Risk	£169
CEF0873		MATRX FLO-TECH IMAGE CUSHION - High Risk	£288
CEF0874		MATRX FLO-TECH SOLUTION XTRA CUSHION - Very High Risk	£398
CEF0876		MATRX LIBRA CUSHION - Very High Risk	£410
SEAT CUSHION DEPTH			
CEF0878		FOR 400 MM SEAT DEPTH	○
CEF0879		FOR 450 MM SEAT DEPTH	○
OTHER SEATING OPTIONS			
CEF0852		SEAT RIGIDIZER	£39
CEF0905		NECKREST BLACK MULTI STRETCH PU COATED FABRIC (NC 0906, 0907, 0908)	£131
CEF0906		HEADREST BLACK MULTI STRETCH PU COATED FABRIC (NC 0907, 0908)	£159
CEF0907		HEADREST WITH CHEEK SUPPORT BLACK MULTI STRETCH PU COATED FABRIC (NC 0908)	£183
CEF0908		MATRX LOXX SOFT PAD HEADREST (SOFT PAD 250MM) WITH FIXED HARDWARE (ONLY FOR MATRX BACK)	£269
SWO3333		MULT ADJUSTABLE LATERALS	£284
CEF0930		CALF STRAP	◆
CEF0935		LONGER FOOTREST TUBE	£30
CEF0936		SHORTER FOOTREST SUPPORT	£42
CEF0937		BODYPOINT ANKLE HUGGERS (PAIR) SIZE M - UNMOUNTED	£114
CEF0938		BODYPOINT ANKLE HUGGERS (PAIR) SIZE L - UNMOUNTED	£114
CEF0955		W/O AUTOLOCK SAFETY DEVICE FOR ARMREST N°1 ONLY	○

COLOUR*			
CEF0961	<input type="checkbox"/>	SLATE GREY (RAL7024 FINE TEXTURE)	○
CEF0964	<input type="checkbox"/>	PEARL WHITE	○
CEF0965	<input type="checkbox"/>	ELECTRIC GREEN	○
CEF0966	<input type="checkbox"/>	MATT BLACK	○
CEF0969	<input type="checkbox"/>	RED MATT	○
CEF0972	<input type="checkbox"/>	WATER BLUE	○
CEA0974	<input type="checkbox"/>	BLUE	○
CEA0976	<input type="checkbox"/>	DARK PURPLE	○
CEF0980	<input type="checkbox"/>	ANTHRACITE BLACK	○
CEF0991	<input type="checkbox"/>	SPECIAL COLOUR ON DEMAND ONLY	£169
SWO4444	<input type="checkbox"/>	DARK RED - SPECIAL	£169
SWO5555	<input type="checkbox"/>	PEARL GREY - SPECIAL	£169
SWO6666	<input type="checkbox"/>	INTENSE YELLOW	£169
SWO7777	<input type="checkbox"/>	VIBRANT ORANGE	£169
SWO8888	<input type="checkbox"/>	POLAR BLUE	£169
SWO9999	<input type="checkbox"/>	CANDY PINK	£169
CASTOR HOUSING*			
CEF0995	<input type="checkbox"/>	STANDARD	◆
CEF0996	<input type="checkbox"/>	LIGHT WEIGHT	£28
FRONT CASTORS*			
CEF1010	<input type="checkbox"/>	8" X 1.1/4" BLACK SOLID TYRE	◆
CEF1030	<input type="checkbox"/>	8" X 2" BLACK SOLID TYRE	£28
REAR WHEEL SIZE*			
CEF1062	<input type="checkbox"/>	12"	◆
REAR WHEELS*			
CEF1073	<input type="checkbox"/>	DRUM BRAKE - TRANSIT 12" (NC 1201, 1270)	○
CEF1074	<input type="checkbox"/>	TRANSIT COMPOSITE (BLACK PLASTIC HUB AND RIM) (NC 1280)	◆
TYRES*			
CEF1082	<input type="checkbox"/>	SOLID BLACK	◆
REAR WHEEL AXLE*			
CEF1170	<input type="checkbox"/>	FIXED AXLE	◆
REAR WHEEL POSITION*			
CEF1200	<input type="checkbox"/>	PASSIVE REAR WHEELS POSITION	◆
CEF1201	<input type="checkbox"/>	DUAL AMPUTEE REAR WHEEL POSITION (NC 1280)	£72
CHASSIS OPTION			
CEF1241	<input type="checkbox"/>	FRONT WHEELS EXTENSION (+ 50 MM OUTWARD)	£22
BRAKES*			
CEF1270	<input type="checkbox"/>	PUSH TO LOCK BRAKES	◆
CEF1280	<input type="checkbox"/>	ATTENDANT HUB BRAKES (INCLUDING PUSH TO LOCK BRAKES)	£231
TRANSPORTATION KIT*			
CEF1345	<input type="checkbox"/>	TRANSPORTATION BRACKETS	◆

SAFETY OPTIONS*			
CEF1360		LAP BELT	(NC 1361, 1362, 1363, 1364) ◆
CEF1361		BODYPOINT 2-PT POSTURE BELT M	(NC 1362, 1363, 1364) £83
CEF1362		BODYPOINT 2-PT POSTURE BELT L	(NC 1363, 1364) £83
CEF1363		BODYPOINT 4-PT POSTURE BELT M	(NC 1364) £98
CEF1364		BODYPOINT 4-PT POSTURE BELT L	£98
CEF1370		REAR ANTI-TIPPERS PAIR	(NC 1375, 1440, 1441) £36
CEF1375		REAR ANTI-TIPPERS PAIR WITH TIPPING AID	(NC 1440, 1441) £72
OTHER OPTIONS			
CEF1430		CANE HOLDER	£23
CEF1440		TIPPING AID RH	£14
CEF1441		TIPPING AID LH	£14
CEF1480		ACRYLIC TRAY	(NC 1483, 1484, 1485) £116
CEF1483		HALF TABLE TRAY RH	(NC 1484) £89
CEF1484		HALF TABLE TRAY LH	£89
CEF1485		COMFORT PADDING FOR HALF TABLE TRAY	£44
SWO10001		ABDUCTED LEGRESTS (PR) 0° TO 30° - INTERNAL REFERENCE DE259618	£107
1595002		ViaGo Push aid. Includes as standard: Drive Unit, Charger, Control unit & Brackets (Note: 10 minute fitting time)	£1,995

TOTAL PRICE	
--------------------	--

Non Conforming Options:

CEF0330	NYLON W/O HOOK AND LOOP STRIPS
CEF0350	PADDED NYLON
CEF0482	360mm BACKREST WIDTH, TO FIT CHEST WIDTH 265mm
CEF0483	380mm BACKREST WIDTH, TO FIT CHEST WIDTH 290mm
CEF0484	410mm BACKREST WIDTH, TO FIT CHEST WIDTH 315mm
CEF0485	430mm BACKREST WIDTH, TO FIT CHEST WIDTH 340mm
CEF0486	460mm BACKREST WIDTH, TO FIT CHEST WIDTH 370mm
CEF0487	480mm BACKREST WIDTH, TO FIT CHEST WIDTH 390mm
CEF0488	510mm BACKREST WIDTH, TO FIT CHEST WIDTH 420mm
CEF0600	FLIP UP RH FOOTPLATE
CEF0610	FLIP UP DEPTH AND ANGLE AJUSTABLE RH FOOTPLATE - NB: IF SELECTED WITH A FIXED FOOTREST - FOOTPLATE WILL BE ANGLE AND DEPTH ADJ AND COMPATIBLE WITH ANKLE HUGGERS. IF SELECTED WITH AN ELEVATING LEGREST - FOOTPLATE WILL BE ANGLE, DEPTH AND WIDTH ADJ - NOT COMPATIBLE WITH ANKLE HUGGERS.
SWO1111	FLIP UP DEPTH AND ANGLE ADJUSTABLE RH FOOTPLATE - NB: COMPATIBLE WITH ANKLE HUGGERS (FOR ELEVATING LEGRESTS ONLY)
CEF0611	FLIP UP ONE PIECE
CEF66666	NO FOOTREST
CEF0620	FLIP UP LH FOOTPLATE
CEF0630	FLIP UP DEPTH AND ANGLE AJUSTABLE LH FOOTPLATE - NB: IF SELECTED WITH A FIXED FOOTREST - FOOTPLATE WILL BE ANGLE AND DEPTH ADJ AND COMPATIBLE WITH ANKLE HUGGERS. IF SELECTED WITH AN ELEVATING LEGREST - FOOTPLATE WILL BE ANGLE, DEPTH AND WIDTH ADJ - NOT COMPATIBLE WITH ANKLE HUGGERS.
SWO2222	FLIP UP DEPTH AND ANGLE ADJUSTABLE RH FOOTPLATE - NB: COMPATIBLE WITH ANKLE HUGGERS (FOR ELEVATING LEGRESTS ONLY)
CEF55555	NO FOOTREST
CEF0740	REMOVABLE LH (N°1) DESK ADJ. 200 TO 280 MM HEIGHT AUTOLOCK SAFETY DEVICE INCLUDED
CEF0750	REMOVABLE LH (N°1) FULL LENGTH ADJ. 200 TO 280 MM HEIGHT AUTOLOCK SAFETY DEVICE INCLUDED
CEF0760	SWING AWAY DETACHABLE LH (N°2) DESK 240 MM HEIGHT
CEF0770	SWING AWAY DETACHABLE LH (N°2) FULL LENGTH 240 MM HEIGHT
CEF0780	SWING AWAY DETACHABLE LH (N°3) DESK ADJ. IN HEIGHT (170-295 MM)
CEF0790	SWING AWAY DETACHABLE LH (N°3) FULL LENGTH ADJ. IN HEIGHT (170-295 MM) AND DEPTH
CEF0795	REMOVABLE LH (N°4) DESK ADJ. IN HEIGHT (210-345 MM) + SKIRTGUARD HEIGHT ADJ. AUTOLOCK SAFETY DEVICE INCLUDED
CEF0796	REMOVABLE LH (N°4) FULL LENGTH ADJ. IN HEIGHT (210-345 MM) AND DEPTH + SKIRTGUARD HEIGHT ADJ. AUTOLOCK SAFETY DEVICE INCLUDED
CEF0878	FOR 400 MM SEAT DEPTH
CEF0879	FOR 450 MM SEAT DEPTH
CEF0852	SEAT RIGIDIZER
CEF0906	HEADREST BLACK MULTI STRETCH PU COATED FABRIC
CEF0907	HEADREST WITH CHEEK SUPPORT BLACK MULTI STRETCH PU COATED FABRIC
CEF0908	MATRIX LOXX SOFT PAD HEADREST (SOFT PAD 250MM) WITH FIXED HARDWARE (ONLY FOR MATRIX BACK)
CEF0930	CALF STRAP
CEF0935	LONGER FOOTREST TUBE
CEF0936	SHORTER FOOTREST SUPPORT
CEF0937	BODYPOINT ANKLE HUGGERS (PAIR) SIZE M - UNMOUNTED
CEF0938	BODYPOINT ANKLE HUGGERS (PAIR) SIZE L - UNMOUNTED
CEF0955	W/O AUTOLOCK SAFETY DEVICE FOR ARMREST N°1 ONLY
CEF1074	TRANSIT COMPOSITE (BLACK PLASTIC HUB AND RIM)
CEF1201	DUAL AMPUTEE REAR WHEEL POSITION

CEF1270	PUSH TO LOCK BRAKES
CEF1280	ATTENDANT HUB BRAKES (INCLUDING PUSH TO LOCK BRAKES)
CEF1361	BODYPOINT 2-PT POSTURE BELT M
CEF1362	BODYPOINT 2-PT POSTURE BELT L
CEF1363	BODYPOINT 4-PT POSTURE BELT M
CEF1364	BODYPOINT 4-PT POSTURE BELT L
CEF1375	REAR ANTI-TIPPERS PAIR WITH TIPPING AID
CEF1440	TIPPING AID RH
CEF1441	TIPPING AID LH
CEF1480	ACRYLIC TRAY
CEF1483	HALF TABLE TRAY RH
CEF1484	HALF TABLE TRAY LH
CEF1485	COMFORT PADDING FOR HALF TABLE TRAY



For a more complex customer requirement please refer to our Invacare Unique special order form. email specialsuk@invacare.com or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: _____ Email: _____

Phone: _____ Fax: _____

Description of special request (please include relevant option code or part number):