

Leo®

up to
38km
range


KG
up to
136 kg
user weight

8
km/h
speed



Main features

- Compact and easy to control
- 360° swivel seat, height and depth adjustable
- Can be dismantled into 6 parts
- Ergonomic and adjustable steering
- Road equipment with headlights and indicators
- Automatic brake light activation



Onyx Blue Ruby Red

Quality
★★★★★

Comfort
↑↑↑

Safety
✓

Orion® Metro

up to
43km
range


KG
up to
136 kg
user weight

12
km/h
speed



Main features

- 360° swivel seat, adjustable in height and depth, with headrest
- 4 suspensions
- Ergonomic and adjustable steering
- Automatic brake light activation
- Automatic speed reduction on bends
- Optional connectivity



Garnet Red Alium Blue Marble White Diamond Silver Jade Green

Quality
★★★★★

Comfort
↑↑↑

Safety
✓

Orion® Pro

up to
54km
range


KG
up to
160 kg
user weight

12
km/h
speed



Main features

- 360° swivel seat, height and depth adjustable
- 4 suspensions
- On-board computer with intuitive LCD screen: display of speed, time, temperature and odometer
- Ergonomic and adjustable steering
- Automatic brake light activation
- Automatic speed reduction when cornering
- Optional connectivity



Garnet Red Alium Blue Marble White Diamond Silver Jade Green

Quality
★★★★★

Comfort
↑↑↑

Safety
✓

Comet® Pro

up to
55km
range


KG
up to
160 kg
user weight

12
km/h
speed



Main features

- 360° swivel seat, height and depth adjustable
- 2 suspensions at the front and 1 central suspension at the rear
- On-board computer with intuitive LCD screen: display of speed, time, temperature and odometer
- Ergonomic and adjustable steering
- Automatic brake light activation
- Automatic speed reduction when cornering
- Optional connectivity



Garnet Red Alium Blue Marble White Diamond Silver Jade Green

Quality
★★★★★

Comfort
↑↑↑


Safety
✓

Comet® Ultra

up to
49km
range


KG
up to
220 kg
user weight

10
km/h
speed



Main features

- Maximum user weight up to 220 kg
- 26" Ultra seat, 360° swivel, height and depth adjustable
- 4 suspensions
- On-board computer with intuitive LCD screen: display of speed, time, temperature and odometer
- Automatic brake light activation
- Reinforced armrests and seat support
- Automatic speed reduction when cornering
- Optional connectivity



Garnet Red Alium Blue Marble White Diamond Silver Jade Green

Quality
★★★★★

Comfort
↑↑↑

Safety
✓

Selection Guide

My priorities	Leo®	Orion® Metro	Orion® Pro	Comet® Pro	Comet® Ultra
Home or shopping centre	●	●	●	●	-
Urban walks	●	●	●	●	●
Countryside	-	●	●	●	●
Transport in the boot of a car?	●	-	-	-	-
My desired travel distance					
1 - 15 km	●	●	●	●	●
15-25 km	●	●	●	●	●
25-35 km	●	●	●	●	●
35-50 km	●	●	●	●	●
My weight					
Less than 115kg	●	●	●	●	●
Less than 136kg	●	●	●	●	●
Less than 160kg	-	-	●	●	●
Less than 220kg	-	-	-	-	●
Desired speed					
Up to 6km/h	-	●	●	●	●
Up to 8km/h	●	●	●	●	●
Up to 10km/h	-	●	●	●	●
Up to 12km/h	-	●	●	●	-

Please note that this selection guide is for reference only. Ask your retailer for advice.

Contact us
for more information
and book a test drive today

Retailer information



Discover more about
our scooter range.
www.invacare.co.uk



Yes, you can.

© 2023 Invacare International GmbH. All rights reserved. All information quoted is believed to be correct at time of print. Invacare reserves the right to alter product specifications without prior consultation. Scooter - UK - 03/2023

Other options



○ Location of accessories

Features and options	Leo [®]	Orion [®] Metro	Orion [®] Pro	Comet [®] Pro	Comet [®] Ultra
Range	38 km	43 km	54 km	55 km	49 km
Speed	8km/h	12km/h	12km/h	12km/h	10km/h
Lighting package	●	●	●	●	●
Automatic stop lights	●	●	●	●	●
Flashing lights	●	●	●	●	●
Load indicator	●	●	●	●	●
Automatic speed reduction for cornering		●	●	●	●
Deluxe seat		●	●	●	●
Premium seat				●	
Ultra seating					●
Easily transports in a car	●				
Headrest	●	●	●	●	●
Suspension		●	●	●	●
Optional connectivity		○	○	○	○

Accessories

External charging kit	○				
Front basket ①	●	●	●	●	●
Back basket ②	○	○	○	○	○
Backpack ③		○	○	○	○
Left/right mirror ④		○	○	○	○
Delux rear view mirror ⑤	○	○	○	○	○
Walking frame holder ⑥	○	○	○	○	○
Rollator holder ⑦		○	○	○	○
Crutch holder ⑧	○	○	○	○	○
Oxygen bottle holder ⑨	○	○	○	○	○
Mobile phone holder ⑩		○	○	○	○
Drink holder ⑪		○	○	○	○
Seat belt ⑫	○	○	○	○	○
Front bumper ⑬		○	○	○	○
Rear bumper ⑭		○	○	○	○
Windscreen ⑮	○	○	○	○	○
Seat suspension ⑯	○			○	○
Protective coverage ⑰	○	○	○	○	○

● Standard ○ Optional

Feel the Freedom
with Invacare Scooters



Yes, you can.