



via GO®

NEW

The compact push aid for wheelchairs



Pushing
the easy way.

via GO – push smarter, not harder.



The powerful, simplified push aid for wheelchairs.

There are times when even a poor pavement surface or a slight incline can become an insurmountable obstacle for everyday tasks. Via Go - the powerful push aid from Alber - provides the right amount of assistance for a family member or carer, helping to handle these situations with confidence. The simplified design and on demand push assistance will allow enjoyment of excursions with peace of mind.

The high quality, compact push aid fits most common wheelchair models, your Invacare (Alber) approved retailer will check whether your wheelchair is compatible and provide expert support for fitting and operation.

via GO – pushing the limits.

www.invacare.co.uk

COMPACT AND CONVENIENT



Experience more

With its impressive versatility, via GO will become an everyday companion – for day-to-day tasks and excursions alike. Thanks to its compact dimensions and lightweight design, this push aid is very easy to handle and can be stowed away quickly to save space.

STRONG AND ENDURING



No effort

Almost every outdoor situation has inclines and ramps. The via GO gives you the right amount of assistance just when you need it, helping you to get your loved ones on the move with no problem at all – for up to 10 miles on a single charge.

INTUITIVE AND FLEXIBLE



Clever thinking

The via GO is so straightforward you can use it immediately, without any training. Large buttons with a clear design make control a breeze. What's more, you don't have to remove the via GO to fold up your wheelchair. The Invacare retailer will provide guidance for using via Go.

SAFE AND RELIABLE



Alber quality

The via GO push aid incorporates the depth and breadth of experience and expertise of Alber, the specialist and market leader in add-on drives for wheelchairs. Alber products regularly set new benchmarks and stand for enduring reliability and TÜV-tested safety.



MADE IN GERMANY 

Technical data for via GO

Dimensions

Drive unit: 6.7 x 9.4 x 10.6 in (W x H x L)

Control unit: 1.4 x 3.7 x 1.1 in (W x H x D)

Speed 0.9–3.4 mph

Range up to 10 miles
integrated rechargeable
lithium-ion battery*

Weight capacity up to 330 lbs (user weight)**

Gradient up to 16% **

Weight

Charger: 1.19 lbs

Drive unit: 12.87 lbs

Control unit: 0.31 lbs

* Depending on user weight, topography and ambient temperature in accordance with ISO 7176-4

** May be lower depending on the wheelchair used and operator ability




Daily life

Both the control and the installation of the components are particularly uncomplicated and user-friendly with via GO.

The small, ergonomic **control unit** attaches simply to the push handles of your wheelchair.

The clever **mounting device** allows your wheelchair to be folded up as normal, even with the push aid still attached.

As a result of its lightweight design and large handle, the **drive unit** is extremely easy to attach to and remove from the wheelchair.



via GO[®]

The compact push aid for wheelchairs



28.0004.6.02.01 Version 03/2023 © Alber GmbH, Albstadt.
Subject to technical changes.

The via GO is CE compliant according to the European Medical Device Regulation (MDR) 2017/745 and has been tested on a voluntary basis by TÜV Süd according to the latest edition of the product standard EN 12184.

Our quality commitment is your safety! Alber is certified under the current quality management standard ISO 13485 for medical devices. Our environmental management system, which is certified under ISO 14001, ensures the environmental responsibility manufacture of our products.

Invacare Uk
Unit 4 Pencoed Technology Park
Bridgend CF35 5AQ
t: 01656 776222
e: uk@invacare.com

www.invacare.co.uk



LET'S via GO!

Experience the practical push aid at your local Invacare Alber retailer. Use our Retailer Locator to find your nearest stockist.

www.invacare.co.uk/retailer-locator

MADE IN GERMANY