



Invacare Ltd
 Unit 4 Pencoed Technology Park
 Pencoed, Bridgend
 CF35 5AQ
 T: 01656 776222
 E: ordersuk@invacare.com
 www.invacare.co.uk

Please fill in all specifications with each order.

Account code

Customer order number

Company name

Address



e-motion Duo Drive	£5775
Included as standard:	22"
• 2 e-motion drive wheels	Total diameter 550mm
• Control unit (bracket included)	24"
• Speed: 6km/h	Total diameter 600mm
• Weight: Only 7.8 kg per wheel	25" *
• Brushed stainless steel push rims	
• Charger	
• Lithium-ion batteries in the Wheel hubs	
• Pneumatic tyres with puncture protection	

Ergonomic Control System (ECS) 1592486 **£409**

- Switches e-motion on and off
- Two remote assistance levels
- Roll back delay
- Battery monitor
- Self diagnosis



Anti-tippers with jack-up function (pair) 1489214 **£281**
 NB Always use anti-tippers (Alber or existing)

Mobility Plus Package

Download from App Store or Playstore as an In-App purchase

- Remote control
- Increase speed to 8.5km/h
- Cruise control function,
- ECS functions
- Push counter - evaluates your driving behaviour during journeys



Wheelchair details

Wheelchair manufacturer (e.g. Kuschall)

Serial Year

Model (e.g. Compact)

Seat width  cm

Seat height (e.g. Hemi)  cm

Armrest & Brake type

Degree of wheel cam-

(Drum brakes must be removed)

* Total diameter with Right Run 625mm, 615mm with Marathon Plus

Brackets required for installation onto wheelchair

1467178 **STD**
 Omit brackets **-£67**

Fitting on wheelchair at Invacare factory & delivery** **£75**
 Wheelchair order number


Collection & fitting to third party wheelchair & delivery*** **£75+£96**

Fitting by trained dealer*** **£0**

** Only when ordered with an Invacare wheelchair.
 *** Mounting kits can only be sent to selected trained dealers.

Other

- Control Unit RH
- Control Unit LH
- Equipped with vehicle retention system
- Rear wheel extension available
- Available wheel position:



- 1=rear
- 2=centre
- 3=front



Accessories

Push rims

Stainless steel (standard)	1592526	£0
22" or 24" push rim cover, Tetra (pair)	1592527	£133
CarboLife Curve L Black Powder coated	1592528	£289
CarboLife Quadro Tetra-Gripp coated	1592520	£350



Tyres (pairs)

Schwalbe Rightrun (Standard) Black	1592530	£0
Schwalbe Marathon Plus Evolution with SmartGuard puncture protection, Black	1565894	£36
Pro Speed puncture-proof solid PU 1" Black (+ 300g per wheel, not available on 22" or 26")	1592532	£97
Puncture proof with inlay 1 3/8" (+ 500g per wheel, not available on 26")	1592533	£97

Push rim covers (retrofit)

Push rim cover 22"	1592425	£143
Push rim cover 24"	1592421	£143
Push rim cover 26"	1592426	£143

Spoke Covers (pairs)

22" spoke cover Design 'Comic'	1592266	£143
24" spoke cover Design 'Sphere'	1592268	£149
26" spoke cover, Design 'Sphere'	1592269	£154
24" spoke cover, Plain	1592553	£149



Extra options

Brake extension*	1467780	£30
Transport axle	1591362	£68
Ergonomic stubaxle	1591500	£131



Do you want to use both your standard and e-motion wheels in the future? If so, please order the retrofit stub axle.

Manual wheel axle retrofit for e-motion * / **

Retrofit, stub axle to 1/2" diameter **£103**

*If the manual wheels are used in combination with e-motion, these must be fitted with 1/2" stub axles, without drum brakes and 1" wide tyres.

**If drum brakes are installed, these must be removed as they are not compatible with e-motion. Standard lever brakes are used instead.

*maybe required depending on wheelchair type

Comments (mounting etc.)