

# Action<sup>®1</sup> R



A feature rich, entry level chair



Yes, you can.®

The **Action<sup>1</sup> R** is an entry level wheelchair that doesn't compromise on functionality. Built using decades of in-house design experience, the **Action<sup>1</sup> R** will be available in two versions; a self-propel model with 24" spoked rearwheels and aluminium handrims and a transit version with small 12" rearwheels.

With six standard sizes to choose from, it is quick and easy to find the right size. It features easy to remove swingaway armrests, adjustable armpads, push to lock brakes and is quick to fold for storage or transport.

The **Action<sup>1</sup> R** offers:

- **Improved function and comfort**
- **Proven quality from a market leader**

### Improved functionality

The new armrests swing away for transfers and offers two armpad positions - desk and full length - which can be removed quickly and easily for transport.



### Proven reliability and quality

With years of experience and using tried-and-tested technology, the **Action<sup>1</sup> R** is a reliable product that can be issued quickly and with confidence to clients.



### Elevating legrests

The optional Eco elevating legrests provide extra comfort for users.



### Push to lock brake

Easy to use push to lock brake can be easily activated by attendants foot.
































### Anti-tippers

These removable anti-tippers, available as an option, provide extra safety on slopes or when backrest is reclined.

## Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please see your local Invacare website.

	Seat width	Seat depth	Seat angle	Seat height	Backrest angle	Backrest height
Action1 R	 380 / 405 / 430 / 455 / 480 / 505 mm	 400 - 430 mm	 3°	 495 mm	 0° / 7°	 430 mm
	Front frame angle	Length of footrests	Armrest height	Driving unit width <sup>4</sup>	Total width, folded	Total length inc. footrests <sup>4</sup>
Self Propulsion	 70°	 330 - 455 mm	 240 mm	 608 - 700 mm (Seat width +195 mm)	 250 mm	 1085 - 1125 mm
Transit	 70°	 330 - 455 mm	 240 mm	 578 - 670 mm (seat width +160 mm)	 260 mm	 955 mm
	Total length without footrests	Total height	Max. user weight	Total weight	Mass of the heaviest part <sup>3</sup>	
Self Propulsion	 835 - 875 mm	 995 mm	 125 kg	 from 17.5 kg	 from 10,2 kg	
Transit	 710 mm	 995 mm	 125 kg	 from 17.5 kg	 from 10,2 kg	

### Frame colour

Anthracite  
Black



### Upholstery

Black Nylon



#### 1. Stability Values

Max. slope parking brake: 7°  
Static stability downhill: 18°  
Static stability uphill: 9°  
Static stability sideways: 19°

#### 3. Folded Dimensions

Folded height: 735 - 950 mm  
Folded depth: 1025 mm  
Folded width: 250 mm

#### 2. Radius

Mini turning radius: 915 mm  
Pivot width: 1280 mm

#### 4. Overall dimensions

In some configurations, overall width or length exceed the authorised limits, access to emergency escape routes is not possible

T: +44 1656 776222  
uk@invacare.com  
www.invacare.co.uk  
www.invacare.ie

Invacare UK & Ireland  
Unit 4 Pencoed  
Technology Park  
Bridgend CF35 5AQ



**Yes, you can.®**