

# Accent<sup>®</sup>

For everyday care



Accent profiling bed  
Trusted by healthcare professionals



# Affordable quality

The Accent bed provides a simple and affordable quality solution for both Care and Nursing home settings.



### Safety

All Invacare beds meet IEC 60601-2-52 safety standards.



### Quality

Our bed designs and specifications adhere to the highest standards.



### Care-led

Made for the caring environment.



### Comfort

Every Invacare bed promotes rest, dignity and independence.



Accent is a cost effective solution that offers great patient care. Simple aesthetics blend well into residential care settings.



Invacare mattress sold separately.

# Designed with care



With straightforward styling, smooth adjustments and thoughtful details, Accent supports quality care every day.

01



**Affordable choice**  
A simple and affordable quality bed.

02



**Quick & easy setup**  
Simple tool-free assembly and disassembly. Easy to handle.

03



**Adjustable positioning**  
Backrest and thigh sections adjust at the touch of a button.

04



**Height adjustable**  
Adjusts between 40 - 80 cm for comfortable working heights and optimal patient care.

05



**Intuitive hand control**  
Clear visual display for ease of use. Optional lockable tilt available.

06



**Range of side rails**  
Choose from the Slim, Scala 2, Verso HC or Aria options.

08



**Easy to clean**  
Suitable to wipe with mild detergents making it ideal for single or multiple use environments. See User Manual.

07




**Compact power unit & cabling**  
Smaller unit with neater and detachable cables gives more space under the bed and makes moving the bed easier.

09

11°

**Supportive tilt**  
For positioning and pressure management.




# Accent® options


Simplifying care with a choice of everyday options.



Lifting pole



**Slim**  
Full length side rail in Pine wood.



**Scala 2**  
Side rail in steel\*.



**Verso HC**  
Collapsible side rail in steel.



**Aria**  
Full length side rail in steel.




**Bed end**  
Optional panel.



**Alternative handset**  
With lockable tilt for peace of mind.

\* available in a range of heights to accommodate foam and dynamic mattress systems



## Easy to transport

### Optional Transportation Kit

Reduces delivery and set up time. Ideal for large bed placements. Reduced cardboard packaging, great for the environment.



Watch the How to video on assembly & disassembly.



### Customer Support & After Sales

Ask our team. We are always here to help.

Tel: 01656 776222 &  
Tel: 01656 776333

# Technical specifications



Width inside  
90 cm



Width outside  
102 cm



Height adjustment  
40 - 80 cm



Length inside  
208 cm



Length outside  
224 cm



Tilt angle  
11°



Max user weight  
145 kg



Total weight  
64 kg



Heaviest part  
16 kg



Mattress support  
80 - 23 - 30  
- 67 cm



Backrest angle  
0 - 70°



Thigh angle  
0 - 33°

Contact your Authorised  
Distributor:

## Colours



Steel



Beech



Complies to the safety standard,  
IEC 60601-2-52.



ISO 13485 certified company

T: +44 1656 776222  
uk@invacare.com  
www.invacare.co.uk  
www.invacare.ie

Invacare UK & Ireland  
Unit 4 Pencoed Technology  
Park  
Bridgend CF35 5AQ