



Yes, you can:

Tel: 01656 776222
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

July 2026
Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

**M10 TDX SP2 HD 900092 - Weekly Price
£56.64**

Standard Features

- Max user weight 180kg (28 stone)
- Crash tested ISO7176 - 19
- Modulite Seating System
- Matrix Flo-tech Cushion
- Flex 3 adjustable Backrest plates with Laguna backrest cushion
- Standard T-armrests (height and depth adjustable), detachable, with long armpad (375mm)
- Manual settable Seat Tilt 0° - 20°
- Manual settable Backrest Recline 0° - 30° in 3.75° steps
- 80° hangers with fixed footplate
- LiNX DLX-REM110 swing away bracket
- Solid Black tyres
- 6km/h (4mph) 4 pole motors
- 4-Point tie-down system (Q'Straint QLK-150 compatible)
- Lapbelt with autolock (adjustable both sides)

Standard Base price: £8,650



Please make your selection with a cross in the box ☐ Free of charge ☒ standard feature POA special request

Customer weight = KG

Code Description

BASE MODEL*

ARJ0000	☐	TDX SP2 Heavy Duty	£8,650
---------	--------------------------	--------------------	--------

USER WEIGHT*

ARJ1918	☐	User weight less 100 kg	☐
ARJ1919	☐	User weight 101 to 125 kg	☐
ARJ1921	☐	User weight 126 to 136 kg	☐
ARJ1922	☐	User weight 137 to 150 kg	☐

SEATING TYPE*

ARJ2100	☐	Modulite Seat Unit	☒
---------	--------------------------	--------------------	----------------------------------

SEAT WIDTH*				
ARJ2121	☐	380 mm Seat Plate. Can be reduced to 330 mm in combination with inrigged armrests.		☐
ARJ2122	☐	430 mm Seat Plate. Can be reduced to 380 mm in combination with inrigged armrests.		☐
ARJ2123	☐	480 mm Seat Plate		☐
ARJ2124	☐	530 mm Seat Plate		☐
SEAT DEPTH RAIL LENGTH*				
ARJ2129	☐	Standard Depth Package (410mm - 510mm. Simple stepless adjustment from the side)	(NC 2138, 2196)	☐
ARJ2130	☐	Extended Depth Package (510mm - 560mm. Simple stepless adjustment from the side)	(NC 2132, 2134)	☐
SEAT DEPTH SETTING*				
ARJ2132	☐	410 mm	(NC 2196)	☐
ARJ2134	☐	460 mm		☐
ARJ2136	☐	510 mm		☐
ARJ2138	☐	560 mm		☐
SEAT CUSHION*				
ARJ2142	☐	Matrx Flo-tech Cushion (Lite - 75mm up to 127kg) or (Contour - 100mm up to 191kg)		◆
BACKREST SUPPORT*				
ARJ2162	☐	Backposts tall (540 mm) For Flex3 backrests/Tension Adjustable Upholstery/Matrx Backs, measured from seat plate to top of upholstery. Width adjustable cast bracket.		◆
PUSH HANDLES FOR BACKPOSTS*				
ARJ2107	☐	Backposts WITH push handles		◆
BACKREST UPHOLSTERY / CUSHION*				
ARJ2173	☐	Flex3 adjustable Backrest Plates with Laguna Backrest cushion		◆
UPHOLSTERY COLOUR*				
ARJ2152	☐	Black Multi Stretch PU Coated Wipe Down Fabric TR26		☐
ARJ2155	☐	Black Spacer Fabric TR36		☐
ARMRESTS*				
ARJ2191	☐	T-Arm detachable, Standard. Adjustment range from seat plate: 245 - 315 mm		◆
ARMPADS*				
ARJ2196	☐	Long Armpad (375mm)	(NC 2132)	◆
HEADREST / NECKREST (Adequate headrest option must be selected)*				
ARJ1495	☐	Non Invacare Headrest - no mounting hardware. With Backbrace/Pushbar dia. 22mm in conjunction with telescopic seat frame/backposts		◆
MANUAL / POWERED SEATING OPTIONS*				
ARJ1091	☐	Manually settable Seat Tilt 0° - 20°		◆
SEAT TO GROUND HEIGHT*				
ARJ0354	☐	450mm (17.75") (from floor to front of seat plate at 0° tilt)		◆
BACKREST RECLINER*				
ARJ0975	☐	Manual Settable Recliner (0° - 30°, in 3.75° steps) Foldable		◆
LEG REST TYPE*				
ARJ0526	☐	Side Mounted: Standard hanger 80° pair. Lower leg length: 350 - 450 mm		◆
FOOTPLATE TYPE*				
ARJ0512	☐	Fixed Footplate pair		◆

LINX CONTROLS*		
ARJ2201	☐	DLX-REM110 (Drive Only) ◆
CONTROL BRACKETRY*		
ARJ1707	☐	Swing Away RH Bracket ○
ARJ1708	☐	Swing Away LH Bracket ○
SHROUD COLOUR*		
ARJ2640	☐	Obsidian Black Frame Shrouds ○
ARJ2641	☐	Ruby Red Frame Shrouds ○
ARJ2642	☐	Purple Passion Frame Shrouds ○
ARJ2643	☐	Peacock Blue Frame Shrouds ○
ARJ2644	☐	Arctic Blue Frame Shrouds ○
RIM COLOUR*		
ARJ2650	☐	Obsidian Black Frame Shrouds ○
ARJ2651	☐	Ruby Red Frame Shrouds ○
ARJ2652	☐	Purple Passion Frame Shrouds ○
ARJ2653	☐	Peacock Blue Frame Shrouds ○
ARJ2654	☐	Arctic Blue Frame Shrouds ○
FORKS*		
ARJ2632	☐	Double sided fork ○
CENTRE WHEELS (DRIVE WHEELS)*		
ARJ0705	☐	DW 14" (355 mm) x 3" (75 mm), Solid inlay, dual density foam insert ◆
MOTOR TYPE*		
ARJ1282	☐	6 km/h (4mph) (4 pole) ◆
BATTERY*		
ARJ1260	☐	73.5 Ah (Gel) ◆
CHARGER*		
ARJ1343	☐	OptiCharge Charger 8A (200 - 250V) ◆
TRANSPORTATION KIT*		
ARJ1494	☐	4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19). Adequate headrest is compulsory. ◆
FIXATION BELT*		
ARJ2113	☐	Lapbelt (adjustable both sides) adjustable mounting on seat rail extrusion ○

TOTAL PRICE	
--------------------	--

Non Conforming Options:

ARJ2132	410 mm
ARJ2134	460 mm
ARJ2138	560 mm
ARJ2196	Long Armpad (375mm)



For a more complex customer requirement,
please email: orders@invacare.com.au

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):