



Yes, you can:

Tel: 01656 776222
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

July 2026
Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

**M3 TDX SP2 NB 735801 - Weekly Price
£69.07**

Standard Features

- Max User Weight 150kg (23.5 Stone)
- Modulte Seating
- Matrix Flo-tech cushion
- Flex3 Adjustable Plates with Laguna back cushion
- Standard T-armrests (height and depth adjustable), detachable, with long armpad (375mm)
- Seat Lifter with Powered Tilt (300mm, 30°) with CoG shift
- Manual settable Backrest Recliner (0° - 30° in 3.75° steps)
- Manual centre mounted hanger, with angle adjustable footplates
- LiNX DLX-REM216 with swing away bracket
- Solid Black tyres
- 10km/h (6mph) 4 pole motors (including indicators, LED headlights and mirror)
- 4 Point tie-down system (Q'Straint QLK-150 compatible)
- Lapbelt with autolock

Standard Base price: £11,005



Please make your selection with a cross in the box ☐ Free of charge ☒ standard feature POA special request

Customer weight = KG

Code	Description	
BASE MODEL*		
ARM0000	☐ TDX SP2 Narrow Base	£11,005
USER WEIGHT*		
ARM1920	☐ User weight up to 125kg	☐
ARM1921	☐ User weight 126kg to 136kg	☐
ARM1922	☐ User weight 137kg to 150kg	☐
SEATING TYPE*		
ARM2100	☐ Modulte Seat Unit	☒

SEAT FRAME*				
ARM2111		Telescopic Seat Frame + Plate with Side Rails. Width 380 mm to 530 mm. Additional 50 mm width by arm rest adjustment.		◆
SEAT WIDTH*				
ARM2121		380 mm Seat Plate. Can be reduced to 330 mm in combination with inrigged armrests.		○
ARM2122		430 mm Seat Plate. Can be reduced to 380 mm in combination with inrigged armrests.		○
ARM2123		480 mm Seat Plate		○
ARM2124		530 mm Seat Plate		○
SEAT DEPTH RAIL LENGTH*				
ARM2129		Standard Depth Package (410mm - 510mm. Simple stepless adjustment from the side)	(NC 2138)	○
ARM2190		Extended Depth Package (510mm - 560mm. Simple stepless adjustment from the side)	(NC 2132, 2134)	○
SEAT DEPTH*				
ARM2132		410 mm	(NC 2196)	○
ARM2134		460 mm		○
ARM2136		510 mm		○
ARM2138		560 mm		○
SEAT CUSHION*				
ARM2142		Matrx Flo-tech Cushion (Lite - 75mm up to 127kg) or (Contour - 100mm up to 191kg)		◆
BACKREST SUPPORT*				
ARM2162		Backposts tall (540 mm) For Flex3 backrests/Tension Adjustable Upholstery/Matrx Backs, measured from seat plate to top of upholstery. Width adjustable cast bracket.		◆
BACKREST UPHOLSTERY / CUSHION*				
ARM2173		Flex3 adjustable Backrest Plates with Laguna Backrest cushion		◆
UPHOLSTERY COLOUR*				
ARM2152		Black Multi Stretch PU Coated Wipe Down Fabric TR26		○
ARM2155		Black Spacer Fabric TR36		○
ARMRESTS*				
ARM2191		T-Arm detachable, Standard. Adjustment range from seat plate: 245 - 315 mm		◆
ARMPADS*				
ARM2196		Long Armpad (375mm)		◆
HEADREST / NECKREST (Adequate headrest option must be selected)*				
ARM2324		Elan Standard Pad 355mm (10") wide with standard cover complete with Elan Headrest Hardware		◆
MANUAL / POWERED SEATING OPTIONS*				
ARM1094		Seat Lifter with Powered Tilt (300mm, 30°) with CoG shift (Seat height 460mm - 490mm)		◆
BACKREST RECLINER*				
ARM0975		Manual Settable Recliner (0° - 30°, in 3.75° steps) Foldable		◆
LEG REST TYPE*				
ARM0556		Centre Mounted: Manual Centre Mounted Legrest, angle settable. Length 280 - 385 mm		◆
FOOTPLATE TYPE*				
ARM0559		Angle & width adjustable footplates M (140 W x 240 mm D / 5.5"x9.5")		◆

LINX CONTROLS*		
ARM2204	☐	DLX-REM216 (Drive, Actuators & Lights) ◆
CONTROL BRACKETRY*		
ARM1707	☐	Swing Away RH Bracket ○
ARM1708	☐	Swing Away LH Bracket ○
SHROUD COLOUR*		
ARM2640	☐	Obsidian Black Frame Shrouds ○
ARM2641	☐	Ruby Red Frame Shrouds ○
ARM2642	☐	Purple Passion Frame Shrouds ○
ARM2643	☐	Peacock Blue Frame Shrouds ○
ARM2644	☐	Arctic Blue Frame Shrouds ○
RIM COLOUR*		
ARM2650	☐	Obsidian Black Rim Inserts ○
ARM2651	☐	Ruby Red Rim Inserts ○
ARM2652	☐	Purple Passion Rim Inserts ○
ARM2653	☐	Peacock Blue Rim Inserts ○
ARM2654	☐	Arctic Blue Rim Inserts ○
FORKS*		
ARM2632	☐	Double sided fork ◆
CENTRE WHEELS (DRIVE WHEELS)*		
ARM0705	☐	DW 14" (355 mm) x 3" (75 mm), Solid inlay, dual density foam insert ◆
MOTOR TYPE*		
ARM1908	☐	10 km/h (6mph) (includes light system and indicators) ◆
BATTERY*		
ARM1257	☐	50 Ah Gel ◆
CHARGER*		
ARM1343	☐	OptiCharge Charger 8A (200 - 250V) ◆
TRANSPORTATION KIT*		
ARM1494	☐	4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19). Adequate headrest is compulsory. ◆
FIXATION BELT*		
ARM1329	☐	Lap Belt with Auto Lock ○
ARM2113	☐	Lap belt with auto lock (flexible mount in side rails) ◆
REAR VIEW MIRROR*		
ARM1415	☐	Rear view Mirror RH - Must be selected with 10 km/h Motor ○
ARM1339	☐	Rear view Mirror LH - Must be selected with 10 km/h Motor ○

TOTAL PRICE	
--------------------	--

Non Conforming Options:

ARM2132	410 mm
ARM2134	460 mm
ARM2138	560 mm
ARM2196	Long Armpad (375mm)



For a more complex customer requirement,
please email: orders@invacare.com.au

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):