



Yes, you can:

Tel: 01656 776222
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

July 2026
Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

**M4 TDX SP2 NB 1001336 - Weekly Price
£56.71**

Standard Features

- Max User Weight 150kg (23.5 Stone)
- Crash tested ISO7176 - 19
- Modulite Seating
- Matrix Flo-tech cushion
- Basic tension adjustable backrest
- Standard T-armrests (height and depth adjustable), detachable, with long armpad (375mm)
- Manual settable Seat Tilt (0° - 20°)
- Manual settable Backrest Recliner (0° - 30° in 3.75° steps)
- 80° hangers with fixed footplate
- LiNX DLX-REM216 with swing away bracket
- Solid Black tyres
- 10km/h (6mph) motors
- 4 Point tie-down brackets (Q'Straint QLK-150 compatible)
- Lap belt with auto lock (flexible mount in side rails)

Standard Base price: £8,010



Please make your selection with a cross in the box ☐ Free of charge ☒ standard feature POA special request

Customer weight = KG

Code	Description	
BASE MODEL*		
ARM0000	☐ TDX SP2 Narrow Base	£8,010
USER WEIGHT*		
ARM1920	☐ User weight up to 125kg	☐
ARM1921	☐ User weight 126kg to 136kg	☐
ARM1922	☐ User weight 137kg to 150kg	☐
SEATING TYPE*		
ARM2100	☐ Modulite Seat Unit	☒

SEAT FRAME*				
ARM2111		Telescopic Seat Frame + Plate with Side Rails. Width 380 mm to 530 mm. Additional 50 mm width by arm rest adjustment.		◆
SEAT WIDTH*				
ARM2121		380 mm Seat Plate. Can be reduced to 330 mm in combination with inrigged armrests.		○
ARM2122		430 mm Seat Plate. Can be reduced to 380 mm in combination with inrigged armrests.		○
ARM2123		480 mm Seat Plate		○
ARM2124		530 mm Seat Plate		○
SEAT DEPTH RAIL LENGTH*				
ARM2129		Standard Depth Package (410mm - 510mm. Simple stepless adjustment from the side)	(NC 2138)	○
ARM2190		Extended Depth Package (510mm - 560mm. Simple stepless adjustment from the side)	(NC 2132, 2134)	○
SEAT DEPTH*				
ARM2132		410 mm	(NC 2196)	○
ARM2134		460 mm		○
ARM2136		510 mm		○
ARM2138		560 mm		○
SEAT CUSHION*				
ARM2142		Matrx Flo-tech Cushion (Lite - 75mm up to 127kg) or (Contour - 100mm up to 191kg)		◆
BACKREST SUPPORT*				
ARM2162		Backposts tall (540 mm) For Flex3 backrests/Tension Adjustable Upholstery/Matrx Backs, measured from seat plate to top of upholstery. Width adjustable cast bracket.		◆
BACKREST UPHOLSTERY / CUSHION*				
ARM2171		Tension Adjustable Upholstery with Nylon overlay. Short = 480 mm / 7 straps. Long = 540 mm / 8 straps. Measured from seat plate to top of upholstery.		◆
UPHOLSTERY COLOUR*				
ARM2151		Black Nylon Fabric TR18		◆
ARMRESTS*				
ARM2191		T-Arm detachable, Standard. Adjustment range from seat plate: 245 - 315 mm		◆
ARMPADS*				
ARM2196		Long Armpad (375mm)		◆
HEADREST / NECKREST (Adequate headrest option must be selected)*				
ARM1495		Non Invacare Headrest - no mounting hardware. With Backbrace/Pushbar dia. 22mm in conjunction with telescopic seat frame/backposts		◆
MANUAL / POWERED SEATING OPTIONS*				
ARM1091		Manually settable Seat Tilt 0° - 20°		◆
SEAT TO GROUND HEIGHT*				
ARM0354		450mm (17.75") (from floor to front of seat plate at 0° tilt)		◆
BACKREST RECLINER*				
ARM0975		Manual Settable Recliner (0° - 30°, in 3.75° steps) Foldable		◆
LEG REST TYPE*				
ARM0526		Side Mounted: Standard hanger 80° pair. Lower leg length: 350 - 450 mm		◆

FOOTPLATE TYPE*			
ARM0512	☐	Fixed Footplate pair	◆
LINX CONTROLS*			
ARM2204	☐	DLX-REM216 (Drive, Actuators & Lights)	◆
CONTROL BRACKETRY*			
ARM1707	☐	Swing Away RH Bracket	○
ARM1708	☐	Swing Away LH Bracket	○
SHROUD COLOUR*			
ARM2640	☐	Obsidian Black Frame Shrouds	○
ARM2641	☐	Ruby Red Frame Shrouds	○
ARM2642	☐	Purple Passion Frame Shrouds	○
ARM2643	☐	Peacock Blue Frame Shrouds	○
ARM2644	☐	Arctic Blue Frame Shrouds	○
RIM COLOUR*			
ARM2650	☐	Obsidian Black Rim Inserts	○
ARM2651	☐	Ruby Red Rim Inserts	○
ARM2652	☐	Purple Passion Rim Inserts	○
ARM2653	☐	Peacock Blue Rim Inserts	○
ARM2654	☐	Arctic Blue Rim Inserts	○
FORKS*			
ARM2632	☐	Double sided fork	◆
CENTRE WHEELS (DRIVE WHEELS)*			
ARM0705	☐	DW 14" (355 mm) x 3" (75 mm), Solid inlay, dual density foam insert	◆
MOTOR TYPE*			
ARM1908	☐	10 km/h (6mph) (includes light system and indicators)	◆
BATTERY*			
ARM1257	☐	50 Ah Gel	◆
CHARGER*			
ARM1343	☐	OptiCharge Charger 8A (200 - 250V)	◆
TRANSPORTATION KIT*			
ARM1494	☐	4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19). Adequate headrest is compulsory.	◆
FIXATION BELT*			
ARM1329	☐	Lap Belt with Auto Lock	○
ARM2113	☐	Lap belt with auto lock (flexible mount in side rails)	◆
REAR VIEW MIRROR			
ARM1415	☐	Rear view Mirror RH - Must be selected with 10 km/h Motor	○
ARM1339	☐	Rear view Mirror LH - Must be selected with 10 km/h Motor	○

TOTAL PRICE			
--------------------	--	--	--

Non Conforming Options:

ARM2132	410 mm
ARM2134	460 mm
ARM2138	560 mm
ARM2196	Long Armpad (375mm)



For a more complex customer requirement,
please email: orders@invacare.com.au

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):