

Invacare Aquatec Sansibar

Product Data Sheet

Product Range Available



- **Aquatec Sansibar**
Product code: **9.10.101**

Standard Features

- Space saving: mounted on the wall and can be folded when not in use
- Large hygiene recess makes personal hygiene access easier
- Height adjustable support legs for individual adjustment and additional stability
- Compensates for floor unevenness through finely adjustable feet

Technical Specification

- Seat width 400mm
- Seat depth 410mm
- Seat height 460 – 530mm
- Total width 400mm
- Total depth 460mm
- Total height 530 – 600mm
- Graduation, height adjustment 25mm
- Depth folded 110mm
- Height folded 540 – 610mm
- Maximum user weight 130kg (19 Stone)
- 3 Year Warranty

Supporting Documentation

- Independent living catalogue PN: 1560945
-

Manufacturers Contact Details

Invacare Ltd, Unit 4 Pencoed Technology Park, Bridgend, CF35 5AQ

Tel: 01656 776222 E-mail: UK@invacare.com

Web: www.invacare.co.uk