

# Invacare Cadiz shower stool

## Product Data Sheet

---

### Product Range Available



- Cadiz Shower Stool  
Product code: 1544791

### Standard Features

- A high quality and cost effective solution for showering and dressing
- Easy to handle seating height adjustment from 360 to 510mm – in 13 steps
- Ergonomic seat surface
- Stable and secure footprint, due to angled rubber feet
- Cadiz stool can be upgraded to a Cadiz shower chair

### Technical Specification

- Seat width 510mm
- Seat depth 310mm
- Seat height 360 - 510mm
- Total width 440mm
- Maximum user weight 143kg (22.5 Stone)
- 2 Year Warranty

### Supporting Documentation

- Independent living catalogue PN: 1560945
- 

### Manufacturers Contact Details

Invacare Ltd, Unit 4 Pencoed Technology Park, Bridgend, CF35 5AQ

Tel: 01656 776222 E-mail: [UK@invacare.com](mailto:UK@invacare.com)

Web: [www.invacare.co.uk](http://www.invacare.co.uk) / [www.thinkindependence.co.uk](http://www.thinkindependence.co.uk)