



Personalise your
wheelchair with our
UNIQUE service
from Invacare



Based on passion and care, **Invacare UNIQUE** has been launched as a dedicated service that provides expertise in customising an individual’s product to cater for a specific need, or even the most complex requirement.

Whether the modification is a more traditional requirement such as; longer armrests, wider seats or customised controls, or something more personal, such as a hot pink frame, our team has the knowledge, experience and capability to deliver a **UNIQUE** product that meets your specifications.

Contents

Introduction	2
Invacare UNIQUE	3
How to Order	4
Request for Quotation (RFQ) form	5
Medium Active Wheelchairs	6
• Positioning and Comfort	6
• Safety Features	8
• Miscellaneous	9
• Propulsion	10
• Backrest and Upholstery	10
• Personalise your wheelchair	11
Passive Wheelchairs	12
• Positioning and Comfort	12
• Miscellaneous	13
• Safety Features	14
• Options and accessories	15
• Propulsion	16
• Personalise your wheelchair	17
Power Wheelchairs	18
• Seat Accessories	18
• Personalise your wheelchair	19
• Modulite seat accessories	21
Quality	24





..... Invacare UNIQUE

A truly personal specification can sometimes be so specific that these adaptations are **UNIQUE** - this is why we offer the **Invacare UNIQUE** service. Our dedicated team of qualified specialists and engineers will assess every new request. This process can take anywhere from a few minutes (i.e. a new colour), to several weeks when there are more complicated modifications to be made. We always keep **UNIQUE** products on file, traceable by the serial number, so we are able to quickly respond to any queries regarding orders through this service.

We offer three levels of service based on complexity;

Level 1: Direct order - our most simple service. When you submit your completed prescription/order form for your requested wheelchair, add the reference number for your special request, which is found in this brochure to the 'special comments' section. We'll take care of the rest.

- **1st response in 2 days, quote in 4 days**
- **Extra lead time, on average, 1 week**

Level 2: Our 'On Demand' service - this usually requires further information than Level 1. We will ask you to specify the additional information required, either on your prescription form or and/or the Request for Quotation (RFQ) form, for example, the seat width/depth or a special dedicated colour (you will need to provide the RAL reference if you request a specific colour).

- **1st response in 2 days, quote in 1 week**
- **Extra lead time, on average, 2 weeks**

Level 3: Our 'Custom Made' service requires further investigation into the technical implications of modifications, resulting in our team carrying out a full and detailed risk assessment. You will have to fill in the RFQ form for your type of wheelchair. You will also need to attach pictures or drawings of your specific requirements. The more detail supplied, the quicker the quotation process will be.

- **1st response in 2 days, quote in 2 weeks**
- **Extra lead time, on average, 3 weeks**

How to order

Simply follow the easy step by step process to order your **Invacare UNIQUE**:

1. Submit a request

This can be submitted on your Prescription Form with adaptations noted in the **Invacare UNIQUE** section, or if necessary, the Request for Quotation form for your type of wheelchair within this document.

This can be either emailed to specialsuk@invacare.com, or, faxed to **01656 776220**.

Please remember to include;

- Sketches with dimensions
- Your contact number and/or email address
- For a modification to an existing chair, please provide the serial number of the chair
- Any additional information could help us to process your request more quickly

2. Repeat process

If you have ordered the special modification before, please include the previous quotation number and state the quotation details along with the price. This will allow you to bypass the quotation phase and speed up your ordering process.

3. Feasibility study

Once the request is understood it will undergo a risk assessment. If it passes, a price and lead time will be issued. The more complex the request, the longer the lead time* and the more in depth the risk assessment must be.

4. Processing the order

Once the quote has been signed, faxed back or emailed and an order raised, the **UNIQUE** modification will be manufactured and shipped within the stated lead time quoted.

Remember - the more information that can be supplied with your initial request, the better. In some instances this could save time and money with simple alternatives.

**Delivery time frames outlined in quotes are valid for a period of 90 days.
Delivery times on repeated specials may also vary due to workload.*

Request for Quotation (RFQ) form

What is an RFQ form?

Here is an example of a form and what information will be required when you come to ordering your **Invacare UNIQUE** Level 3.

Customised Powered Wheelchair - Required Information
Date 23.06.2016 Revision 1.0 Page 5 of 7

Provider's information

Company
 Customer Account Number
 Address
 Phone
 Email
 Date
 Model


Wheelchair usage:
 Mainly indoor Mainly outdoor Both indoor and outdoor

Required chair dimensions

Position	Description	Dimension	Unit (please state)
A	Overall width	<input type="text"/>	mm / inches
B	Back support angle	<input type="text"/>	deg.
C	Overall height (w/o head support)	<input type="text"/>	mm / inches
D	Overall length	<input type="text"/>	mm / inches
E	Seat to floor height	<input type="text"/>	mm / inches
F	Leg support angle	<input type="text"/>	deg.

Clients information

Name / Reference
 Age
 Disability



Trusted Solutions, Passionate People



Customised Manual Wheelchair - Required Information
Date 23.06.2016 Revision 1.0 Page 2 of 7

Provider's information

Company
 Customer Account Number
 Address
 Phone
 Email
 Date
 Model


Wheelchair usage:
 Mainly indoor Mainly outdoor Both indoor and outdoor
 Self Propelling Attendant Secured in car
 Used as a seat in Transport

Required chair dimensions

Position	Description	Dimension	Unit (please state)
A	Seat width	<input type="text"/>	mm / inches
B	Seat depth	<input type="text"/>	mm / inches
C	Seat to floor height	<input type="text"/>	mm / inches
D	Backrest height	<input type="text"/>	mm / inches
E	Overall length	<input type="text"/>	mm / inches
F	Knee to Heel length	<input type="text"/>	mm / inches

Client's information

Name / Reference
 Age
 Disability



Trusted Solutions, Passionate People



Medium Active Wheelchairs

The following section shows the adaptations available to our Medium Active Wheelchairs and the corresponding **UNIQUE** Levels.



Positioning and comfort

Positioning



Trunk supports

- ▶ Pair of Dartex trunk supports
- ▶ Adjustable height and width
- ▶ Swing away and detachable
- ▶ Swing away armrests will no longer swing away
- ▶ Serial Number needed for spare part order

Available on

- ▶ Action³_{NG}
- ▶ Action⁴_{NG}

Level

- ▶ Level 1

Non Compatibilities

- ▶ Not available with height adjustable, fixed 7° backrest with fixed or height adjustable push handles.

Positioning



Pair of adjustable foot straps

- ▶ Serial Number needed for spare part order DE181226 (+HD)

Available on

- ▶ Action³_{NG}
- ▶ Action⁴_{NG}

Level

- ▶ Level 1

Positioning and comfort

Comfort



Knee pads / adduction pads

- ▶ Pair of vinyl knee pads
- ▶ Adaptable with a clamp on all standard legrests

Available on

- ▶ Action^{®3}_{NG}
- ▶ Action^{®4}_{NG}
- ▶ Action^{®3} Junior
- ▶ Esprit Action^{®4}_{NG}

Level

- ▶ Level 1

Comfort



Pair of amputee legrests for Seat Width 380 and 405mm

- ▶ Smaller size pads
- ▶ Size of pad: 167x142mm

Available on

- ▶ Action^{®3}_{NG}
- ▶ Action^{®4}_{NG}
- ▶ Esprit Action^{®4}_{NG}
- ▶ Action^{®5}

Level

- ▶ Level 1

Comfort



Elevating legrests with calf pads for Seat Width 255, 280 or 305mm

- ▶ Smaller pads to be ordered for one side or both
- ▶ Size of pad: 70mm x 170mm or 82mm x 170mm or 95mm x 170mm

Available on

- ▶ Action^{®3} Junior

Level

- ▶ Level 1

Comfort



Anatomic seat

- ▶ Wooden preformed seat with foam and black vinyl cover
- ▶ Locking system
- ▶ Removable
- ▶ Seat Width 380mm, 405mm and 480mm. Others, except Heavy Duty are available on the Prescription Form
- ▶ Depth adjustable 400mm or 450mm

Available on

- ▶ Action^{®3}_{NG}
- ▶ Action^{®4}_{NG}

Level

- ▶ Level 1

Positioning and comfort

Comfort



Dynamic harness 4 points

- ▶ Available in 3 sizes (small/medium/large)

Available on

- ▶ Action^{®3} Junior

Level

- ▶ Level 1

	Small	Medium	Large
A	11 cm	14,5 cm	17,5 cm
B	18 cm	23 cm	26 cm
C	16 cm	20 cm	23 cm
D	4,5 cm	5 cm	6 cm

Comfort



Abduction pommel

- ▶ Foam and Vinyl
- ▶ Height adjustable
- ▶ Removable

Available on

- ▶ Action^{®3}_{NG}
- ▶ Action^{®4}_{NG}
- ▶ Action^{®3} Junior

Level

- ▶ Level 1

Safety features

Anti-theft



Anti-theft bar

- ▶ Enables the wheelchair to become non-foldable
- ▶ Removable with tools only

Available on

- ▶ Action^{®3}_{NG}
- ▶ Action^{®4}_{NG} (incl. HD)

Level

- ▶ Level 1

Anti-theft



Anti-theft device for legrests

- ▶ The legrests are no longer removable
- ▶ Swing-in and out functions remain available
- ▶ Removable with tools only

Available on

- ▶ Action^{®3}_{NG}
- ▶ Action^{®4}_{NG} (incl. HD)

Level

- ▶ Level 1

Miscellaneous



Gas spring lever on armrests

- ▶ Allows the person to recline the backrest without an attendant

Available on

- ▶ Action^{®3}_{NG}
- ▶ Action^{®4}_{NG}

Level

- ▶ Level 1



Oxygen bottle holder

- ▶ Multi-adjustable
- ▶ Steel
- ▶ To be placed under the seat
- ▶ Max. bottle dimensions: 150mm & Length 500mm

- ▶ 2 sizes: Adjustable from SW380mm to 605mm or from SW305mm to 355mm



Available on

- ▶ Action^{®3}_{NG}
- ▶ Action^{®4}_{NG}
- ▶ Action^{®3} Junior

Level

- ▶ Level 1



Swing away and detachable armrests with additional handle for transfer

Available on

- ▶ Action^{®3}_{NG}
- ▶ Action^{®4}_{NG}

Level

- ▶ Level 1

Propulsion

Propulsion



Dual handrim propulsion

- ▶ Only with 560mm (22") rear wheels. Available for SW280mm, 305mm, 355mm, 380mm

Available on

- ▶ Action^{®3} Junior

Level

- ▶ Level 2

Backrests and upholsteries

Backrest



Recliner and Hub brakes

- ▶ Gas spring reclining backrest with attendant hub brakes
- ▶ Hub brakes lever and reclining lever on the same handle

Available on

- ▶ Action³ Junior
DE173789

Level

- ▶ Level 1

Upholstery



Tension adjustable upholstery for 510mm folding backrest

- ▶ Available with 510mm folding backrest

Available on

- ▶ Action^{®3}_{NG}
- ▶ Action^{®4}_{NG}

Level

- ▶ Level 1

Backrest



Folding backrest 0°

- ▶ Backrest height 430 or 510mm
- ▶ Serial number needed for spare part order

Available on

- ▶ Action^{®3}_{NG}
- ▶ Action^{®4}_{NG}
- ▶ Action^{®3} Junior

Level

- ▶ Level 2

Personalise your wheelchair

Colour



Frame colour (RAL)

- ▶ List of RAL colours available on request
- ▶ Tick special colour option on prescription form

Available on

- ▶ Action^{®3}_{NG}
- ▶ Action^{®4}_{NG}
- ▶ Action^{®3} Junior[®]
- ▶ Esprit Action^{®4}_{NG}

Level

- ▶ Level 2

Colour



Frame colour (from other Action chairs)

- ▶ Any colour that is not on the prescription form but exists on another Action wheelchair, or as spare parts
- ▶ Tick special colour option on prescription form

Available on

- ▶ Action^{®3}_{NG}
- ▶ Action^{®4}_{NG}
- ▶ Action^{®3} Junior
- ▶ Esprit Action^{®4}_{NG}

Level

- ▶ Level 2

Colour



Colour touch

- ▶ Painting on rear wheels support, rims and castor forks
- ▶ List of RAL colours available on request

Available on

- ▶ Action^{®3}_{NG}
- ▶ Action^{®4}_{NG}
- ▶ Action^{®3} Junior
- ▶ Esprit Action^{®4}_{NG}

Level

- ▶ Level 2

Personalise



Personalised backrest upholstery

- ▶ Silk screen print on black nylon upholstery (rear)
- ▶ Cost according to size, colours and quantity

Available on

- ▶ Action^{®3}_{NG}
- ▶ Action^{®4}_{NG}
- ▶ Action^{®3} Junior
- ▶ Esprit Action^{®4}_{NG}

Level

- ▶ Level 3

Passive Wheelchairs

The following section shows the adaptations available to our Passive Wheelchairs and the corresponding **UNIQUE** Levels



Positioning and comfort

Positioning



Shorter seat depth

- ▶ Move the backrest forward

Available on

- ▶ Rea® Azalea®
- ▶ Rea® Clematis®

Level

- ▶ Level 2

Ref. no.

- ▶ UN0014D

Comfort



Higher mounted armrests

Available on

- ▶ Rea® Azalea®

Level

- ▶ Level 3

Ref. no.

- ▶ UN0022D

Miscellaneous

Sandblasting



XLT

- ▶ Sandblasting of XLT Swing. A pure clean Titanium surface without painting

Available on

- ▶ Invacare® XLT

Level

- ▶ Level 1

Ref. no.

- ▶ UN0010D

Options



Options

- ▶ Use of options from different chairs.

Available on

- ▶ Rea® Azalea®
- ▶ Rea® Dahlia®
- ▶ Rea® Clematis®
- ▶ Invacare® XLT

Level

- ▶ Level 1 - 2

Frame Sides



Flat frame sides

Available on

- ▶ Invacare® XLT

Level

- ▶ Level 3

Ref. no.

- ▶ UN0023D

Foot plates



Separate foot plates on a XLT 75°

Available on

- ▶ Invacare® XLT

Level

- ▶ Level 3

Ref. no.

- ▶ UN0024D

Miscellaneous

Camber



6° camber on XLT

Available on

- ▶ Invacare® XLT

Level

- ▶ Level 3

Ref. no.

- ▶ UN0025D

Footbox



Complete footbox

Available on

- ▶ Rea® Azalea® (Std)
- ▶ Rea® Dahlia®
- ▶ Rea® Clematis®

Level

- ▶ Level 3

Ref. no.

- ▶ UN0026D

Safety features

Anti-theft



Anti-theft device for legrests

- ▶ The legrests are not removable however, they can swing in and out

Available on

- ▶ Invacare® XLT

Level

- ▶ Level 2

Ref. no.

- ▶ UN0020D

Options and accessories

Options



Pushbar

- ▶ For chairs with external push handles

Available on

- ▶ Invacare® XLT

Level

- ▶ Level 2

Ref. no.

- ▶ UN0016D

Options



Pushbar

- ▶ For Azalea® Minor

Available on

- ▶ Rea® Azalea®

Level

- ▶ Level 2

Ref. no.

- ▶ UN0017D

Options



Ventilator tray

- ▶ With specific measurements from the customer

Available on

- ▶ Rea® Azalea®
- ▶ Rea® Clematis®

Level

- ▶ Level 2

Ref. no.

- ▶ UN0018D

Propulsion

Propulsion



One arm brake

- ▶ Azalea® and Dahlia®

Available on

- ▶ Rea® Azalea®
- ▶ Rea® Dahlia®
- ▶ Rea® Clematis®
- ▶ Invacare® XLT (Std)

Level

- ▶ Level 2

Ref. no.

- ▶ UN0005D

Propulsion



One arm drive

- ▶ Azalea® (Drive Wheels included)

Available on

- ▶ Rea® Azalea®
- ▶ Rea® Dahlia®
- ▶ Invacare® XLT (Std)


Level

- ▶ Level 2

Ref. no.

- ▶ UN0006D

Propulsion



Shorter spacers for handrim

- ▶ Handrim closer to drive wheel, a more narrow total width of the chair

Available on

- ▶ Rea® Azalea®
- ▶ Rea® Dahlia®
- ▶ Rea® Clematis®
- ▶ Invacare® XLT

Level

- ▶ Level 2

Ref. no.

- ▶ UN0013D

Propulsion



Drive wheels

- ▶ Different combinations of spokes, tyres and handrims

Available on

- ▶ Rea® Azalea®
- ▶ Rea® Dahlia®
- ▶ Rea® Clematis®
- ▶ Invacare® XLT

Level

- ▶ Level 2

Personalise your wheelchair

Colour



Frame colour

- ▶ Special powder from our stock
- ▶ Standard colours from Azalea®, Clematis®, Dahlia® and XLT

Available on

- ▶ Rea® Azalea®
- ▶ Rea® Dahlia®
- ▶ Rea® Clematis®
- ▶ Invacare® XLT

Level

- ▶ Level 2

Ref. no.

- ▶ UN0001D

Colour



Frame colour

- ▶ Special powder not in stock

Available on

- ▶ Rea® Azalea®
- ▶ Rea® Dahlia®
- ▶ Rea® Clematis®
- ▶ Invacare® XLT

Level

- ▶ Level 2

- ▶ UN0002D

Ref. no.

Personalise



Print on Dartex

- ▶ On the backrest upholstery

Available on

- ▶ Rea® Azalea®
- ▶ Rea® Dahlia®
- ▶ Rea® Clematis®
- ▶ Invacare® XLT

Level

- ▶ Level 2

Ref. no.

- ▶ UN0003D

Power Wheelchairs

The following section shows the adaptations available to our Power Wheelchairs and the corresponding **UNIQUE** Levels.



Seat accessories

Seating



Yes, you can.®

Amputee support RH or LH Vari A Base

- ▶ Adding an amputee support instead of the standard legrest (Vari F and fix footplate)

Available on

- ▶ Kite®
- ▶ TDX® SP2
- ▶ Storm®4
- ▶ Spectra® XTR2
- ▶ Pronto™ M41

Level

- ▶ Level 1

Non

Compatibilities

- ▶ Not available with central legrest

Seating



Yes, you can.®

Amputee support RH and LH Vari A Base

- ▶ Adding a pair of amputee supports instead of the standard legrest (Vari F and fix footplates)

Available on

- ▶ Kite®
- ▶ TDX® SP2
- ▶ Storm®4
- ▶ Spectra® XTR2
- ▶ Pronto™ M41

Level

- ▶ Level 1

Non

Compatibilities

- ▶ Not available with central legrest, seat width 380mm not possible, please check centre of gravity with the patient inside the chair

Seat accessories

Seating



Amputee support RH or LH Action

- ▶ Adding an amputee support instead of the standard legrest

Available on

- ▶ Kite®
- ▶ TDx® SP2
- ▶ Storm®4
- ▶ Spectra® XTR2
- ▶ Pronto™ M41

Level

- ▶ Level 1

Non

Compatibilities

- ▶ Not available with central legrest; seat width 380mm might cause problems; standard legrest type may not fit

Seating



Amputee support RH and LH Action

- ▶ Adding a pair of amputee supports instead of the standard legrest

Available on

- ▶ Kite®
- ▶ TDx® SP2
- ▶ Storm®4
- ▶ Spectra® XTR2
- ▶ Pronto™ M41

Level

- ▶ Level 1

Non

Compatibilities

- ▶ Not available with central legrest; seat width 380mm not possible; please check centre of gravity with the patient inside the chair

Personalise your wheelchair

Colour



Shroud colour

- ▶ Customise the shroud of the chair by choosing from any available RAL colour

Available on

- ▶ Kite®
- ▶ TDx® SP2
- ▶ Storm®4
- ▶ Spectra® XTR2
- ▶ Pronto™ M41
- ▶ Fox™

Level

- ▶ Level 1

Non

Compatibilities

- ▶ Other colour codes on request

Colour



Frame colour

- ▶ Customise the main frame of the chair by choosing from any available RAL colour

Available on

- ▶ Kite®
- ▶ Spectra® XTR2
- ▶ Pronto™ M41
- ▶ Fox™

Level

- ▶ Level 1

Non

Compatibilities

- ▶ Other colour codes on request

Personalise your wheelchair

Colour



Rear rim colour

▶ Customise the rear wheel rims of the chair by choosing from any available RAL colour

Available on

- ▶ Kite®
- ▶ Storm®4
- ▶ Spectra® XTR2
- ▶ Pronto™ M41
- ▶ Fox™

Level

- ▶ Level 1

Non Compatibilities

- ▶ Other colour codes on request

Colour



Front rim colour

▶ Customise the front wheel rims of the chair by choosing from any available RAL colour

Available on

- ▶ Kite®
- ▶ Storm®4

Level

- ▶ Level 1

Non Compatibilities

- ▶ Other colour codes on request

Colour



Centre rim colour

▶ Customise the centre wheel rims of the chair by choosing from any available RAL colour

Available on

- ▶ TDX® SP2

Level

- ▶ Level 1

Non Compatibilities

- ▶ Other colour codes on request

Colour



Front sticker colour

▶ For those seating models with front sticker sets, these can be customised by choosing from any available RAL colour

Available on

- ▶ Storm®4

Level

- ▶ Level 1

Non Compatibilities

- ▶ Other colour codes on request

Personalise your wheelchair

Personalise



Special programming

- ▶ Indoor programming of 6 km/h chairs, changed to 3 km/h

Available on

- ▶ Kite®
- ▶ TDx® SP2 / NB
- ▶ Storm®4
- ▶ Spectra® XTR2

Level

- ▶ Level 1

Personalise



Special programming

- ▶ Special customised programming on request

Available on

- ▶ Kite®
- ▶ TDx® SP2
- ▶ Storm®4
- ▶ Spectra® XTR2
- ▶ Pronto™ M41
- ▶ Fox™

Level

- ▶ Level 1

Colour



Special colour combinations

- ▶ Choose from a combination of different colours on the chair, shrouds mixed, or any other customised request

Available on

- ▶ Kite®
- ▶ TDx® SP2
- ▶ Storm®4
- ▶ Spectra® XTR2
- ▶ Pronto™ M41
- ▶ Fox™

Level

- ▶ Level 2

Non Compatibilities

- ▶ On request

Modulite seat accessories

Modulite



Armrest height

- ▶ Adding up to 100mm to the standard armrest height

Available on

- ▶ Kite®
- ▶ TDx® SP2
- ▶ Storm®4
- ▶ Spectra® XTR2
- ▶ Pronto™ M41

Level

- ▶ Level 1

Non Compatibilities

- ▶ Only for T-armrest, due to possible restriction with lateral supports

Modulite seat accessories

Modulite



Pommel

- ▶ Adding a pommel at the frame of the Modulite seat

Available on

- ▶ Kite®
- ▶ TDx® SP2
- ▶ Storm®4
- ▶ Spectra® XTR2
- ▶ Pronto™ M41

Level

- ▶ Level 2

Non Compatibilities

- ▶ Not available with central legrest, small seat width and calf pads

Modulite



Hemiplegic armrest pad fixed (one side)

- ▶ Adding a hemiplegic armrest pad fixed instead of the existing armrest

Available on

- ▶ Kite®
- ▶ TDx® SP2
- ▶ Storm®4
- ▶ Spectra® XTR2
- ▶ Pronto™ M41

Level

- ▶ Level 2

Modulite



Hemiplegic armrest pad swivel style (one side)

- ▶ Adding a hemiplegic armrest pad with swivel mounting bracket instead of the existing armrest

Available on

- ▶ Kite®
- ▶ TDx® SP2
- ▶ Storm®4
- ▶ Spectra® XTR2
- ▶ Pronto™ M41

Level

- ▶ Level 2

Modulite



Hemiplegic armrest pad fixed combined with remote (one side)

- ▶ Adding a hemiplegic armrest pad fixed instead of the existing armrest on same side as the remote

Available on

- ▶ Kite®
- ▶ TDx® SP2
- ▶ Storm®4
- ▶ Spectra® XTR2
- ▶ Pronto™ M41

Level

- ▶ Level 2

Modulite seat accessories

Modulite



Hemiplegic armrest pad swivel style, combined with remote (one side)

- ▶ Adding a hemiplegic armrest pad with swivel mounting bracket instead of the existing armrest on same side as the remote

Available on

- ▶ Kite®
- ▶ TDX® SP2
- ▶ Storm®4
- ▶ Spectra® XTR2
- ▶ Pronto™ M41

Level

- ▶ Level 2

Modulite



Lateral supports on Flex3 tension back

- ▶ Adding laterals from the REA product range to the Flex3 tension back

Available on

- ▶ Kite®
- ▶ TDX® SP2
- ▶ Storm®4
- ▶ Spectra® XTR2
- ▶ Pronto™ M41

Level

- ▶ Level 2

Non Compatibilities

- ▶ Depending on chair configuration, there may be restrictions, however these will be given at quote stage

Modulite



LNx central manual legrest

- ▶ Adding the LNx manual legrest to the Modulite seating

Available on

- ▶ Kite®
- ▶ Storm®4
- ▶ Spectra® XTR2
- ▶ Pronto™ M41

Level

- ▶ Level 2

Non Compatibilities

- ▶ Possible restriction with the footplates, patient information required

Modulite



LNx central electrical legrest

- ▶ Adding the LNx electrical legrest to the Modulite seating

Available on

- ▶ Kite®
- ▶ Storm®4
- ▶ Spectra® XTR2
- ▶ Pronto™ M41

Level


- ▶ Level 2

Non Compatibilities

- ▶ Possible restriction with the footplates, patient information required

Quality

All **UNIQUE** products that have been designed and built by our **Invacare UNIQUE** team are thoroughly tested by experts and comply with all relevant legislation, for your extra piece of mind.



In accordance with regulatory requirements Invacare undertake normal lifetime and durability test methods, to ensure our products meet EN-standard, ISO standards and the Medical Device Directive (MDD*) requirements.

All our products are Class 1 medical devices according to the MDD**. This is also shown on the Declaration of Conformity belonging to your product. Our manufacturing sites are certified according to ISO 13485 - Medical devices - Quality management systems - Requirements for regulatory purposes.

* MDD 93/42/EEC, Annex 1 (Essential Requirements)

** Annex VII

Invacare Limited
Pencoed Technology Park - Pencoed -
Bridgend CF35 5AQ - United Kingdom
Tel: +44 1 656 776 222
Fax: +44 1656 776 220
E-mail: uk@invacare.com
Sales Order E-mail:
ordersuk@invacare.com
www.invacare.co.uk

© 2016 Invacare International Sàrl.
All rights reserved.

All information quoted is believed to be correct
at time of print. Invacare reserves the right
to alter product specifications without prior
consultation.

Project Unique - UK - 06/2016



Follow us on Twitter:
[@Think_Mobility](https://twitter.com/Think_Mobility)



Follow us on YouTube:
Invacare Ltd - UK



Like us on Facebook:
Invacare UK

Trusted Solutions,
Passionate People



Yes, you can.®