

Tel: 01656 776222 Fax: 01656 776220
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

Motability Order Form
Q2 April – June 2020

Customer Ref : Contact No: Delivery Address :	Account No: Order No : Order Qty : Order Date :
-----------------------------------------------------	----------------------------------------------------------

Invacare® Pronto M41 Standard - M1
MOTABILITY
Max User Weight 136kgs (21.4 stone)

Standard Features



- Standard Seating System
- Seat width of seat frame 440mm - 480mm
- Seat cushion (width 430mm, 50mm)
- 460mm seat depth (adjustable to 410mm or 510mm)
- Black Nylon Tension adjustable backrest upholstery
- Detachable armrests, height adjustable (290mm - 360mm without seat cushion) with extended armpad
- Manual settable backrest Recline
- Adjustable flip-up centre footboard
- **LiNX** Control REM110 with standard Bracket
- Candy Red shroud
- 255mm x 80mm (10" x 3") solid insert / Puncture proof Grey drive wheels
- 6km/h (4mph) motors and 35Ah AGM batteries
- 8A (200 - 250V) charger
- Metal buckle lapbelt (adjustable from both sides)

M1 Standard Base Price : £2052



Motability Weekly Rental : £24.81



* Configurations may vary from the image shown

Code	X	Description	◆ Standard feature	○ Free of Charge	Motability
BACKREST HEIGHT (without Seat Cushion)					
AQB0448		480mm (19")	◆		◆
AQB0454		540mm (21.25")		○	○
SEAT HEIGHT					
AQB0325		Seat height 455mm (18")		○	○
AQB0339		Seat height 480mm (19")		○	○
AQB0387		Seat height 505mm (20")		○	○
AQB0338		Seat height 530mm (21")		○	○
AQB0389		Seat height 555mm (22")		○	○

REMOTE & BRACKETRY			
AQB2201	X	REM110 Drive only	◆
AQB1207		Standard Bracket RH	○
AQB1208		Standard Bracket LH	○
AQB1707		Swing Away Bracket RH	£105
AQB1708		Swing Away Bracket LH	£105

Need an option that's not included on this Motability Fixed Specification form?
Please enter Option code, Description and retail price from our standard spec form.
Please note this will be subject to compatibility, this may increase the total price.

Option Code	Description	Option Price