

Customer Ref :	Account No :
Contact No :	Order No :
Delivery Address :	Order Qty :
	Order Date :

## Invacare® Bora Modulite - M12 - M17 - M20 - M60 - M37

MOTABILITY

Max User Weight 130Kg (20.5 Stone)

### Standard Package



### Standard Features

- **Modulite** highly adjustable seating system
- Telescopic adjustable seat widths from 380mm - 530mm (15" - 21") with an additional 50mm, (2") by arms
- Seat depths 410mm - 510mm stepless adjustment (16" - 20")
- Matrix Flo-tech Lite - pressure relieving seat cushion
- Flex3 adjustable backrest plates with Laguna backrest cushion
- Standard T-armrests, detachable with short armpad
- Manual settable seat Tilt (0° - 20°)
- Manual backrest Recline (0° - 30°, in 3.75° steps)
- 80° hangers with fixed footplate
- **LiNX** Control REM110 with standard bracket
- Choice of 2 chassis colours, 7 Body shroud colours
- Silver Rim colour
- 200mm x 45mm (8" x 1.175") Puncture proof PU-foam sealed, Black castors
- 318mm x 57mm (12.5" x 2.25") Puncture proof solid inlay Black drive wheels
- Rear suspension (standard Silver spring)
- 6km/h (4mph) 2 pole motors and 50Ah batteries
- 8A (200 - 250V) charger
- Bolt on or quick release anti-tippers
- 4 Point tie-down system (*Q'Straint QLK-150 compatible*)
- Lapbelt with autolock (adjustable both sides)

**M12 Standard Base Price : £4789**



**Weekly Allowance: £38.95**

Invacare® **Modulite**™  
*An evolution in seating designed around you.*



\* Configurations may vary from the image shown

**Q'STRAINT**  
**QLK-150**  
Compatible

Packages below include the above standard features with a variation to the list, any chargeable options are included in the revised price.



### Changes/ additions to standard package

- White, Silver or Black Rim colour
- LiNX control REM216 with lights and actuators
- 9" Puncture Proof Solid Inlay Black front wheels
- 10km/h (6mph) motors (includes LED lights, fenders and mirror)
- Central kerbrider

**M17 Standard Base Price : £5,508**



**Motability Weekly Rental : £42.41**







### Changes/ additions to standard package

- Headrest, height, depth and sideways adjustable
- Electric Tilt (0° - 20°) seat height 450mm
- Electric Recliner (90° - 120°)
- Central kerbrider

**M20 Standard Base Price : £5,889**



**Motability Weekly Rental : £44.12**







### Changes/ additions to standard package

- Headrest, height, depth and sideways adjustable
- Electric Tilt (0° - 20°) seat height 450mm
- Electric Recliner (90° - 120°)
- Swing away bracket
- Central kerbrider

**M60 Standard Base Price : £5,997**



**Motability Weekly Rental : £44.49**




## Changes/ additions to standard package

- Headrest, height, depth and sideways adjustable
- Seat lifter with electrical Tilt (300mm, 30°) with CoG shift seat height 450mm
- Electric Recliner (90° - 120°)
- Rear suspension (strong Silver spring)
- G-Trac

**M37 Standard Base Price : £7,274**



**Motability Weekly Rental : £51.83**



X		Please make your selection with a cross in the box				◆ Standard feature	○ Free of Charge
Code	Description						Motability
<b>MODULITE SEAT SIZE SETTING CHART:</b>							
		SEAT WIDTHS				<p>Invacare <b>Modulite</b> An evolution in seating designed around you.</p>	
		ASM2121 380mm (15")	ASM2122 430mm (17")	ASM2123 480mm (19")	ASM2124 530mm (21")		
SEAT DEPTHS	ASM2132 410mm (16")	○	○	○	○		
	ASM2134 460mm (18")	○	○	○	○		
	ASM2136 510mm (20")	○	○	○	○		
<b>UPHOLSTERY TYPE</b>							
ASM2151		Black fabric					○
ASM2152		Black Dartex					○
ASM2155		Black spacer fabric					○

REMOTE OPTIONS			
<b>BRACKET OPTION SELECT FOR M12 - M17 - M20 - M29</b>			
ASM1207		Standard Bracket RH	<input type="radio"/>
ASM1208		Standard Bracket LH	<input type="radio"/>
<b>BRACKET OPTION SELECT FOR M60 - M37 ONLY</b>			
ASM1707		Swing Away RH Bracket	<input type="radio"/>
ASM1708		Swing Away LH Bracket	<input type="radio"/>
DRIVE UNIT OPTIONS			
<b>SHROUD COLOUR</b>			
ASM0228		Black	<input type="radio"/>
ASM0252		Arctic White	<input type="radio"/>
<b>BODY SHROUD COLOURS</b>			
ASM0242		Pearl Night Blue	<input type="radio"/>
ASM0243		Metallic Silver	<input type="radio"/>
ASM0224		Black	<input type="radio"/>
ASM0204		Transparent Cherry Red	<input type="radio"/>
ASM0222		Arctic White	<input type="radio"/>
ASM0248		Grasshopper Green	<input type="radio"/>
ASM0279		Sunshine Yellow	<input type="radio"/>
<b>RIM COLOUR SELECT FOR M17 and M29 ONLY</b>			
ASM0290		Rim colour White	<input type="radio"/>
ASM0292		Rim colour Silver	<input type="radio"/>
ASM0291		Rim colour Black	<input type="radio"/>
<b>CASTORS SELECT FOR M17 and M29 ONLY</b>			
ASM0836		9" Castor, Puncture proof solid inlay, Black (non-marking)	<input type="radio"/>
<b>ANTI-TIPPERS</b>			
ASM1234		Bolted Anti-tippers	<input type="checkbox"/>
ASM1201		Quick release	<input type="radio"/>
<b>REAR VIEW MIRROR - Must be selected with 10 km/h Motors (for M17 and M29)</b>			
ASM1415		Rear view Mirror RH	<input type="radio"/>
ASM1339		Rear view Mirror LH	<input type="radio"/>
<p><b>Need an option that's not included on this Motability Fixed Specification form?</b>  <b>Please enter Option code, Description and retail price from our standard spec form.</b>  <b>Please note this will be subject to compatibility, this may increase the total price.</b></p>			
Option Code	Description	Option Price	