

Customer Ref :	Account No :
Contact No :	Order No :
Delivery Address :	Order Qty :
	Order Date :

Invacare® TDX SP2 Narrow Base - M16 - M24 - M51 - M61 - M50
MOTABILITY
Max User Weight 150Kg (23.5 Stone)

Standard Package



Standard Features

- One Piece seat plate width 430mm (adjustable by arms 440mm - 480mm)
- Seat depth 410mm - 560mm stepless adjustment (16" - 22")
- Standard foam cushion 50mm thick
- Backrest frame, detachable, height 480mm
- Tension adjustable backrest
- Standard T-armrests, detachable with short armpad
- Push bar dia 22mm in conjunction with telescopic seat frame / backposts, No headrest - No mounting hardware
- Manual settable seat Tilt (0° - 20°)
- Manual settable backrest Recline (0° - 30° in 3.75° steps)
- 80° hangers with fixed footplate
- LINX control REM110 with standard bracket
- Selection of 10 Stunning shroud and Rim colours
- Centre (driving wheels) 355mm x 76mm Black tyres
- 6km/h (4mph) 4 pole motors and 50Ah (Gel) battery
- 8A (200 - 250V) charger
- 4 Point tie-down system (*Q'Straint QLK-150 compatible*)
- Lapbelt with autolock (adjustable both sides)



* Configurations may vary from the image shown

Q'STRAINT
QLK-150
Compatible

M16 Standard Base Price : £4990

Motability weekly rate £39.66

Packages below include the above standard features with a variation to the list, any chargeable options are included in the revised price.

Changes / additions to standard package

- Modulite seating unit
- Telescopic seat frame with metal seat plate (seat width 380mm - 530mm)
- Matrix Flo-tech Lite cushion
- Backposts long 540mm adjustable backrest plates
- Flex3 adjustable backrest plates with Laguna backrest cushion

M24 Standard Base Price : £5443

Motability Weekly Rental : £42.46








Changes / additions to standard package

- Modulite seating unit
- Telescopic seat frame with metal seat plate (seat width 380mm - 530mm)
- Matrix Flo-tech cushion
- Backposts long 540mm adjustable backrest plates
- Flex3 adjustable backrest plates with Laguna backrest cushion
- **Elan Headrest Standard Pad 255mm (10")**
- Electric Tilt (0°-30°) with CoG shift
- Angle and depth adjustable footplates (adapt)
- LiNX control REM216 with lights and actuators
- 10km/h (6mph) (includes light system with indicators, LED headlights and Mirror)

M51 Standard Base Price : £7113



Motability Weekly Rental : £50.24







Changes / additions to standard package

- Modulite seating unit
- Telescopic seat frame with metal seat plate (seat width 380mm - 530mm)
- Matrix Libra cushion
- Backposts long 540mm adjustable backrest plates
- Flex3 adjustable backrest plates with Laguna backrest cushion
- Without push handles
- Flip-up armrests
- **Elan Headrest Standard Pad 255mm (10")**
- Manual centre mounted hanger with angle adjustable footplate
- Electric Tilt 20° with electric Recliner (0° - 30°)
- Swing away brackets

M61 Standard Base Price : £7185



Motability Weekly Rental : £51.33




Changes / additions to standard package

- Modulite seating unit
- Telescopic seat frame with metal seat plate (seat width 380mm - 530mm)
- Matrix Flo-tech cushion
- Backposts long 540mm adjustable backrest plates
- Flex3 adjustable backrest plates with Laguna backrest cushion
- Seat lifter with electrical Tilt (300mm, 30°) with CoG shift
- Manual centre mounted hanger with angle adjustable footplate
- LiNX control REM216 with lights and actuators
- 10km/h (6mph) (includes light system with indicators, LED headlights and Mirror)

M50 Standard Base Price : £7820



Motability Weekly Rental : £57.26



	X	Please make your selection with a cross in the box		◆ Standard feature	○ Free of Charge		
Code		Description					Motability
MODULITE SEAT SIZE SETTING CHART:							
			SEAT WIDTHS				<p><i>An evolution in seating designed around you.</i></p>
			ARM2121 380mm (15")	ARM2122 430mm (17")	ARM2123 480mm (19")	ARM2124 530mm (21")	
Seat Depth ARM2129	Standard Seat Rail 410mm - 510mm	ARM2132 410mm (16")	○	○	○	○	
		ARM2134 460mm (18")	○	○	○	○	
		ARM2136 510mm (20")	○	○	○	○	
Seat Depth ARM2130	Extended Seat Rail 510mm - 560mm	ARM2136 510mm (20")	○	○	○	○	
		ARM2138 560mm (22")	○	○	○	○	
Seat depth 410mm: cannot combine T-Armrests (2191) with flexible lapbelt mount (2113)							
UPHOLSTERY TYPE (Select for M24, M51, M61 & M50)							
ARM2151		Black fabric					
ARM2152		Black Dartex					○
ARM2155		Black Spacer fabric					○
SEAT HEIGHT (for M61 Only)							
ARM0354		Seat Height setting 450mm					○
ARM0364		Seat Height setting 475mm					○
REMOTE OPTIONS							
BRACKET OPTIONS (Select for M16 - M24 - M51 - M50)							
ARM1207		Standard Bracket RH					○
ARM1208		Standard Bracket LH					○

BRACKET OPTIONS (Select for M61 only)			
ARM1707		Swing Away RH Bracket	<input type="radio"/>
ARM1708		Swing Away LH Bracket	<input type="radio"/>
DRIVE UNIT OPTIONS			
SHROUD COLOURS (ALL COLOURED PARTS ARE DETACHABLE AND CHANGEABLE PLASTIC SHROUDS)			
ARM0228		Wet Black	<input type="radio"/>
ARM2601		SpaceStation Silver	<input type="radio"/>
ARM2602		Island Blue	<input type="radio"/>
ARM2603		Grape JellyBean (Purple)	<input type="radio"/>
ARM2604		Monster Green	<input type="radio"/>
ARM2605		Tangy Orange	<input type="radio"/>
ARM2606		WhiteOut	<input type="radio"/>
ARM2607		Rockstar Red	<input type="radio"/>
ARM2608		Lights Out Black - Matt	<input type="radio"/>
ARM2609		Invacare® Blue	<input type="radio"/>
RIM INSERT COLOURS			
ARM2611		SpaceStation Silver	<input type="radio"/>
ARM2612		Island Blue	<input type="radio"/>
ARM2613		Grape JellyBean (Purple)	<input type="radio"/>
ARM2614		Monster Green	<input type="radio"/>
ARM2615		Tangy Orange	<input type="radio"/>
ARM2616		WhiteOut	<input type="radio"/>
ARM2617		Rockstar Red	<input type="radio"/>
ARM2618		Lights Out Black - Matt	<input type="radio"/>
ARM2619		Invacare® Blue	<input type="radio"/>
ARM2620		Wet Black	<input type="radio"/>
ARM2621		Without Rim Insert	<input type="radio"/>
FIXATION BELT			
ARM1329		Lap Belt (adjustable both sides) fitted to base of backrest	<input type="radio"/>
ARM2113		Lap belt (adjustable both sides) adjustable mounting on seat rail extrusion	Not available with 410mm seat depth <input type="checkbox"/>
REAR VIEW MIRROR - Must be selected with 10km/h Motor (for M50 & M51)			
ARM1415		Rear View Mirror RH	<input type="radio"/>
ARM1339		Rear View Mirror LH	<input type="radio"/>
Need an option that's not included on this Motability Fixed Specification form? Please enter Option code, Description and retail price from our standard spec form. Please note this will be subject to compatibility, this may increase the total price.			
Option Code	Description	Option Price	