

Tel: 01656 776222 Fax: 01656 776220
Email: ordersuk@invacare.com
Online Spares Available at: www.invacare.co.uk

Motability Order Form
Q2 April – June 2020

Customer Ref :	Account No :
Contact No :	Order No :
Delivery Address :	Order Qty :
	Order Date :

Invacare® Kite Modulite - M23 - M38 - M78
MOTABILITY
Max User Weight 160Kg (25 Stone)

Standard Features



- One Piece Seat Plate 430mm (adjustable by arms 440mm - 480mm)
- Standard foam cushion 50mm thick
- Backrest height 480mm
- Upholstery colour Black fabric
- Tension adjustable backrest in Black Nylon fabric
- Standard T-armrests, detachable with short armpad
- Manual settable seat Tilt (0° - 12°)
- Manual backrest Recline (90° - 120°)
- Standard 80° hangers with fixed footplate
- **LINX** control REM216 with lights & actuators with standard bracket 2 Frame colours, 7 body and side shroud colours
- Silver Rim
- 10" Puncture proof solid inlay Grey front wheels
- 14" Puncture proof solid inlay Grey rear wheels
- Dual swing technology suspension system
- Rear suspension with hard shock absorbers
- 10km/h (6mph) motors including LED lights, indicators and Mirror
- 60Ah battery
- 8A (200 - 250 V) charger
- Standard anti-tipper
- 4-Point tie-down system
- Fixation belt with auto-lock



M23 Standard Base Price : £5649

Motability Weekly Rental : £45.05

Invacare® Modulite™
An evolution in seating designed around you.



* Configurations may vary from the image shown

The Package below includes the above standard features with a variation to the list, any chargeable options are included in the revised price.


Changes / additions to standard package

- Modulite seating unit
- Telescopic seat frame with metal seat plate (seat width 380mm - 530mm)
- Matrix Flo-tech cushion
- Headrest, height, depth and sideways adjustable
- Powered Tilt (0° - 20°)

M38 Standard Base Price : £6266



Motability Weekly Rental : £45.85




Changes / additions to standard package

- Modulite seating unit
- Telescopic seat frame with metal seat plate (seat width 380mm - 530mm)
- Matrix Flo-tech cushion
- Backposts 540mm (Tall)
- Flex3 adjustable backrest plates with Laguna backrest plates
- Black upholstery
- Long armpad
- Matrix Flo-tech cushion
- Electrical Tilt (30°) with CoG shift
- Powered Recliner (90° -120°)
- 12km/h (8mph) motors

M78 Standard Base Price : £7371



Motability Weekly Rental : £51.76



	X	Please make your selection with a cross in the box	◆ Standard feature	○ Free of Charge
Code		Description		Motability

MODULITE SEAT SIZE SETTING CHART:

		SEAT WIDTHS				
		AVB2121 380mm (15")	AVB2122 430mm (17")	AVB2123 480mm (19")	AVB2124 530mm (21")	
SEAT DEPTHS	Seat Rail 410mm - 510mm	AVB2132 410mm (16")	○	○	○	○
		AVB2134 460mm (18")	○	○	○	○
		AVB2136 510mm (20")	○	○	○	○

Invacare **Modulite**
An evolution in seating designed around you.



UPHOLSTERY TYPE SELECT FOR M34 AMD M44

AVB2151		Black fabric	<input type="radio"/>
AVB2152		Black Dartex	<input type="radio"/>
AVB2155		Black Spacer fabric	<input type="radio"/>

LINX REMOTE & BRACKETRY

AVB2204	<input checked="" type="checkbox"/>	Remote REM216 Drive, Actuators & Lights	◆
AVB1207		Standard Bracket RH	<input type="radio"/>
AVB1208		Standard Bracket LH	<input type="radio"/>
AVB1707		Swing Away RH Bracket	£105
AVB1708		Swing Away LH Bracket	£105

DRIVE UNIT OPTIONS

FRAME COLOUR

AVB0228		Black	<input type="radio"/>
AVB0252		Arctic White	<input type="radio"/>

BODY SHROUD COLOURS

AVB0242		Body shroud colour Pearl Night Blue	<input type="radio"/>
AVB0243		Body shroud colour Metallic Silver	<input type="radio"/>
AVB0224		Body shroud colour Black	<input type="radio"/>
AVB0204		Body shroud colour Transparent Cherry Red	<input type="radio"/>
AVB0222		Body shroud colour Arctic White	<input type="radio"/>
AVB0248		Body shroud colour Grasshopper Green	<input type="radio"/>
AVB0279		Body shroud Sunshine Yellow	<input type="radio"/>

SIDE SHROUD COLOURS

AVB0284		Side shroud colour Pearl Night Blue	<input type="radio"/>
AVB0285		Side shroud colour Metallic Silver	<input type="radio"/>
AVB0265		Side shroud colour Black	<input type="radio"/>
AVB0261		Side shroud colour Transparent Cherry Red	<input type="radio"/>
AVB0253		Side shroud colour Arctic White	<input type="radio"/>
AVB0250		Side shroud colour Grasshopper Green	<input type="radio"/>
AVB0254		Side shroud colour Sunshine Yellow	<input type="radio"/>

REAR VIEW MIRROR - Must be selected with 10km/h / 12km/h Motor

AVB1415		Rear View Mirror RH	<input type="radio"/>
AVB1339		Rear View Mirror LH	<input type="radio"/>

**Need an option that's not included on this Motability Fixed Specification form?
Please enter Option code, Description and retail price from our standard spec form.
Please note this will be subject to compatibility, this may increase the total price.**

Option Code	Description	Option Price