



Customer Ref :	Account No :
Contact No :	Order No :
Delivery Address :	Order Qty :
	Order Date :

Invacare® Esprit Action⁴ NG Adult - M48
MOTABILITY
Max User Weight : 125Kg (19.7 stone) Double Crossbar

Standard Package

Standard Features

- Alber E-FIX 12" x 2.1/4" propulsion, 12Ah battery
- Seat widths 380mm to 605mm (15" - 24")
- Seat depths 400mm, 450mm or 500mm (16", 18" or 20")
- Seat heights 460mm or 485mm (18" or 19")
- Nylon seat upholstery
- Half folding 7° 430mm backrest
- Contour adjustable Nylon backrest canvas
- 80° swing in / swing out legrest hangers
- Flip up footplates includes heel strap
- Height adjustable, removable full armpad armrests
- Choice of 6 colours
- Anti-tippers with chassis jack and buckle belt

M48 Standard Base Price : £3343

Motability Weekly Rental : £31.81



* Configurations may vary from the image shown

Code	X	Description	◆ Standard feature	○ Free of Charge	Motability
CHASSIS					
SFD0210	X	Double crossbar (max user weight 125kg) adult			○
CONTROL UNIT					
<i>USER REMOTE</i>					
SFD0400		User remote control RH (swing away)	◆		
SFD0410		User remote control LH (swing away)		○	
BATTERY					
SFD0600		Additional battery 12Ah			£240
SEATING					
SFD1040		380mm (15") seat width (double crossbrace, including backbrace)	Not available with 20" seat depth or adjustable seat depth 18"- 20"		○
SFD1050		405mm (16") seat width (double crossbrace, including backbrace)	Not available with 20" seat depth or adjustable seat depth 18"- 20"		○
SFD1060		430mm (17") seat width (double crossbrace, including backbrace)			○
SFD1070		455mm (18") seat width (double crossbrace, including backbrace)			○
SFD1080		480mm (19") seat width (double crossbrace, including backbrace)			○
SFD1090		505mm (20") seat width (double crossbrace, including backbrace)			○
SFD1100		555mm (22") seat width (double crossbrace, including backbrace reinforced seat upholstery & adjustable heel strap)	Not available with 16" seat depth, mistral cushion, foam cushion 400mm / 450mm		£86
SFD1110		605mm (24") seat width (double crossbrace, including backbrace reinforced seat upholstery & adjustable heel strap)	Not available with 16" seat depth, mistral cushion, foam cushion 400mm / 450mm		£91

SEAT DEPTH			
SFD1210		400mm (16") seat depth	Not available with 450mm foam cushion <input type="radio"/>
SFB1220		400mm to 450mm (16" - 18") adjustable seat depth	<input type="radio"/>
SFD1230		450mm (18") seat depth	Not available with 400mm foam cushion <input type="radio"/>
SFB1240		450mm to 500mm (18" - 20") adjustable seat depth	Not available with 400mm foam cushion £27
SFD1250		500mm (20") seat depth	Not available with 400mm foam cushion <input type="radio"/>
SEAT HEIGHT			
SFD1300		460mm (18") seat height	<input type="radio"/>
SFD1310		485mm (19") seat height	<input type="radio"/>
BACKREST UPHOLSTERY TYPE			
SFD2105		Lumbar support (+ strap) padded shaped cover	Not available with 22" or 24" seat width <input type="radio"/>
SFD2110		Contour adjustable (4 + 1 straps) padded shaped cover	<input type="radio"/>
SFD2120		Contour adjustable (10 + 1 straps) padded shaped cover	<input type="radio"/>
SFD2125		Lumbar support (+ strap) BI material padded cover includes back pocket	Not available with 22" or 24" seat width <input type="radio"/>
SFD2126		Contour adjustable (3 + 1Straps) BI material padded cover includes back pocket	Not available with 22" or 24" seat width £47
SFB2130		Contour adjustable (4 + 1 Straps) BI material cover includes back pocket	Not available with 22" or 24" seat width <input checked="" type="radio"/>
SFB2140		Contour adjustable (10 + 1 straps) BI material cover includes back pocket	Not available with 22" or 24" seat width <input type="radio"/>
SEATING OPTIONS / ACCESSORIES			
SFD5022		Foam cushion 50mm Nylon & anti-slip cover	<input type="radio"/>
SEAT CUSHION DEPTH			
SFD5071		For seat depth 400mm	<input type="radio"/>
SFD5072		For seat depth 450mm	<input type="radio"/>
SFD5073		For seat depth 500mm	<input type="radio"/>
CHASSIS			
COLOURS			
SFD6010		Slate Grey	<input type="radio"/>
SFD6030		Pearl White	<input type="radio"/>
SFD6031		Blue Pop	<input type="radio"/>
SFD6032		Hypnotic Purple	<input type="radio"/>
SFD6033		Red Matt	<input type="radio"/>
SFD6034		Matt Black	<input type="radio"/>
Need an option that's not included on this Motability Fixed Specification form? Please enter Option code, Description and retail price from our standard spec form. Please note this will be subject to compatibility, this may increase the total price.			
Option Code	Description	Option Price	