

Tel: 01656 776222 Fax: 01656 776220
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

Motability Order Form
Q2 April – June 2020

Customer Ref :	Account No :
Contact No :	Order No :
Delivery Address :	Order Qty :
	Order Date :

Invacare® Pronto M41 with CAB Seat - M49
MOTABILITY
Max User Weight 136kgs (21.4 stone)

Standard Features



- CAB seating system (M61 style)
- Seat width of seat frame 406mm or 457mm (16" or 18")
- Height and width adjustable flip up armrests
- Headrest height adjustable
- Backrest height 480mm (19")
- Manual settable reclining backrest (0° to +25°)
- Adjustable flip up footboard
- Standard remote Shark2 REM D01 with standard bracket
- Candy Red shroud
- Centre wheel 10" x 3" (255mm x 80mm) solid insert / Puncture proof Grey
- 6km/h (4mph) motor
- Battery 35Ah AGM pair
- Charger 8A (200-250V)
- Metal buckle belt (adjustable from both sides)



* Configurations may vary from the image shown

M49 Standard Base Price : £1820

Motability Weekly Rental : £23.87

Code	X	Please make your selection with a cross in the box	◆ Standard feature	○ Free of Charge	Motability
AQB0000	X	Pronto M41 Basic Model			○
SEAT SIZES					
AQB0341		Seat width 406mm (16"), with seat depth 406mm (16")			○
AQB0346		Seat width 457mm (18"), with seat depth 457mm (18")			○
SEAT HEIGHT					
AQB0364		Seat height 470mm			○
AQB0322		Seat height 495mm			○
AQB0321		Seat height 520mm			○
AQB0323		Seat height 545mm			○
AQB0359		Seat height 570mm			○
REMOTE & BRACKETRY					
AQB1801	X	Standard Remote Shark2 REM D01			◆
AQB1207		Standard Bracket RH			○
AQB1208		Standard Bracket LH			○

Need an option that's not included on this Motability Fixed Specification form?
Please enter Option code, Description and retail price from our standard spec form.
Please note this will be subject to compatibility, this may increase the total price.

Option Code	Description	Option Price