

Tel: 01656 776222 Fax: 01656 776220
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

Motability Order Form

Q2 April – June 2020

Customer Ref :	Account No :
Contact No :	Order No :
Delivery Address :	Order Qty :
	Order Date :

Invacare® Fox in a Box - M54

MOTABILITY
Maximum User Weight 127kg (20 stone)

Standard Features



With Simply Smart
LINX

- Seat width : 400mm (16") or 450mm (18")
- Seat depth : Set at 460mm (18") adjustable to 370mm (14.5")
- Seating : Black Nylon sling seat
- Backrest: Black Nylon sling upholstery
- Backrest height : 490mm (19")
- Detachable T-armrest with width adjustment
- Long Ergonomic armpad (397mm)
- Low seat to floor height 480mm adjustable to 510mm
- Manually settable seat Tilt (0° / 2° / 4° / 6° / 8°) factory setting: 4°
- Backrest Recliner, settable in steps (85° / 95° / 105° / 115°)
- Legrest hangers : Standard hangers 80°
- Footrests : Fixed angle, flip-up footplates
- Control system : LiNX Drive only simply SMART
- Shroud colours: Pearly Night Blue or Cherry Red
- Castor Grey 200mm x 50mm (8" x 2") solid inlay, Silver Rim
- Drive wheel Puncture proof solid inlay, (12 1/2" x 2 1/4"), Silver Rim
- Speed setting : 6km/h (4mph)
- 34Ah battery (range 27km)
- Lap belt with auto-lock (adjustable both sides)
- Manual wheel locks, Charger 8A (200 - 250V)
- 4-Point tie-down (tested version for using wheelchair as a vehicle seat according to ISO 7176-19:2001)



* Configurations may vary from the image shown

Standard Base Price : £2408

Motability Weekly Rental : £26.11

Part number:

Qty Req:

1640339 400mm x 460mm (16" x 18") RH Controller - Red

1640345 450mm x 460mm (18" x 18") RH Controller - Red 48HR DELIVERY

1640346 400mm x 460mm (16" x 18") RH Controller - Blue

1640347 450mm x 460mm (18" x 18") RH Controller - Blue 48HR DELIVERY
