



Yes, you can.®



Tel: 01656 776222 Fax: 01656 776220
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

Motability Order Form

Q2 April – June 2020

Customer Ref :

Contact No :

Delivery Address :

Account No :

Order No :

Order Qty :

Order Date :

Invacare® Pronto M41 Modulite - M57

MOTABILITY

Max user weight 136kg (21.4 stone)

Invacare® Modulite™
An evolution in seating designed around you.

Standard Features



- Modulite seating unit
- Compact, manoeuvrable, **MWD Pronto M41 Base**
- Telescopic adjustable seat widths from 380mm - 530mm (15" - 21") with an additional 50mm (2") by arms
- Seat depths 410mm - 560mm stepless adjustment (16" - 22")
- Matrix Flo-tech Lite seat cushion
- Backposts 540mm (tall)
- Flex3 adjustable backrest plates with Laguna backrest cushion
- Standard T-armrests, detachable with short armpad
- Headrest, height, depth and sideways adjustable
- Seat Lifter with electric Tilt (COG Shift)
- Manual settable Recline (0° - 30°)
- Centre mounted, angle settable legrests
- **LiNX** Control REM211 with standard Bracket
- Candy Red shroud colour
- 255mm x 80mm (10" x 3") solid insert / Puncture proof Grey drive wheels
- 6 km/h (4mph) in-line motors and 35Ah batteries
- 8A (200 - 250V) charger
- 4-Point tie-down system
- Lapbelt with autolock (adjustable both sides)

M57 Standard Base Price : £5422





Motability Weekly Rental : £41.30

Powered Positioning



* Configurations may vary from the image shown

		X	Please make your selection with a cross in the box ◆ Standard feature ○ Free of Charge				
Code	Description					Motability	
MODULITE SEAT SIZE SETTING CHART:							
			SEAT WIDTHS				 
			AQD2121 380mm (15")	AQD2122 430mm (17")	AQD2123 480mm (19")	AQD2124 530mm (21")	
Seat Depth AQD2129	Standard Seat Rail 410mm - 510mm	AQD2132 410mm (16")	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	
		AQD2134 460mm (18")	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	
		AQD2136 510mm (20")	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	
Seat Depth AQD2130	Extended Seat Rail 510mm - 560mm	AQD2136 510mm (20")	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	
		AQD2138 560mm (22")	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	
SEAT CUSHION (sizes according to seat frame settings)							
AQD2142	X	Matrix Flo-tech Lite - pressure relieving seat cushion				◆	
UPHOLSTERY TYPE							
AQD2151		Black fabric				○	
AQD2152		Black Dartex				○	
AQD2155		Black Spacer fabric				○	
DX2 REMOTE & BRACKETRY							
AQD2203	X	REM211 Drive & Actuators				◆	
AQD1207		Standard Bracket RH				○	
AQD1208		Standard Bracket LH				○	
AQD1707		Swing away Bracket RH				£105	
AQD1708		Swing away Bracket LH				£105	
Need an option that's not included on this Motability Fixed Specification form? Please enter Option code, Description and retail price from our standard spec form. Please note this will be subject to compatibility, this may increase the total price.							
Option Code		Description				Option Price	