

Tel: 01656 776222 Fax: 01656 776220
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

Motability Order Form

Q2 April – June 2020

Customer Ref :	Account No :
Contact No :	Order No :
Delivery Address :	Order Qty :
	Order Date :



Standard Package



- Modulite seating system
- Compact, manoeuvrable, **CWD TDX SP2**
- Telescopic adjustable seat widths from 380mm - 530mm (15" - 21") with an additional 50mm, (2") by arms
- Seat depths 410mm - 560mm stepless adjustment (16" - 22")
- Matrix Flo-tech Lite cushion
- Tension adjustable backrest
- Standard T-armrests, detachable with short armpad
- Seat Lifter with powered Tilt (300mm, 30°) CoG Shift
- Choice of seat height either 420mm or 403mm
- Manual settable Recliner 0° - 30° in 3.75° steps
- Manual centre mounted legrest with angle adjustable footplate
- **LiNX** control REM211 with swing away bracket
- Selection of 10 stunning shroud and Rim colours
- Centre (driving wheels) 14" X 3" Black tyres
- Suspension and sure step technology
- 6km/h (4mph) 4 Pole motors
- 8A (200 - 250V) charger
- 4-Point tie-down system
- Lapbelt with autolock (adjustable both sides)

Low!!
403mm or 420mm
seat to ground height

M59 Standard Base Price : £7234



Motability Weekly Rental : £52.38

Invacare **Modulite**
An evolution in seating designed around you.



* Configurations may vary from the image shown


Standard Package

- Backposts WITHOUT push handles
- Elan Headrest Standard Pad
- Flex3 adjustable backrest plates with Laguna backrest cushion
- Standard 80° hangers with angle and depth adjustable footplates (adapt)

M75 Standard Base Price : £7640



Motability Weekly Rental : £53.34



	X	Please make your selection with a cross in the box	◆ Standard feature	○ Free of Charge	
Code		Description			Motability
BASE MODELS					
ARQ0071		TDX SP2 Basic Model Low Rider 73Ah Batteries 420mm seat to ground height			○
ARQ0072		TDX SP2 Basic Model Low Rider 60Ah Batteries 403mm seat to ground height			○

SEATING TYPE : MODULITE MULTI ADJUSTABLE SEATING SYSTEM

MODULITE SEAT SIZE SETTING CHART:

		SEAT WIDTH			
		ARQ2122 430mm (17")	ARQ2123 480mm (19")	ARQ2124 530mm (21")	
Seat Depth ARQ2129	Standard Seat Rail 410mm - 510mm	ARQ2132 410mm (16")	○	○	○
		ARQ2134 460mm (18")	○	○	○
		ARQ2136 510mm (20")	○	○	○
Seat Depth ARQ2130	Extended Seat Rail 510mm - 560mm	ARQ2136 510mm (20")	○	○	○
		ARQ2138 560mm (22")	○	○	○

Invacare® **Modulite**™
An evolution in seating designed around you.



Seat depth 410mm: cannot combine T-Armrests (2191) with flexible lapbelt mount (2113)

UPHOLSTERY TYPE

ARQ2151		Black fabric	◆
ARQ2152		Black Dartex	○
ARQ2155		Black Spacer fabric	○

LINX REMOTE & BRACKETRY

ARQ2203	X	Remote REM211 Drive & Actuator	◆
ARQ1707		Swing away RH Bracket	○
ARQ1708		Swing away LH Bracket	○

SHROUD COLOURS (ALL COLOURED PARTS ARE DETACHABLE AND CHANGEABLE PLASTIC SHROUDS)

ARQ0228		Wet Black	○
ARQ2601		SpaceStation Silver	○
ARQ2602		Island Blue	○
ARQ2603		Grape JellyBean (Purple)	○
ARQ2604		Monster Green	○
ARQ2605		Tangy Orange	○
ARQ2606		WhiteOut	○
ARQ2607		Rockstar Red	○
ARQ2608		Lights Out Black - Matt	○
ARQ2609		Invacare® Blue	○

RIM INSERT COLOURS			
ARQ2611		SpaceStation Silver	○
ARQ2612		Island Blue	○
ARQ2613		Grape JellyBean (Purple)	○
ARQ2614		Monster Green	○
ARQ2615		Tangy Orange	○
ARQ2616		WhiteOut	○
ARQ2617		Rockstar Red	○
ARQ2618		Lights Out Black - Matt	○
ARQ2619		Invacare® Blue	○
ARQ2620		Wet Black	○
ARQ2621		Without Rim insert	○
FIXATION BELT			
ARQ1329		Lapbelt (adjustable both sides) fitted to base of backrest	○
ARQ2113		Lapbelt (adjustable both sides) adjustable mounting on seat rail extrusion	◆
Not available with 410mm seat depth			
Need an option that's not included on this Motability Fixed Specification form? Please enter Option code, Description and retail price from our standard spec form. Please note this will be subject to compatibility, this may increase the total price.			
Option Code	Description	Option Price	