



Yes, you can.®



Tel: 01656 776222 Fax: 01656 776220
Email: ordersuk@invacare.com
Online Spares Available at: www.invacare.co.uk

Motability Order Form

Q2 April – June 2020

Customer Ref :

Contact No :

Delivery Address :

Account No :

Order No :

Order Qty :

Order Date :

Invacare® Fox Modulite - M62 - M63 - M64

MOTABILITY

Max User Weight 127Kg (20 Stone)

Standard Package



Standard Features

- **Modulite** highly adjustable seating system
- Telescopic adjustable seat widths from 380mm - 530mm (15" - 21")
- Seat depths 410mm - 510mm stepless adjustment (16" - 20")
- Tension Adjustable backrest
- Matrix Flo-tech cushion
- Standard T-armrest, detachable with short armpads
- Manual settable seat Tilt (from 0° to 8°)
- Manual settable backrest Recline (0° - 30°)
- Standard 80° hangers with fixed flip-up footplates
- **LiNX** control REM110 with standard bracket
- Choice of 5 colour options
- Solid Grey castors 200mm x 50mm (8" x 2"), Silver Rim
- Grey drive wheels, Puncture proof solid inlay, (12 1/2" x 2 1/4"), Silver Rim (aluminium)
- 6km/h (4mph) motors with 40Ah Gel battery
- 8A (200 - 250V) charger
- 4-Point tie-down system
- Lapbelt with autolock (adjustable both sides)
- Wheel lock / parking brake

M62 Standard Base Price : £2973

Weekly Allowance: £28.68



Invacare **Modulite**
An evolution in seating designed around you.




The Package below include the above standard features with a variation to the list, any chargeable options are included in the revised price.

Changes / additions to standard package

- Flex 3 adjustable backrest plates with Laguna backrest cushion
- Central Kerbrider

M63 Standard Base Price : £3477

Motability Weekly Rental : £31.36




Changes / additions to standard package

- Flex 3 adjustable backrest plates with Laguna backrest cushion
- Headrest, height, depth and sideways adjustable.
- Lateral supports with extended arm.
- Powered Recliner.
- LiNX REM211 Drive & Actuators

M64 Standard Base Price : £4039



Motability Weekly Rental : £34.31



X		Please make your selection with a cross in the box				◆ Standard feature	○ Free of Charge
Code	Description						Motability
MODULITE SEAT SIZE SETTING CHART:							
		SEAT WIDTHS					
		AXB2121 380mm (15")	AXB2122 430mm (17")	AXB2123 480mm (19")	AXB2124 530mm (21")		
SEAT DEPTHS	AXB2132 410mm (16")	○	○	○	○		
	AXB2134 460mm (18")	○	○	○	○		
	AXB2136 510mm (20")	○	○	○	○		
UPHOLSTERY TYPE FOR M63 AMD M64 ONLY							
AXB2151		Black fabric					○
AXB2152		Black Dartex					○
AXB2155		Black Spacer fabric					○
SEAT ADJUSTMENTS							
SEAT HEIGHT SETTING							
AXB0312		Low seat to floor height 480mm at 4° factory setting (-30mm / + 30mm depending on Tilt settings)					○
AXB0344		High seat to floor height 510mm at 4° factory setting (- 30mm / + 30mm depending on Tilt settings)					○
LiNX REMOTE & BRACKETRY							
AXB2201	X	Remote REM110 Drive only					◆
AXB1207		Standard Bracket RH					○
AXB1208		Standard Bracket LH					○
AXB1707		Swing away RH Bracket					£105
AXB1708		Swing away LH Bracket					£105

DRIVE UNIT OPTIONS			
SHROUD COLOUR			
AXB0201		Pearl Night Blue	<input type="radio"/>
AXB0233		Metallic Silver	<input type="radio"/>
AXB0228		Silky Matt Black	<input type="radio"/>
AXB0204		Transparent Cherry Red	<input type="radio"/>
AXB0222		White	<input type="radio"/>
<p>Need an option that's not included on this Motability Fixed Specification form? Please enter Option code, Description and retail price from our standard spec form. Please note this will be subject to compatibility, this may increase the total price.</p>			
Option Code		Description	Option Price