

Tel: 01656 776222 Fax: 01656 776220
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

Motability Order Form
Q2 April – June 2020

Customer Ref :	Account No :
Contact No :	Order No :
Delivery Address :	Order Qty :
	Order Date :

Invacare® Bora "Express" Standard - M68 - M70

MOTABILITY
Max User Weight 130Kg (20.5 Stone)

Standard Features (Spec 1)



- One piece seat plate width 430mm (17") (adjustable by arms 440mm - 480mm, 17" - 19") with standard foam cushion (50mm, 2")
- Seat depth set at 460mm (adjustable 410mm - 510mm)
- Tension adjustable backrest, standard 50mm cushion
- Standard T-armrests, detachable with long armpad
- Adjustable seat Tilt (0° - 20°) & backrest recline (90° - 120°)
- Swing away legrest hangers, fixed footplates
- LiNX REM110 control unit, standard RH bracket (can be swapped to LH)
- 9" Black castors, 12.5" x 2.25" solid Black drive wheels
- Rear suspension
- Black body shrouds, Black Rims
- 6km/h (4mph) motors and 50Ah batteries, 8Ah charger
- Bolt on anti-tippers, 4-Point tie-down system, lapbelt



Configurations may vary to the image shown

Express shipment!

Orders will be shipped **same day** if ordered before 2pm, next day if ordered after 2pm!

Changes / additions to SPEC 1

M68

- LiNX control REM211
- Powered 20° seat Tilt
- 6mk/h (4mph)

Express SPEC 2 RRP : £3,026

Motability Weekly Rental : £29.37

Changes / additions to SPEC 1

M70

- LiNX control REM216
- Powered 20° seat Tilt
- 10km/h (6mph) (includes lights, indicators & mirror)

Express SPEC 4 RRP : £3,725

Motability Weekly Rental : £32.22

Item: 1647648