

Tel: 01656 776222 Fax: 01656 776220
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

Motability Order Form
Q2 April – June 2020

Customer Ref :	Account No :
Contact No :	Order No :
Delivery Address :	Order Qty :
	Order Date :

Invacare® TDX SP2 Narrow Base Express Standard - M71 - M72
MOTABILITY
Max User Weight 150Kg (23.5 Stone)

Standard Features



- One Piece Seat Plate width 430mm (17") (adjustable 440mm - 480mm)
- Seat depth set at 450mm (adjustable 410mm - 510mm)
- Basic Tension adjustable backrest, standard 50mm cushion
- Standard T-armrests, detachable with long armpad
- Manual adjustable seat Tilt (0°-20°)
- Manual backrest Recline (0°- 30° in 3.75° steps)
- Standard 80° with fixed Footplates
- LiNX REM110 control unit, Standard RH Bracket (can be swapped to LH)
- 6" Black castors, 14" solid Black drive wheels
- Sure Step and Stability Lock
- Matt Black body shrouds, Black rims
- 6 km/h (4mph) motors and 50Ah Batteries, 8Ah charger
- 4 Point tie down system, lapbelt



Configurations may vary to the image shown

Express shipment!

Orders will be shipped **same day** if ordered before 2pm, next day if ordered after 2pm!

6km/h (4mph)

10km/h (6mph)

Please see Standard Features Above **M71**

Express SPEC 1 RRP : £3,982

Motability Weekly Rental : £34.37

Item: 1650309

Changes / additions to SPEC 1 **M72**

- LiNX control REM216
- 10km/h (6mph) (includes lights, indicators & mirror)

Express SPEC 2 RRP : £4,529

Motability Weekly Rental : £37.19

Item: 1650311