

Tel: 01656 776222 Fax: 01656 776220  
Email: ordersuk@invacare.com  
Online Spares available at: www.invacare.co.uk

Motability Order Form  
Q2 April – June 2020

Customer Ref :	Account No :
Contact No :	Order No :
Delivery Address :	Order Qty :
	Order Date :

**Invacare® TDX SP2 Narrow Base Express Modulite - M73 - M74**  
**MOTABILITY**  
**Max User Weight 150Kg (23.5 Stone)**

**Standard Features**



- Telescopic adjustable seat width 430mm (17") (adjustable 440mm - 480mm)
- Seat depth set at 460mm (adjustable 410mm - 510mm)
- Flex3 adjustable backrest, Matrix Flo-tech Lite cushion
- Standard T-armrests, detachable with long armpad
- Powered Tilt (0°-20°)
- Manual settable backrest Recline (0°- 30° in 3.75° steps)
- Standard 80° hangers with fixed Footplates
- LiNX REM110 control unit, Standard RH Bracket (can be swapped to LH)
- 6" Black castors, 14" solid Black drive wheels
- Sure Step and Stability Lock
- Matt Black body shrouds, Black rims
- 6 km/h (4mph) motors and 50Ah Batteries, 8Ah charger
- 4 Point tie down system, lapbelt



Configurations may vary to the image shown

**Express shipment!**

Orders will be shipped **same day** if ordered before 2pm, next day if ordered after 2pm!

**6km/h (4mph)**

**10km/h (6mph)**

Please see Standard Features Above **M73**

**Express SPEC 3 RRP : £ 5786**

**Motability Weekly Rental : £43.65**

Item: 1650310

Changes / additions to SPEC 3 **M74**

- LiNX control REM216
- 10km/h (6mph) (includes lights, indicators & mirror)

**Express SPEC 4 RRP : £6332**

**Motability Weekly Rental : £46.50**

Item: 1650312