

Tel: 01656 776222 Fax: 01656 776220
 Email: ordersuk@invacare.com
 Online Spares available at: www.invacare.co.uk

Motability Order Form

Q2 April – June 2020

Customer Ref :	Account No :
Contact No :	Order No :
Delivery Address :	Order Qty :
	Order Date :

Invacare® Storm 4 Modulite® - M76
 MOTABILITY
 Max User Weight 150Kg (23.5 Stone)

Standard Features



- Modulite highly adjustable seating system
- Telescopic adjustable Seat Widths from 330mm - 530mm (13" - 21") with an additional 50mm, (2") by arms
- Seat Depths 410mm - 510mm Stepless Adjustment (16" - 20")
- Matrix Flo-tech Lite Cushion
- Tension Adjustable Backrest
- Standard T-Armrests, detachable with Short armpad
- Manual settable seat Tilt 0° - 20°
- Manual settable Recliner (0° - 30° in 3.75° steps)
- VARI F hangers with Fixed Footplate
- **LiNX** Control REM216 with standard bracket
- Choice of 6 Colours
- Puncture Proof Castors
- Puncture Proof Drive Wheels
- 73.5 Ah Batteries
- 12km/h (8mph) Motors
- 8A (200 - 250V) Charger
- LED light system with indicators and Mirror
- 4 Point tie down system
- Lapbelt with autolock (adjustable both sides)
- Anti-Tippers
- Rear wheel Suspension

M76 Standard Base Price : £6679

Motability Weekly Rental : £47.67



Invacare® **Modulite™**
 An evolution in seating designed around you.



* Configurations may vary from the image shown

Q'STRAIT.
QLK-150™
 Compatible

	X	Please make your selection with a cross in the box	◆ Standard feature	○ Free of Charge	
Code		Description			Motability

MODEL

ATB0000	X	Storm 4 Basic Model			◆
---------	---	---------------------	--	--	---

SEAT SIZES

MODULITE SEAT SIZE SETTING CHART

		SEAT WIDTHS					
		330mm* (13")	ATB2121 380mm (15")	ATB2122 430mm (17")	ATB2123 480mm (19")	ATB2124 530mm (21")	
SEAT DEPTHS	Standard Seat Rail 410 - 510mm	ATB2132 (16") 410mm	○	○	○	○	○
		ATB2134 (18") 460mm	○	○	○	○	○
		ATB2136 (20") 510mm	○	○	○	○	○



Invacare® **Modulite™**
 An evolution in seating designed around you.

HEADREST/NECKREST (Dartex) - Adequate headrest option must be selected		
ATB1432		Headrest sidewise adjustable £179
ATB1495		Push bar dia. 22mm in conjunction with telescopic seat frame / backposts, No Headrest - no mounting hardware. ○
SEAT HEIGHT SETTINGS		
ATB1023		Fixed Seat Height 450mm (17.75") ◆
ATB1025		Fixed Seat Height 480mm (19") ○
LINX REMOTE & BRACKETRY		
ATB2204	X	REM216 Drive, Actuators & Lights ◆
ATB1207		Standard Bracket RH ○
ATB1208		Standard Bracket LH ○
ATA1707		Swing Away RH Bracket £105
ATA1708		Swing Away LH Bracket £105
CHASSIS OPTIONS		
COLOURS		
ATB0202		Blue ○
ATB0228		Black ○
ATB0220		Champagne ○
ATB0204		Transparent Cherry Red ○
ATB0222		White ○
ATB0216		Olive Green ○
REAR SUSPENSION		
ATB1240		Suspension low (up to 80kg user weight) ○
ATB1241		Suspension middle (from 70kg up to 100kg user weight) ○
ATB1242		Suspension hard (from 90kg up to 150kg user weight) ○
REAR VIEW MIRROR - Must be selected with 12km/h motors		
ATB1415		Rear view Mirror RH ○
ATB1339		Rear view Mirror LH ○
Need an option that's not included on this Motability Fixed Specification form? Please enter Option code, Description and retail price from our standard spec form. Please note this will be subject to compatibility, this may increase the total price.		
Option Code	Description	Option Price