



Yes, you can:

Powerchairs Retail Prescription Form

Tel: 01656 776222 Fax: 01656 776220
 Email: ordersuk@invacare.com
 Online Spares available at: www.invacare.co.uk

Apr - Jun 2020
 Customer Ref:

| | |
|--------------------------|--------------------|
| Customer Ref* | Account No* |
| Contact No* | Order No* |
| Delivery Address* | Order Qty* |
| | Order date* |

M83 AVIVA RX20 Modulite

Standard Features

- Max user weight 136kg (21 Stone)
- Crash tested ISO 7176 -Part 19
- Modulite highly adjustable seating system
- Telescopic adjustable seat widths from 330mm -530mm (13" -21") with an additional 50mm, (2") by arms
- Seat depths 410mm - 510mm stepless adjustment (16" - 20")
- Flex3 adjustable Backrest Plates with Laguna Backrest cushion
- Standard T-Armrests, detachable with short armpad
- Headrest sideways adjustable
- Powered "fixed pivot" Tilt 20°
- Choice of seat height's
- Powered Recliner 30°
- 80° hangers with fixed footplate
- Remote LiNX DLX-REM216 with standard bracket
- Choice of 5 Shroud colours, 5 Rim insert colours
- 8" (200 mm) x 2.50-4 (65 mm) Castor, puncture proof solid tyre
- 14" (355 mm) x 3" (75 mm), Solid inlay black drive wheels
- CTC Suspension Control - Traction - Comfort
- 10 km/h 4 pole Durawatt motors Package combined with LED Lights
- 52Ah batteries
- 8A (200 - 250V) charger
- Anti tippers
- 4 Point tie down system
- Lapbelt with autolock (adjustable both sides)



Standard Base price: £7,446

Motability Weekly Rental: £55.13

Please make your selection with a cross in the box Free of charge standard feature POA special request

Customer weight = KG

| Code | Description | |
|-----------------------|---|-------------------------------------|
| BASE MODEL* | | |
| JBC0000 | AVIVA RX 20 Base Model Max user weight up to 136 kg | £7,446 |
| JBX1300 | Special Extended lead time: + 10 days | <input type="radio"/> |
| CHASSIS WIDTH* | | |
| JBL2810 | Chassis Narrow (560 mm/22") | <input checked="" type="checkbox"/> |

| | | | |
|---|--------------------------|---|---|
| SEATING TYPE* | | | |
| JBC2100 | <input type="checkbox"/> | Invacare® Modulite Seat Unit | ◆ |
| SEAT FRAME* | | | |
| JBC2111 | <input type="checkbox"/> | Telescopic Seat Frame and Plate with Side Rails | ◆ |
| SEAT WIDTH* | | | |
| JBC2121 | <input type="checkbox"/> | 380 mm Seat Plate | ○ |
| JBC2122 | <input type="checkbox"/> | 430 mm Seat Plate | ○ |
| JBC2123 | <input type="checkbox"/> | 480 mm Seat Plate | ○ |
| JBC2124 | <input type="checkbox"/> | 530 mm Seat Plate (NC 2106) | ○ |
| Body Fit Kit | | | |
| JBC2106 | <input type="checkbox"/> | Body fit Kit | ○ |
| SEAT DEPTH* | | | |
| JBC2132 | <input type="checkbox"/> | Seat Depth 410 mm | ○ |
| JBC2134 | <input type="checkbox"/> | Seat Depth 460 mm | ○ |
| JBC2136 | <input type="checkbox"/> | Seat Depth 510 mm | ○ |
| SEAT CUSHION* | | | |
| JBC2142 | <input type="checkbox"/> | Matrx Flo-tech Cushion Flo-tech Lite | ◆ |
| BACKREST SUPPORT* | | | |
| JBC2162 | <input type="checkbox"/> | Backposts tall (540 mm) | ◆ |
| PUSH HANDLES FOR BACKPOSTS | | | |
| JBC2108 | <input type="checkbox"/> | Backposts without push handles | ◆ |
| BACKREST UPHOLSTERY / CUSHION* | | | |
| JBC2173 | <input type="checkbox"/> | Flex3 adjustable Backrest Plates with Laguna Backrest cushion | ◆ |
| UPHOLSTERY TYPE* | | | |
| JBC2151 | <input type="checkbox"/> | Black Fabric Upholstery | ○ |
| JBC2152 | <input type="checkbox"/> | Black Dartex Upholstery | ○ |
| JBC2155 | <input type="checkbox"/> | Black Spacer Fabric Upholstery | ○ |
| ARMREST (Pair)* | | | |
| JBC2191 | <input type="checkbox"/> | T-Arm detachable, Standard | ◆ |
| ARMPAD (Pair)* | | | |
| JBC2195 | <input type="checkbox"/> | Short Armpad (325 mm) | ◆ |
| HEADREST / NECKREST - Adequate headrest option must be selected* | | | |
| JBC1432 | <input type="checkbox"/> | Headrest sideways adjustable | ◆ |
| SEAT LIFTER / SEAT TILT* | | | |
| JBC1092 | <input type="checkbox"/> | Powered "fixed pivot" Tilt 20° | ◆ |
| SEAT TO GROUND HEIGHT* | | | |
| JBC2857 | <input type="checkbox"/> | Seat-to-Floor Height 435 mm | ○ |
| JBC2858 | <input type="checkbox"/> | Seat-to-Floor Height 460 mm | ○ |
| JBC2859 | <input type="checkbox"/> | Seat-to-Floor Height 485 mm | ○ |
| BACKREST RECLINER* | | | |
| JBC0977 | <input type="checkbox"/> | Powered Recliner 30° | ◆ |
| FOOTREST HANGER TYPE* | | | |
| JBC0526 | <input type="checkbox"/> | Standard hanger 80° pair Lower leg length: 350 - 450 mm | ◆ |
| FOOTPLATE TYPE* | | | |
| JBC0512 | <input type="checkbox"/> | Fixed Footplate pair | ◆ |

| CONTROLS* | | | |
|----------------------------|--------------------------|--|---|
| JBC2204 | <input type="checkbox"/> | Remote LiNX DLX-REM216 = Drive, Actuators & Lights 3 Drive functions and power module incl. 2 Actuators & Lights | ◆ |
| CONTROL BRACKETRY* | | | |
| JBC1207 | <input type="checkbox"/> | Standard Bracket RH | ○ |
| JBC1208 | <input type="checkbox"/> | Standard Bracket LH | ○ |
| SHROUD COLOUR* | | | |
| JBC2871 | <input type="checkbox"/> | Mid range - Black/Anthracite Shroud | ○ |
| JBC2872 | <input type="checkbox"/> | Mid range - Silver Shroud | ○ |
| JBC2873 | <input type="checkbox"/> | Mid range - Bengal Red Shroud | ○ |
| JBC2874 | <input type="checkbox"/> | Mid range - Mercury Blue Shroud | ○ |
| JBC2875 | <input type="checkbox"/> | Mid range - Midnight Purple Shroud | ○ |
| RIM COLOUR* | | | |
| JBC2882 | <input type="checkbox"/> | Mid range - Black Rim Inserts | ○ |
| JBC2883 | <input type="checkbox"/> | Mid range - Silver Rim Inserts | ○ |
| JBC2884 | <input type="checkbox"/> | Mid range - Bengal Red Rim Inserts | ○ |
| JBC2885 | <input type="checkbox"/> | Mid range - Mercury Blue Rim Inserts | ○ |
| JBC2886 | <input type="checkbox"/> | Mid range - Midnight Purple Rim Inserts | ○ |
| FORKS* | | | |
| JBC2632 | <input type="checkbox"/> | Double sided fork Without castor rim insert (enlarges wide base to 656 mm) | ◆ |
| FRONT WHEEL* | | | |
| JBC0828 | <input type="checkbox"/> | 8" (200 mm) x 2.50-4 (65 mm) Castor, puncture proof solid tyre | ◆ |
| DRIVE WHEEL* | | | |
| JBC0705 | <input type="checkbox"/> | DW 14" (355 mm) x 3" (75 mm), Solid inlay Dual density foam insert | ◆ |
| C.T.C SUSPENSION* | | | |
| JBC2870 | <input type="checkbox"/> | CTC Suspension Control - Traction - Comfort | ◆ |
| SUSPENSION SETTING* | | | |
| JBC1910 | <input type="checkbox"/> | Suspension Soft | ◆ |
| JBC1911 | <input type="checkbox"/> | Suspension Medium | ○ |
| JBC1915 | <input type="checkbox"/> | Suspension Hard | ○ |
| MOTOR TYPE* | | | |
| JBC1908 | <input type="checkbox"/> | 10 km/h Package Always combined with LED 1437 | ◆ |
| BATTERY TYPE* | | | |
| JBC2860 | <input type="checkbox"/> | 52 Ah AGM (Top-Gel) | ◆ |
| CHARGER* | | | |
| JBC1343 | <input type="checkbox"/> | Opti Charge Charger 8A (IP44) | ◆ |
| ANTI-TIPPERS* | | | |
| JBC1205 | <input type="checkbox"/> | Anti-tipper | ◆ |
| KERB RIDER | | | |
| JBC1244 | <input type="checkbox"/> | Kerb Rider With park position feature | ◆ |
| TRANSPORTATION KIT* | | | |
| JBC1494 | <input type="checkbox"/> | 4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19) Adequate headrest is compulsory | ◆ |
| LAP BELT* | | | |
| JBC2113 | <input type="checkbox"/> | Lap belt with Auto Lock Flexible mount in side rails | ◆ |

REAR VIEW MIRROR*

| | | | |
|---------|--------------------------|---------------------|-----------------------|
| JBC1415 | <input type="checkbox"/> | Rear view mirror RH | <input type="radio"/> |
| JBC1339 | <input type="checkbox"/> | Rear view mirror LH | <input type="radio"/> |

| | | | |
|--------------------|--|--|----------------------|
| TOTAL PRICE | | | <input type="text"/> |
|--------------------|--|--|----------------------|

Non Conforming Options:

JBC2106

Body fit Kit



For a more complex customer requirement please refer to our Invacare Unique special order form. email specialsuk@invacare.com or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: _____ Email: _____

Phone: _____ Fax: _____

Description of special request (please include relevant option code or part number):