



Yes, you can:

Tel: 01656 776222
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

July 2026
Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

Storm 4 Modulite

Standard Features

- Max User Weight 150kg (23.5 stone)
- Crash tested ISO7176 - 19
- Highly adjustable Modulite Seating System
- Telescopic adjustable seat widths from 330mm - 530mm (13" - 21") with an additional 50mm, (2") by arms
- Seat depths from 410mm - 510mm (16" - 20"), stepless adjustment
- Tension adjustable backrest
- Standard T-Armrests, height and depth adjustable, detachable, with short armpads
- Manual settable Seat Tilt
- Manual settable Backrest Recline (0° - 30° in 3.75° steps)
- VARI F leg rests with fixed footplates
- LiNX DLX-REM110 with standard bracket
- Choice of 3 colours
- Puncture Proof Castors
- Puncture Proof Drive Wheels
- 73.5 Ah Batteries
- 6 km/h (4mph) Motors
- 8A (200 - 250V) Charger
- 4 Point tie down system
- Lapbelt with autolock (adjustable both sides)
- Anti-Tippers
- Rear wheel suspension

Standard Base price: £7,250



Please make your selection with a cross in the box ○ Free of charge ◆ standard feature POA special request

Customer weight = KG

Code	Description	
BASE MODEL*		
ATB0000	Storm 4 Modulite	£7,250
ATZ1300	Special Storm 4 Modulite (Extended Lead Time). Please provide details in comments section.	POA
SEATING TYPE*		
ATB2100	Modulite Seat Unit	◆

SEAT FRAME*			
ATB2111		Telescopic Seat Frame + Plate with Side Rails. Width 380 mm to 550 mm. Additional 50 mm width by arm rest adjustment.	◆
ATB2112		Telescopic Seat Frame + Sling Seat with Side Rails. Width 380 mm to 530 mm. +30 mm seat height. Additional 50 mm width by arm rest adjustment.	○
SEAT WIDTH*			
ATZ1001		Seat width 330mm (13") seat width achieved on a 380mm seat frame (2121) with inrigged armrests (2183)	○
ATB2121		380 mm Seat Plate. Can be reduced to 330 mm in combination with inrigged armrests. (NC 2703, 2705)	○
ATB2122		430 mm Seat Plate. Can be reduced to 380 mm in combination with inrigged armrests. (NC 2703, 2705)	○
ATB2123		480 mm Seat Plate (NC 2704, 2706)	○
ATB2124		530 mm Seat Plate (NC 2175, 2176, 2178, 2702, 2704, 2706)	○
ATB2125		550 mm Seat Plate (NC 2134, 2142, 2143, 2147, 2144, 1980, 1981, 1982, 1983, 1980, 1981, 1982, 1983, 1985, 2179, 2172, 2173, 2174, 2175, 2176, 2178, 2168, 2169, 2380, 2381, 2382, 2383, 2385, 2386, 2387, 2389, 2390, 2391, 2393, 2394, 2395, 2151, 2152, 2183, 2197, 2184, 2186, 2188, 2185, 2187, 2189, 1496, 2701, 2702, 2703, 2704, 2705, 1404, 1405)	○
SEAT DEPTH*			
ATB2132		410 mm (NC 2701, 2702, 2703, 2704, 2705, 2706)	○
ATB2134		460 mm	○
ATB2136		510 mm	○
SEAT OPTION			
ATBSWO		"H" Frame Interface Kit	£495
SEAT CUSHION*			
ATB2142		Matrx Flo-tech Lite Cushion	£120
ATB2143		Matrx Flo-tech Contour Visco; (2-way stretch material) (NC 1981, 1982, 1983, 1980, 1981, 1982, 1983, 1985)	£169
ATB2147		Matrx Libra Cushion, 90 mm thick (Extended lead time)	£410
ATB2144		Matrx Vi Cushion (Mesh fabric) (NC 1981, 1982, 1983, 1980, 1981, 1982, 1983)	£400
ATB2148		Standard Cushion (90mm high) - For 550mm Seat Plate Only (NC 1981, 1982, 1983, 1980, 1981, 1982, 1983, 1985)	○
ATB2140		Without seat cushion (NC 1981, 1982, 1983, 1980, 1981, 1982, 1983, 1985)	○
1st POSITIONING ACCESORIES FOR MATRX LIBRA			
ATB1980		Tapered Wedges (Pair)	£30
ATB1981		Abductor	£23
ATB1982		Lateral Wedge	£20
ATB1983		Anterior Seat Wedges	£27
2nd POSITIONING ACCESORIES FOR MATRX LIBRA			
ATB1980		Tapered Wedges (Pair)	£30
ATB1981		Abductor	£23
ATB1982		Lateral Wedge	£20
ATB1983		Anterior Seat Wedges	£27
POSITIONING ACCESSORIES FOR MATRX VI			
ATB1985		Pelvic Obliquity Kit (Contains 2x positioning accessories)	£54
BACKREST SUPPORT*			
ATB2162		Backposts tall (540 mm) For Flex3 backrests/Tension Adjustable Upholstery/Matrx Backs, measured from seat plate to top of upholstery. Width adjustable cast bracket. Push Handles included.	◆
ATB2163		Backposts short (480 mm). For Tension Adjustable Upholstery/Matrx Backs, measured from seat plate to top of upholstery. Width adjustable cast bracket. Push Handles included. (NC 2173)	○

PUSH HANDLES FOR BACKPOSTS*			
ATB2107		Backposts WITH push handles	◆
ATB2108		Backposts WITHOUT push handles, Push handles will be cut down on the vertical section. (NC 2219, 2227, 2228, 2245)	○
BACKREST UPHOLSTERY / CUSHION*			
ATB2171		Tension Adjustable Upholstery with Nylon overlay. Short = 480 mm / 7 straps. Long = 540 mm / 8 straps. Measured from seat plate to top of upholstery. (NC 2155, 1496, 2230, 2247, 2324, 2325)	◆
ATB2179		Flex3 Tension Adjustable (540 mm or 480 mm) (NC 2155, 2230, 2247, 2324, 2325)	£155
ATB2172		Flex3 adjustable Backrest Plates (without Backrest Cushion) (NC 2128)	£250
ATB2173		Flex3 adjustable Backrest Plates with Laguna Backrest cushion (NC 2128)	£400
ATB2174		Matrx Posture Back, height 400 mm (Mesh fabric) Corresponds to SW setting chosen	£353
ATB2175		Matrx Posture Back, height 500 mm (Mesh fabric) Corresponds to SW setting chosen	£353
ATB2176		Matrx Elite E2 Back , height 500 mm (Mesh fabric) Corresponds to SW setting chosen (NC 2121)	£494
ATBS177		Matrx Elite TR Back, height 500 mm (Mesh fabric) Corresponds to SW setting chosen (NC 2121)	£558
ATB2178		Matrx Elite E2 Deep Back, height 500 mm (Mesh fabric) Corresponds to SW setting chosen	£546
ATB2168		Matrx MAC Back - Foam Cushion [Mesh fabric] (NC 1445)	£475
ATB2169		Matrx MAC Back - Hug Stamina Fibre Cushion [Mesh fabric] (NC 1445)	£575
ATB2170		Without Backrest Upholstery (Back Posts/Back Frame only) (NC 2151, 2152, 2155, 2230, 2247, 2151, 2152, 2155)	○
MAC BACK BACKREST HEIGHT (Only applicable for MAC Back options)			
ATB2380		405 mm (16 inch) MAC BH (NC 1445)	○
ATB2381		455 mm (18 inch) MAC BH (NC 1445)	○
ATB2382		505 mm (20 inch) MAC BH (NC 1445)	○
ATB2383		560 mm (22 inch) MAC BH (NC 1445)	○
MAC BACK LATERAL CONTOUR DEPTH LH (Only applicable for MAC Back options)			
ATB2385		90 mm (3.5 inch) MAC Contour Depth LH (NC 1445)	£100
ATB2386		150 mm (6 inch) MAC Contour Depth LH (NC 1445)	£100
ATB2387		230 mm (9 inch) MAC Contour Depth LH (NC 1445)	£100
MAC BACK LATERAL CONTOUR DEPTH RH (Only applicable for MAC Back options)			
ATB2389		90 mm (3.5 inch) MAC Contour Depth RH (NC 1445)	£100
ATB2390		150 mm (6 inch) MAC Contour Depth RH (NC 1445)	£100
ATB2391		230 mm (9 inch) MAC Contour Depth RH (NC 1445)	£100
MAC BACK BACKREST WIDTH (Only applicable for MAC Back options)			
ATB2393		Backrest -25 mm (-1 inch) to seat width (NC 1445)	○
ATB2394		Backrest matches seat width (NC 1445)	○
ATB2395		Backrest +25 mm (+1 inch) to seat width (NC 1445)	○
UPHOLSTERY COLOUR*			
ATB2151		Black Nylon Fabric TR18 (NC 2170)	◆
ATB2152		Black Multi Stretch PU Coated Wipe Down Fabric TR26 (NC 2170)	○
ATB2155		Black Spacer Fabric TR36 (NC 2171, 2179, 2143, 2147)	○
ATX3002		No Upholstery Colour required - No backrest upholstery required (NC 2179, 2173, 2174, 2175, 2176, 2178)	○

ARMRESTS*			
ATB2191		T-Arm detachable, Standard. Adjustment range from seat plate: 245 - 315 mm	◆
ATB2181		T-Arm detachable, Extended height. Adjustment range from seat plate: 295 - 365 mm	£40
ATB2183		Inrigged Armrests. Reduces SW by 25 mm each side	○
ATB2193		Flip-up Armrest. Height and depth adjustable (height adjustable from 230mm - 360mm).	£65
ATB2192		Following Armrest (height adjustable from 230mm - 300mm).	£185
ATB2194		Following Armrest (height adjustable from 300mm - 360mm).	£185
ARMPADS*			
ATB2195		Short Armpad (325mm)	◆
ATB2196		Long Armpad (375mm)	○
ATB2197		Long Ergonomic Armpad (397mm)	(NC 2230, 2247) ○
ATB2199		Soft Armpad (REA Azalea Wide Pad)	(NC 2179) £95
HIP SUPPORTS LH			
ATB2184		Hip Support LH, pad size 140 x 100 mm	(NC 2702, 2703, 2704, 2705, 2706) £250
ATB2186		Hip Support LH, pad size 230 x 100 mm	(NC 2702, 2703, 2704, 2705, 2706) £250
ATB2188		Hip Support LH, pad size 330 x 100 mm	(NC 2702, 2703, 2704, 2705, 2706) £250
HIP SUPPORTS RH			
ATB2185		Hip Support RH, pad size 140 x 100 mm	(NC 2702, 2703, 2704, 2705, 2706) £250
ATB2187		Hip Support RH, pad size 230 x 100 mm	(NC 2702, 2703, 2704, 2705, 2706) £250
ATB2189		Hip Support RH, pad size 330 x 100 mm	(NC 2702, 2703, 2704, 2705, 2706) £250
MOUNTING HARDWARE FOR HIP SUPPORTS			
ATB2164		Mounting hardware for hip supports (6.1" / 155 mm high)	○
ATB2165		Mounting hardware for hip supports (8.1" / 206mm high)	○
HEADREST / NECKREST (Adequate headrest option must be selected)*			
ATB1350		Neckrest depth adjustable	(NC 2247) £80
ATB1432		Headrest sideways adjustable	(NC 2247) £225
ATB1427		Headrest with cheek support	(NC 2247) £225
ATB2324		Elan Standard Pad 255 mm / 10" wide w/standard cover (ESP26). Complete with Elan Headrest Hardware	(NC 2230, 2247) £325
ATB2325		Elan Occipital Pad 230x180 mm / 9"x7" w/standard cover (EOP23). Complete with Elan Headrest Hardware	(NC 2230, 2247) £325
ATB1496		Non Invacare Headrest but profile mount fitting Flex3 hard back	(NC 2171) ○
ATB1495		Non Invacare Headrest - no mounting hardware. With Backbrace/Pushbar dia. 22mm in conjunction with telescopic seat frame/backposts	○
LATERAL SUPPORTS			
ATB1445		Lateral Support with extended arm, multi-adjustable pad	(NC 2171) £165
SEAT HEIGHT AND SEAT RISER*			
ATB1023		Fixed Seat Height 450mm (17.75")	○
ATB1025		Fixed Seat Height 480mm (19")	○
ATB1030		Electrical Seat Riser (400mm - 650mm) (16" - 25.5")	£1,200
ATB1040		Electrical Seat Riser (440mm - 690mm) incl. adapter (17.25" - 27")	£1,275

SEAT TILT*			
ATB1002		Manual Seat Tilt by spindle	◆
ATB1014		Electrical Tilt 0° - 25°	£530
BACKREST RECLINER*			
ATB0975		Manual Settable Recliner (0° - 30°, in 3.75° steps)	◆
ATB0977		Powered Recliner (0° - 30°)	£500
LEGREST FIXATION POSITION*			
ATB0590		Position for seat cushion up to 70mm thick	◆
LEG REST TYPE*			
ATB0503		Side Mounted: Vari F pair. Lower leg length: 290 - 460 mm. Angle adjustment: 70 - 0° with tools	○
ATB0533		Side Mounted: Vari A pair. Lower leg length: 290 - 460 mm. (NC 0501) Angle adjustment: 70 - 0° with handwheel and calf pad	£235
ATB0531		Side Mounted: ADE pair. Lower leg length: 290 - 460 mm. (NC 0501) Powered angle adjustment: 80 - 0° with calf pad	£660
ATB0501		Without Leg Rests (NC 0504, 0513)	○
FOOTPLATE TYPE*			
ATB0512		Fixed Footplate (NC 0533, 0531)	○
ATB0504		Angle Adjustable Footplate	£75
ATB0513		Angle and depth adjustable Footplate (Adapt)	£120
ATB0501		Without footplates for option 0501 without Leg Rests (NC 0533, 0531)	○
LINX CONTROLS*			
ATB2201		DLX-REM110 (Drive Only) (NC 0977, 0531, 1437, 2209, 2210, 2230, 2247, 2257, 2233, 2240, 2270, 2219, 2245, 1908, 1909)	◆
ATB2203		DLX-REM211 (Drive & Actuators) (NC 2209, 2210, 2230, 2247, 2257, 2233, 2240, 2270, 1908, 1909)	○
ATB2204		DLX-REM216 (Drive, Actuators & Lights) (NC 2210, 2230, 2247, 2257, 2233, 2240, 2270)	£200
ATB2205		DLX-REM400 (Touch Screen Interface. Drive, Act. & Lights. Configurable jack sockets) (NC 2210)	£425
ATB2206		DLX-REM500 (Touch Screen Interface Display Unit - requires a user input/alternative joystick). (NC 2228, 2270, 1709, 1710)	£765
REM400 TOGGLE OPTION			
ATB2270		Assignable Toggle Switch Module for REM400 (DLX-SW400TGL) (NC 2203, 2204, 2206)	£315
CONTROL BRACKETRY*			
ATB1207		Standard Bracket RH (NC 2206)	○
ATB1208		Standard Bracket LH (NC 2206)	○
ATB1707		Swing Away RH Bracket (NC 2206)	£135
ATB1708		Swing Away LH Bracket (NC 2206)	£135
ATB1709		Matrx Resolve Swing in/out Bracket RH (NC 2213, 2235, 1404)	£245
ATB1710		Matrx Resolve Swing in/out Bracket LH (NC 2214, 2236, 1405)	£245
ATB2213		Centre Mount Swing Away to RH Nucleus Bracket (NC 1405)	£495
ATB2214		Centre Mount Swing Away to LH Nucleus Bracket (NC 1405)	£495
ATB2227		Attendant Mounted RH (NC 1404, 1405, 2128)	£65
ATB2228		Attendant Mounted LH (NC 1404, 1405, 2128)	£65
ATB2209		REM500 Swing Away RH Multi-Adjustable Bracket (NC 2203, 2204, 2205, 1404)	○
ATB2210		REM500 Swing Away LH Multi-Adjustable Bracket (NC 2203, 2204, 2205)	○
ATZ2239		Individual mounting option per request (Special). Requires additional information/specification.	POA

DUAL CONTROLS				
ATB2220		LiNX DLX-REM050 (Dual Control used as secondary remote)	(NC 1222, 2227, 2228, 2245, 2128, 2230, 2247, 2257, 2233, 2240)	£395
ATB2219		LiNX DLX-ACU200 (Compact Dual Control)	(NC 1222, 2227, 2228, 2245, 2128, 2201)	£395
ATB2245		LiNX DLX-IDC (Intuitive Dual Control)	(NC 2219, 2227, 2228, 2201)	£745
UNZ9259		mo-vis Scoot Control (Only available via Special)	(NC 2108, 0000, 1221, 1222)	£2,315
DUAL CONTROL BRACKETRY (Only applicable if Dual Control is selected)				
ATB1207		Dual Control Standard Bracket RH	(NC 1222, 2227, 2228, 2245)	○
ATB1208		Dual Control Standard Bracket LH	(NC 1222, 2227, 2228, 2245)	○
ATB1221		IDC Dual Control Detachable Bracket RH	(NC 2219, 2227, 2228)	○
ATB1222		IDC Dual Control Detachable Bracket LH	(NC 2219, 2227, 2228)	○
UNZ9260		Individual mounting option per request (Special). Requires additional information/specification.		POA
COMPACT REMOTES				
ATB2231		LiNX DLX-CR400 (Compact Remote)	(NC 2203, 2204, 2235, 2236, 2241, 2242, 2220)	£825
ATB2232		LiNX DLX-CR400LF (Compact Remote Low Force)	(NC 2203, 2204, 2235, 2236, 2241, 2242, 2220)	£1,075
UNZ9250		mo-vis Micro LiNX (Only available via Special)	(NC 2203, 2204, 0000, 2207, 2208, 2213, 2214, 2235, 2236, 2241, 2242, 2220)	£2,600
UNZ9251		mo-vis Micro Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204, 0000, 2207, 2208, 2213, 2214, 2235, 2236, 2241, 2242, 2220)	£2,730
UNZ9252		mo-vis Multi LiNX (Only available via Special)	(NC 2203, 2204, 0000, 2207, 2208, 2213, 2214, 2235, 2236, 2241, 2242, 2220)	£2,320
UNZ9253		mo-vis Multi Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204, 0000, 2207, 2208, 2213, 2214, 2235, 2236, 2241, 2242, 2220)	£2,520
UNZ9254		mo-vis All Round Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204, 0000, 2207, 2208, 2213, 2214, 2235, 2236, 2241, 2242, 2220)	£1,395
UNZ9255		mo-vis All Round Low Force Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204, 0000, 2207, 2208, 2213, 2214, 2235, 2236, 2241, 2242, 2220)	£2,015
UNZ9256		mo-vis Heavy Duty Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204, 0000, 2207, 2208, 2213, 2214, 2235, 2236, 2241, 2242, 2220)	£3,190
ATB2233		ASL MEC (Proportional MEC Joystick ASL130)	(NC 2203, 2204, 2220, 2230, 2247, 2257, 2240, 2128, 2207, 2208)	£4,100
ATZ2240		ASL Heavy Duty (Proportional Joystick w mode switch ASL133)	(NC 2203, 2204, 2220, 2230, 2247, 2257, 2233, 2128, 2207, 2208)	£2,285
COMPACT REMOTE BRACKETRY (Only applicable if Compact Remote is selected)				
ATB2207		RH Mount	(NC 2203, 2204, 2217, 1404)	○
ATB2208		LH Mount	(NC 2203, 2204, 2218, 1405)	○
ATB2213		Centre Mount Swing Away to RH (Nucleus)	(NC 2203, 2204, 1404, 1405)	£495
ATB2214		Centre Mount Swing Away to LH (Nucleus)	(NC 2203, 2204, 1404, 1405)	£495
ATB2235		Tray Mount RH (ASL603, 100mm x 150mm)	(NC 1415)	○
ATB2236		Tray Mount LH (ASL603, 100mm x 150mm)	(NC 1339)	○
ATB2241		Centre Mount, Swing Away to RH, Nucleus (Inc Tray ASL603, 100mm x 150mm)	(NC 1405, 1415)	£495
ATB2242		Centre Mount, Swing Away to LH, Nucleus (Inc Tray ASL603, 100mm x 150mm)	(NC 1405, 1339)	£495
ATZ2215		Individual mounting option per request (Special). Requires additional information/specification.	(NC 2203, 2204)	POA
CHIN CONTROLS				
UNZ9250		mo-vis Micro LiNX (Only available via Special)	(NC 2201, 2203, 2204, 2220)	£2,600
UNZ9251		mo-vis Micro Omni (Only available via Special. DLX-IN500 module needed)	(NC 2201, 2203, 2204, 2220)	£2,730
UNZ9252		mo-vis Multi LiNX (Only available via Special)	(NC 2201, 2203, 2204, 2220)	£2,320
UNZ9253		mo-vis Multi Omni (Only available via Special. DLX-IN500 module needed)	(NC 2201, 2203, 2204, 2220)	£2,520
UNZ9254		mo-vis All Round Omni (Only available via Special. DLX-IN500 module needed)	(NC 2201, 2203, 2204, 2220)	£1,395
UNZ9255		mo-vis All Round Low Force Omni (Only available via Special. DLX-IN500 module needed)	(NC 2201, 2203, 2204, 2220)	£2,015

CHIN CONTROL BRACKETRY (Only applicable if Chin Control is selected)			
UNZ9261		mo-vis Chin Control Harness (One Size. Special Request Only)	£960
UNZ9262		mo-vis Manual Swing Away (Special Mounting Only. Please specify what side mounting required)	£1,380
UNZ9263		mo-vis Powered Swing Away (Special Mounting Only. Please specify what side mounting required)	£2,280
HEAD CONTROLS			
ATB2230		ASL104 ATOM Head Control with ultrasonic proximity sensors with standard wings	(NC 2179, 2197, 2201, 2203, 2204, 2220, 2233, 2324, 2325, 1350, 1432, 1427, 2257, 2240) £4,220
ATB2247		ASL104P PROTON Head Control with ultrasonic proximity sensors with PROTON wings	(NC 2179, 2197, 2201, 2203, 2204, 2220, 2233, 2324, 2325, 1350, 1432, 1427, 2257, 2240) £4,945
UNZ9257		mo-vis Proportional Head Control Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204, 2233, 2324, 2325, 2257, 2240, 2175, 2171, 2170, 2197, 1350, 1432, 1427, 0000, 2220) £4,775
UNZ9258		mo-vis Switched Head Control Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204, 2233, 2324, 2325, 2257, 2240, 2175, 2171, 2170, 2197, 1350, 1432, 1427, 0000, 2220) £3,315
FOOT CONTROL			
UNZ9264		mo-vis Foot Control LiNX (Only available via Special)	(NC 2201, 2203, 2204, 2220) £2,880
UNZ9265		mo-vis Foot Control Omni (Only available via Special. DLX-IN500 module needed)	(NC 2201, 2203, 2204, 2220) £2,740
FOOT CONTROL BRACKETRY (Only applicable if Foot Control is selected)			
UNZ9266		Individual mounting option per request (Special). Requires additional information/specification.	POA
TRAY CONTROL			
ATZ2257		Four-Switch Proximity Array (ASL106)	(NC 2203, 2204, 2220, 2230, 2247, 2233, 2240, 2128) £6,535
TRAY CONTROL BRACKETRY (Only applicable if Tray Control is selected)			
ATZ2275		Centre Mount, swing-away to RH, nucleus (includes Eclipse Tray ASL619, 360mm x 180mm)	(NC 1405, 1415) £495
ATZ2276		Centre Mount, swing-away to LH, nucleus (includes Eclipse Tray ASL619, 360mm x 180mm)	(NC 1405, 1339) £495
ATZ2215		Individual mounting option per request (Special). Requires additional information/specification.	POA
WIRELESS ACCESSORIES FOR ASL CONTROLS			
ATB2249		Wireless Triple Switch Receiver (ASL557-3)	£340
ATB2250		Wireless Mouse Emulator (ASL558)	£810
LINX 10-WAY SWITCH			
ATB2251		10-Way Switch DLX-SW10-A for direct Powered Seating Access	£500
LINX 10-WAY SWITCH BRACKETRY (Only applicable if 10-Way Switch is selected)			
ATB2217		Bracket RH, mounted for user	(NC 2213, 2275, 2235, 2241) ○
ATB2218		Bracket LH, mounted for user	(NC 2214, 2276, 2236, 2242) ○
ATB2227		Bracket RH, mounted for attendant	£65
ATB2228		Bracket LH, mounted for attendant	£65
LINX ACCESSORY CABLE			
ATB2252		Loom providing jack sockets for Drive and Mode (includes connection for 10-way switch)	£100
LINX USB CHARGER PORT			
ATB2253		USB Charger Port DLX-USB02-A	(NC 2702, 2703, 2704, 2705, 2706) £160
LINX USB CHARGER PORT BRACKETRY (Only applicable if USB Charger is selected)			
ATB2207		RH Mount	○
ATB2208		LH Mount	○
FURTHER LINX COMPONENTS / ACCESSORIES			
ATB2302		Not Necessary	(NC 1686, 1687, 1690, 2254, 1711, 1712, 1713, 1714, 1715, 1716) ◆
ATB2300		To be enclosed only	○
ATZ2301		To be connected and programmed (Special)	POA

BUDDY BUTTONS			
ATB1685		Red Buddy Button	£90
ATB1686		Yellow Buddy Button	£90
ATB1687		Green Buddy Button	£90
ATB1690		Blue Buddy Button	£90
EGG SWITCH			
ATB2254		Egg Switch, black	£90
PIKO BUTTONS			
ATB1711		Red Piko Button	£160
ATB1712		Yellow Piko Button	£160
ATB1713		Green Piko Button	£160
ATB1714		Orange Piko Button	£160
ATB1715		Blue Piko Button	£160
ATB1716		Black Piko Button	£160
TOGGLES / BUTTONS			
ATZ2261		Stereo Toggle - direct control for Tilt only	£270
ATZ2262		Stereo Button - direct control for Tilt only	£255
ATZ2263		Four-Way Toggle - direct control for Tilt & Recline	£330
ATZ2264		Four-Way Button - direct control for Tilt & Recline	£355
ATZ2207		Toggle / Button RH Mount	○
ATZ2208		Toggle / Button LH Mount	○
ATZ2225		Without Mounting	○
BODYPOINT JOYSTICK KNOBS			
ATB2721		3" U-Shaped, Joystick Handle, Bodypoint PC101A	£75
ATB2722		4" U-Shaped, Joystick Handle, Bodypoint PC102A	£75
ATB2723		Domed Rubber, Joystick Handle, Bodypoint PC107A	£75
ATB2724		Mushroom Grooved, Joystick Handle, Bodypoint PC110A	£75
SHROUD COLOUR*			
ATB0202		Blue	○
ATB0228		Black	○
ATB0204		Transparent Cherry Red	○
ATZ0299		Special shroud colour (RAL code required)	POA
CASTORS*			
ATB0851		3.00 - 6 Grey Pneumatic	○
ATB0853		3.00 - 6 Grey Puncture Proof (sealed)	○
ATB0855		3.00 - 6 Grey Puncture Proof (solid inlay)	£105
DRIVE WHEELS*			
ATB0771		3.00 - 8 Grey Pneumatic	○
ATB0773		3.00 - 8 Grey Puncture Proof (sealed)	○
ATB0750		3.00 - 8 Grey Puncture Proof (solid inlay)	£105
MOTOR TYPE*			
ATB1277		6 km/h (4mph)	○
ATB1908		10 km/h (6mph) (includes LED lights system with indicators)	£515
ATB1909		12 km/h (8mph) (includes LED light system with indicators)	£600

BATTERY*				
ATB1364	<input type="checkbox"/>	73.5 Ah Gel		◆
SUSPENSION*				
ATB1240	<input type="checkbox"/>	Suspension low (up to 80kg user weight)		○
ATB1242	<input type="checkbox"/>	Suspension hard (from 90kg up to 150kg user weight)		○
G TRAC				
ATB1215	<input type="checkbox"/>	G-Trac Module		£400
CHARGER*				
ATB1343	<input type="checkbox"/>	OptiCharge Charger 8A (200 - 250V)		◆
ATB1313	<input type="checkbox"/>	OptiCharge Charger 10A (200 - 250V) Charging time will be 20% faster		£155
LIGHTS*				
ATB1437	<input type="checkbox"/>	LED Light system with indicators (30 LUX)	(NC 1909)	£450
ATB1254	<input type="checkbox"/>	Without Light system		◆
ATX1006	<input type="checkbox"/>	Lights Included with Motor Package 10 & 12kmh		○
ANTI TIPPER*				
ATB1205	<input type="checkbox"/>	Anti-tipper (standard version)		○
KERB RIDER				
ATB1244	<input type="checkbox"/>	Kerb rider		£215
TRANSPORTATION KIT*				
ATB1450	<input type="checkbox"/>	Preparation for DAHL-Docking station (Adequate headrest is compulsory). Docking Station for car implementation has to be ordered separately with DAHL Engineering.		£580
ATB1494	<input type="checkbox"/>	4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19). Adequate headrest is compulsory.		◆
FIXATION BELT*				
ATB1329	<input type="checkbox"/>	Lap Belt with Auto Lock	(NC 2702, 2703, 2704, 2705, 2706)	○
ATB2113	<input type="checkbox"/>	Lap belt with auto lock (flexible mount in side rails)	(NC 2702, 2703, 2704, 2705, 2706)	◆
ATB1332	<input type="checkbox"/>	H-Belt		£75
ATB2701	<input type="checkbox"/>	Bodypoint Evoflex Pelvic, Stabilizer, Rehab Latch size L (padded length 620mm)	(NC 2121, 2122, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 1329, 2113)	£95
ATB2702	<input type="checkbox"/>	Bodypoint Evoflex Pelvic, Stabilizer, Rehab Latch size M (padded length 460mm)	(NC 2121, 2122, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 1329, 2113)	£95
ATB2703	<input type="checkbox"/>	Bodypoint 2-Point Hip Belt, Centre pull - push button, Cinch mount size L (padded length 620mm)	(NC 2121, 2122, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 1329, 2113)	£90
ATB2704	<input type="checkbox"/>	Bodypoint 2-Point Hip Belt, Centre pull - push button, Cinch mount size M (padded length 460mm)	(NC 2121, 2122, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 1329, 2113)	£90
ATB2705	<input type="checkbox"/>	Bodypoint 2-Point Hip Belt, Dual pull Side release button, Cinch mount size L (padded length 620mm)	(NC 2121, 2122, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 1329, 2113)	£90
ATB2706	<input type="checkbox"/>	Bodypoint 2-Point Hip Belt, Dual pull Side release button, Cinch mount size M (padded length 460mm)	(NC 2121, 2122, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 1329, 2113)	£90
REAR VIEW MIRROR				
ATB1415	<input type="checkbox"/>	Rear view Mirror RH - Must be selected with 10 km/h Motor	(NC 2213, 2275, 2235, 2241)	£45
ATB1339	<input type="checkbox"/>	Rear view Mirror LH - Must be selected with 10 km/h Motor	(NC 2214, 2276, 2236, 2242)	£45
BAG HOOKS				
ATB2128	<input type="checkbox"/>	Bag Hooks on Modulite Seat Max load up to 6 kg = impact on max user weight	(NC 2172, 2173, 2227, 2228, 2220, 2219, 2245, 2257, 2233, 2240)	£80
CANE HOLDER				
ATB1312	<input type="checkbox"/>	Cane holder RH		£45
ATB1363	<input type="checkbox"/>	Cane holder LH		£45

TRAY			
ATB1404		Tray Swing Away RH	(NC 1707, 2213, 2214, 2209, 2235, 2236, 2241, 2242, 2207, 2275, 2276) £200
ATB1405		Tray Swing Away LH	(NC 1708, 2213, 2214, 2210, 2235, 2236, 2241, 2242, 2208, 2275, 2276) £200
WHEEL FENDERS			
ATB1416		Fender for Castor Wheels	£110
ATB1417		Fender for Drive Wheels	£110

TOTAL PRICE	
--------------------	--

Non Conforming Options:

ATB0000	Storm 4 Modulite
ATB2121	380 mm Seat Plate. Can be reduced to 330 mm in combination with inrigged armrests.
ATB2122	430 mm Seat Plate. Can be reduced to 380 mm in combination with inrigged armrests.
ATB2132	410 mm
ATB2134	460 mm
ATB2142	Matrx Flo-tech Lite Cushion
ATB2143	Matrx Flo-tech Contour Visco; (2-way stretch material)
ATB2147	Matrx Libra Cushion, 90 mm thick (Extended lead time)
ATB2144	Matrx Vi Cushion (Mesh fabric)
ATB1980	Tapered Wedges (Pair)
ATB1981	Abductor
ATB1982	Lateral Wedge
ATB1983	Anterior Seat Wedges
ATB1980	Tapered Wedges (Pair)
ATB1981	Abductor
ATB1982	Lateral Wedge
ATB1983	Anterior Seat Wedges
ATB1985	Pelvic Obliquity Kit (Contains 2x positioning accessories)
ATB2108	Backposts WITHOUT push handles, Push handles will be cut down on the vertical section.
ATB2171	Tension Adjustable Upholstery with Nylon overlay. Short = 480 mm / 7 straps. Long = 540 mm / 8 straps. Measured from seat plate to top of upholstery.
ATB2179	Flex3 Tension Adjustable (540 mm or 480 mm)
ATB2172	Flex3 adjustable Backrest Plates (without Backrest Cushion)
ATB2173	Flex3 adjustable Backrest Plates with Laguna Backrest cushion
ATB2174	Matrx Posture Back, height 400 mm (Mesh fabric) Corresponds to SW setting chosen
ATB2175	Matrx Posture Back, height 500 mm (Mesh fabric) Corresponds to SW setting chosen
ATB2176	Matrx Elite E2 Back , height 500 mm (Mesh fabric) Corresponds to SW setting chosen
ATB2178	Matrx Elite E2 Deep Back, height 500 mm (Mesh fabric) Corresponds to SW setting chosen
ATB2168	Matrx MAC Back - Foam Cushion [Mesh fabric]
ATB2169	Matrx MAC Back - Hug Stamina Fibre Cushion [Mesh fabric]
ATB2170	Without Backrest Upholstery (Back Posts/Back Frame only)
ATB2380	405 mm (16 inch) MAC BH
ATB2381	455 mm (18 inch) MAC BH
ATB2382	505 mm (20 inch) MAC BH
ATB2383	560 mm (22 inch) MAC BH
ATB2385	90 mm (3.5 inch) MAC Contour Depth LH
ATB2386	150 mm (6 inch) MAC Contour Depth LH
ATB2387	230 mm (9 inch) MAC Contour Depth LH
ATB2389	90 mm (3.5 inch) MAC Contour Depth RH
ATB2390	150 mm (6 inch) MAC Contour Depth RH

ATB2391	230 mm (9 inch) MAC Contour Depth RH
ATB2393	Backrest -25 mm (-1 inch) to seat width
ATB2394	Backrest matches seat width
ATB2395	Backrest +25 mm (+1 inch) to seat width
ATB2151	Black Nylon Fabric TR18
ATB2152	Black Multi Stretch PU Coated Wipe Down Fabric TR26
ATB2155	Black Spacer Fabric TR36
ATB2183	Inrigged Armrests. Reduces SW by 25 mm each side
ATB2197	Long Ergonomic Armpad (397mm)
ATB2184	Hip Support LH, pad size 140 x 100 mm
ATB2186	Hip Support LH, pad size 230 x 100 mm
ATB2188	Hip Support LH, pad size 330 x 100 mm
ATB2185	Hip Support RH, pad size 140 x 100 mm
ATB2187	Hip Support RH, pad size 230 x 100 mm
ATB2189	Hip Support RH, pad size 330 x 100 mm
ATB1350	Neckrest depth adjustable
ATB1432	Headrest sideways adjustable
ATB1427	Headrest with cheek support
ATB2324	Elan Standard Pad 255 mm / 10" wide w/standard cover (ESP26). Complete with Elan Headrest Hardware
ATB2325	Elan Occipital Pad 230x180 mm / 9"x7" w/standard cover (EOP23). Complete with Elan Headrest Hardware
ATB1496	Non Invacare Headrest but profile mount fitting Flex3 hard back
ATB1445	Lateral Support with extended arm, multi-adjustable pad
ATB0977	Powered Recliner (0° - 30°)
ATB0533	Side Mounted: Vari A pair. Lower leg length: 290 - 460 mm. Angle adjustment: 70 - 0° with handwheel and calf pad
ATB0531	Side Mounted: ADE pair. Lower leg length: 290 - 460 mm. Powered angle adjustment: 80 - 0° with calf pad
ATB0501	Without Leg Rests
ATB0504	Angle Adjustable Footplate
ATB0513	Angle and depth adjustable Footplate (Adapt)
ATB0501	Without footplates for option 0501 without Leg Rests
ATB2201	DLX-REM110 (Drive Only)
ATB2203	DLX-REM211 (Drive & Actuators)
ATB2204	DLX-REM216 (Drive, Actuators & Lights)
ATB2205	DLX-REM400 (Touch Screen Interface. Drive, Act. & Lights. Configurable jack sockets)
ATB2206	DLX-REM500 (Touch Screen Interface Display Unit - requires a user input/alternative joystick).
ATB2270	Assignable Toggle Switch Module for REM400 (DLX-SW400TGL)
ATB1707	Swing Away RH Bracket
ATB1708	Swing Away LH Bracket
ATB1709	Matrx Resolve Swing in/out Bracket RH
ATB1710	Matrx Resolve Swing in/out Bracket LH
ATB2213	Centre Mount Swing Away to RH Nucleus Bracket
ATB2214	Centre Mount Swing Away to LH Nucleus Bracket
ATB2227	Attendant Mounted RH
ATB2228	Attendant Mounted LH

ATB2209	REM500 Swing Away RH Multi-Adjustable Bracket
ATB2210	REM500 Swing Away LH Multi-Adjustable Bracket
ATB2220	LiNX DLX-REM050 (Dual Control used as secondary remote)
ATB2219	LiNX DLX-ACU200 (Compact Dual Control)
ATB2245	LiNX DLX-IDC (Intuitive Dual Control)
ATB1221	IDC Dual Control Detachable Bracket RH
ATB1222	IDC Dual Control Detachable Bracket LH
ATB2233	ASL MEC (Proportional MEC Joystick ASL130)
ATZ2240	ASL Heavy Duty (Proportional Joystick w mode switch ASL133)
ATB2207	RH Mount
ATB2208	LH Mount
ATB2213	Centre Mount Swing Away to RH (Nucleus)
ATB2214	Centre Mount Swing Away to LH (Nucleus)
ATB2235	Tray Mount RH (ASL603, 100mm x 150mm)
ATB2236	Tray Mount LH (ASL603, 100mm x 150mm)
ATB2241	Centre Mount, Swing Away to RH, Nucleus (Inc Tray ASL603, 100mm x 150mm)
ATB2242	Centre Mount, Swing Away to LH, Nucleus (Inc Tray ASL603, 100mm x 150mm)
ATB2230	ASL104 ATOM Head Control with ultrasonic proximity sensors with standard wings
ATB2247	ASL104P PROTON Head Control with ultrasonic proximity sensors with PROTON wings
ATZ2257	Four-Switch Proximity Array (ASL106)
ATZ2275	Centre Mount, swing-away to RH, nucleus (includes Eclipse Tray ASL619, 360mm x 180mm)
ATZ2276	Centre Mount, swing-away to LH, nucleus (includes Eclipse Tray ASL619, 360mm x 180mm)
ATB2217	Bracket RH, mounted for user
ATB2218	Bracket LH, mounted for user
ATB2227	Bracket RH, mounted for attendant
ATB2228	Bracket LH, mounted for attendant
ATB2253	USB Charger Port DLX-USB02-A
ATB2207	RH Mount
ATB2208	LH Mount
ATB1686	Yellow Buddy Button
ATB1687	Green Buddy Button
ATB1690	Blue Buddy Button
ATB2254	Egg Switch, black
ATB1711	Red Piko Button
ATB1712	Yellow Piko Button
ATB1713	Green Piko Button
ATB1714	Orange Piko Button
ATB1715	Blue Piko Button
ATB1716	Black Piko Button
ATZ2207	Toggle / Button RH Mount
ATZ2208	Toggle / Button LH Mount
ATB1908	10 km/h (6mph) (includes LED lights system with indicators)
ATB1909	12 km/h (8mph) (includes LED light system with indicators)
ATB1437	LED Light system with indicators (30 LUX)

ATB1329	Lap Belt with Auto Lock
ATB2113	Lap belt with auto lock (flexible mount in side rails)
ATB2701	Bodypoint Evoflex Pelvic, Stabilizer, Rehab Latch size L (padded length 620mm)
ATB2702	Bodypoint Evoflex Pelvic, Stabilizer, Rehab Latch size M (padded length 460mm)
ATB2703	Bodypoint 2-Point Hip Belt, Centre pull - push button, Cinch mount size L (padded length 620mm)
ATB2704	Bodypoint 2-Point Hip Belt, Centre pull - push button, Cinch mount size M (padded length 460mm)
ATB2705	Bodypoint 2-Point Hip Belt, Dual pull Side release button, Cinch mount size L (padded length 620mm)
ATB2706	Bodypoint 2-Point Hip Belt, Dual pull Side release button, Cinch mount size M (padded length 460mm)
ATB1415	Rear view Mirror RH - Must be selected with 10 km/h Motor
ATB1339	Rear view Mirror LH - Must be selected with 10 km/h Motor
ATB2128	Bag Hooks on Modulite Seat Max load up to 6 kg = impact on max user weight
ATB1404	Tray Swing Away RH
ATB1405	Tray Swing Away LH



For a more complex customer requirement,
please email: orders@invacare.com.au

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):