



Yes, you can:

Tel: 01656 776222
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

Power Retail Prescription Form

July 2026
Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

TDX SP2 HD Modulite

Standard Features

- Max User Weight up to 180kg (28 stone)
- Crash tested ISO7176 - 19
- Highly adjustable Modulite Seating System
- Telescopic adjustable seat widths from 330mm - 530mm (13" - 21") with an additional 50mm, (2") by arms
- Seat depths 410mm - 560mm stepless adjustment (16" - 22")
- Tension Adjustable backrest
- Standard T-Armrests, height and depth adjustable, detachable, with short armpads
- Manual settable Seat Tilt (0° - 20°)
- Manual settable Backrest Recline (0° - 30° in 3.75° steps)
- Standard 80° hangers with fixed footplate
- LiNX DLX-REM110 with standard bracket
- Selection of 5 shroud and rim colours
- Centre (Driving Wheels) 355mm x 76mm (14" x 3") Black Tyres
- Suspension and Sure Step Stability Technology
- 6 km/h (4mph) 4 Pole Motors
- 8A (200 - 250V) Charger
- 73.5 Ah Batteries
- 4 Point tie down system
- Lapbelt with autolock (adjustable both sides)

Standard Base price: £7,995



Please make your selection with a cross in the box ☐ Free of charge ☒ standard feature POA special request

Customer weight = KG

Code	Description	
BASE MODEL*		
ARJ0000	TDX SP2 Heavy Duty	£7,995
ARZ1300	Special TDX SP2 Heavy Duty (Extended Lead Time). Please provide details in comments section.	POA

USER WEIGHT*			
ARJ1918		User weight less 100 kg	○
ARJ1919		User weight 101 to 125 kg	○
ARJ1921		User weight 126 to 136 kg	○
ARJ1922		User weight 137 to 150 kg (NC 1092, 2052, 2053)	○
ARJ1923		User weight 151 to 160 kg (NC 2192, 2193, 2194, 1092, 2052, 2053)	○
ARJ1924		User weight 161 to 180 kg (NC 2192, 2193, 2194, 1092, 1093, 1094, 1909, 2052, 2053)	○
SEATING TYPE*			
ARJ2100		Modulite Seat Unit	◆
SEAT FRAME*			
ARJ2111		Telescopic Seat Frame + Plate with Side Rails. Width 380 mm to 530 mm. Additional 50 mm width by arm rest adjustment.	◆
SEAT WIDTH*			
ARZ1001		Seat width 330 mm. Seat width achieved on a 380 mm seat frame (2121) with inrigged armrests (2183)	○
ARJ2121		380 mm Seat Plate. Can be reduced to 330 mm in combination with inrigged armrests. (NC 2143, 2701, 2703, 2705)	○
ARJ2122		430 mm Seat Plate. Can be reduced to 380 mm in combination with inrigged armrests. (NC 2703, 2705)	○
ARJ2123		480 mm Seat Plate (NC 2702, 2704, 2706)	○
ARJ2124		530 mm Seat Plate (NC 2174, 2175, 2176, 2177, 2178, 2143, 2702, 2704, 2706)	○
BODY FIT KIT			
UNZ2106		Body Fit Kit - Outrigging Backrest Brackets and Pivot Point (outrig: 1" /25mm per side) (NC 2124, 2192, 2194)	○
SEAT DEPTH RAIL LENGTH*			
ARJ2129		Standard Depth Package (410mm - 510mm. Simple stepless adjustment from the side) (NC 2138)	○
ARJ2130		Extended Depth Package (510mm - 560mm. Simple stepless adjustment from the side) (NC 2132, 2134)	○
SEAT DEPTH*			
ARJ2132		410 mm (NC 2196, 2197, 2701, 2702, 2703, 2704, 2705, 2706)	○
ARJ2134		460 mm	○
ARJ2136		510 mm	○
ARJ2138		560 mm (NC 2143)	○
SEAT OPTION			
ARZSWO		"H" Frame Interface Kit	£495
SEAT CUSHION*			
ARJ2142		Matrx Flo-tech Cushion (Lite - 75mm up to 127kg) or (Contour - 100mm up to 191kg) (NC 1981, 1982, 1983, 1980, 1981, 1982, 1983, 1985)	£120
ARJ2143		Matrx Flo-tech Contour Visco; (2-way stretch material, 430 mm x 430 mm, up to 125kg) (NC 1981, 1982, 1983, 1980, 1981, 1982, 1983, 1985)	£169
ARJ2147		Matrx Libra (90mm - standard up to 136kg or HD up to 227kg). Extended lead time.	£410
ARJ2144		Matrx Vi Cushion (Mesh fabric) (NC 1981, 1982, 1983, 1980, 1981, 1982, 1983)	£400
ARJ2140		Without seat cushion (NC 1981, 1982, 1983, 1980, 1981, 1982, 1983, 1985)	○
1st POSITIONING ACCESORIES FOR MATRX LIBRA			
ARJ1980		Tapered Wedges (Pair)	£30
ARJ1981		Abductor	£23
ARJ1982		Lateral Wedge	£20
ARJ1983		Anterior Seat Wedges	£27

2nd POSITIONING ACCESORIES FOR MATRX LIBRA			
ARJ1980		Tapered Wedges (Pair)	£30
ARJ1981		Abductor	£23
ARJ1982		Lateral Wedge	£20
ARJ1983		Anterior Seat Wedges	£27
POSITIONING ACCESSORIES FOR MATRX VI			
ARJ1985		Pelvic Obliquity Kit (Contains 2x positioning accessories)	£54
BACKREST SUPPORT*			
ARJ2163		Backposts short (480 mm). For Tension Adjustable Upholstery/Matrx Backs, measured from seat plate to top of upholstery. Width adjustable cast bracket.	(NC 2173, 2052, 2053) ○
ARJ2162		Backposts tall (540 mm) For Flex3 backrests/Tension Adjustable Upholstery/Matrx Backs, measured from seat plate to top of upholstery. Width adjustable cast bracket.	◆
PUSH HANDLES FOR BACKPOSTS*			
ARJ2107		Backposts WITH push handles	◆
ARJ2108		Backposts WITHOUT push handles, Push handles will be cut down on the vertical section.	(NC 1332, 2220, 2219, 2245, 2227, 2228) ○
BACKREST UPHOLSTERY / CUSHION*			
ARJ2171		Tension Adjustable Upholstery with Nylon overlay. Short = 480 mm / 7 straps. Long = 540 mm / 8 straps. Measured from seat plate to top of upholstery.	(NC 1445, 1496, 2155, 2230, 2247, 2052, 2053) ◆
ARJ2179		Flex3 Tension Adjustable (540 mm or 480 mm)	(NC 2151, 2192, 2193, 2194, 2199, 2155, 2230, 2247, 2128) £155
ARJ2172		Flex3 adjustable Backrest Plates (without Backrest Cushion)	(NC 2163, 2128) £250
ARJ2173		Flex3 adjustable Backrest Plates with Laguna Backrest cushion	(NC 2163, 2128) £400
ARJ2174		Matrx Posture Back, height 400 mm (Mesh fabric) Corresponds to SW setting chosen	(NC 1496, 2124) £353
ARJ2175		Matrx Posture Back, height 500 mm (Mesh fabric) Corresponds to SW setting chosen	(NC 1496, 2124) £353
ARJ2176		Matrx Elite E2 Back , height 500 mm (Mesh fabric) Corresponds to SW setting chosen	(NC 1496, 2124) £494
ARJ2177		Matrx Elite TR Back, height 500 mm (Mesh fabric) Corresponds to SW setting chosen	(NC 1496, 2124) £496
ARJ2178		Matrx Elite E2 Deep Back, height 500 mm (Mesh fabric) Corresponds to SW setting chosen	(NC 1496, 2124) £546
ARJ2168		Matrx MAC Back - Foam Cushion [Mesh fabric]	(NC 1445) £475
ARJ2169		Matrx MAC Back - Hug Stamina Fibre Cushion [Mesh fabric]	(NC 1445) £575
ARJ2170		Without Backrest Upholstery (Back Posts/Back Frame only)	(NC 1496, 2230, 2247, 2151, 2152, 2155) ○
MAC BACK BACKREST HEIGHT (Only applicable for MAC Back options)			
ARJ2380		405 mm (16 inch) MAC BH	(NC 1445) ○
ARJ2381		455 mm (18 inch) MAC BH	(NC 1445) ○
ARJ2382		505 mm (20 inch) MAC BH	(NC 1445) ○
ARJ2383		560 mm (22 inch) MAC BH	(NC 1445) ○
MAC BACK LATERAL CONTOUR DEPTH LH (Only applicable for MAC Back options)			
ARJ2385		90 mm (3.5 inch) MAC Contour Depth LH	(NC 1445) £100
ARJ2386		150 mm (6 inch) MAC Contour Depth LH	(NC 1445) £100
ARJ2387		230 mm (9 inch) MAC Contour Depth LH	(NC 1445) £100
MAC BACK LATERAL CONTOUR DEPTH RH (Only applicable for MAC Back options)			
ARJ2389		90 mm (3.5 inch) MAC Contour Depth RH	(NC 1445) £100
ARJ2390		150 mm (6 inch) MAC Contour Depth RH	(NC 1445) £100
ARJ2391		230 mm (9 inch) MAC Contour Depth RH	(NC 1445) £100

MAC BACK BACKREST WIDTH (Only applicable for MAC Back options)				
ARJ2393		Backrest -25 mm (-1 inch) to seat width	(NC 1445)	○
ARJ2394		Backrest matches seat width	(NC 1445)	○
ARJ2395		Backrest +25 mm (+1 inch) to seat width	(NC 1445)	○
UPHOLSTERY COLOUR*				
ARJ2151		Black Nylon Fabric TR18		◆
ARJ2152		Black Multi Stretch PU Coated Wipe Down Fabric TR26	(NC 2171)	○
ARJ2155		Black Spacer Fabric TR36	(NC 2171, 2179)	○
ARX1002		No Upholstery Colour required - No backrest upholstery required	(NC 2179, 2172, 2173, 2174, 2175, 2176, 2177, 2178)	○
ARMRESTS*				
ARJ2191		T-Arm detachable, Standard. Adjustment range from seat plate: 245 - 315 mm		◆
ARJ2181		T-Arm detachable, Extended height. Adjustment range from seat plate: 295 - 365 mm		£40
ARJ2183		Inrigged Armrests. Reduces SW by 25 mm each side	(NC 2185, 2186, 2187, 2188, 2189)	£30
ARJ2193		Flip-up Armrest. Height and depth adjustable (height adjustable from 230mm - 360mm). 150 kg max user weight.	(NC 2176, 2199)	£65
ARJ2192		Following Armrest (height adjustable from 230mm - 300mm). 150 kg max user weight	(NC 1445, 2179)	£185
ARJ2194		Following Armrest (height adjustable from 300mm - 360mm). 150 kg max user weight	(NC 1445, 2179)	£185
ARMPADS*				
ARJ2195		Short Armpad (325mm)		◆
ARJ2196		Long Armpad (375mm)		○
ARJ2197		Long Ergonomic Armpad (397mm)	(NC 2247)	○
ARJ2199		Soft Armpad (REA Azalea Wide Pad)		£95
HIP SUPPORTS LH				
ARJ2184		Hip Support LH, pad size 140 x 100 mm	(NC 2702, 2703, 2704, 2705, 2706)	£250
ARJ2186		Hip Support LH, pad size 230 x 100 mm	(NC 2702, 2703, 2704, 2705, 2706)	£250
ARJ2188		Hip Support LH, pad size 330 x 100 mm	(NC 2702, 2703, 2704, 2705, 2706)	£250
HIP SUPPORTS RH				
ARJ2185		Hip Support RH, pad size 140 x 100 mm	(NC 2702, 2703, 2704, 2705, 2706)	£250
ARJ2187		Hip Support RH, pad size 230 x 100 mm	(NC 2702, 2703, 2704, 2705, 2706)	£250
ARJ2189		Hip Support RH, pad size 330 x 100 mm	(NC 2702, 2703, 2704, 2705, 2706)	£250
MOUNTING HARDWARE FOR HIP SUPPORTS				
ARJ2164		Mounting hardware for hip supports (6.1" / 155 mm high)		○
ARJ2165		Mounting hardware for hip supports (8.1" / 206mm high)		○

HEADREST / NECKREST (Adequate headrest option must be selected)*				
ARJ1350		Neckrest depth adjustable	(NC 2247)	£80
ARJ1432		Headrest sideways adjustable	(NC 2247)	◆
ARJ1427		Headrest with cheek support	(NC 2247)	£225
ARJ2324		Elan Standard Pad 255 mm / 10" wide w/standard cover (ESP26). Complete with Elan Headrest Hardware	(NC 2247)	£325
ARJ2325		Elan Occipital Pad 230x180 mm / 9"x7" w/standard cover (EOP23). Complete with Elan Headrest Hardware	(NC 2247)	£325
ARJ2327		Elan 4 Point Headrest Pad 280x255 mm / 11"x10" w/standard cover. Complete with Elan Headrest Hardware	(NC 2247)	£440
ARJ1496		Non Invacare Headrest but profile mount fitting Flex3 hard back		○
ARJ1495		Non Invacare Headrest - no mounting hardware. With Backbrace/Pushbar dia. 22mm in conjunction with telescopic seat frame/backposts		○
LATERAL SUPPORT				
ARJ1445		Lateral Support with extended arm, multi-adjustable pad	(NC 2171)	£165
MANUAL / POWERED SEATING OPTIONS*				
ARJ1091		Manually settable Seat Tilt 0° - 20°		◆
ARJ1092		Powered "fixed pivot" Tilt 20° (136kg max user weight)		£400
ARJ1095		Powered "fixed pivot" Tilt 20°, heavy duty (180kg max user weight)		£400
ARJ1093		Powered Tilt (0° - 30°) with CoG shift (160kg max user weight)	(NC 2201)	£1,100
ARJ2052		Powered Tilt (0°-50°) with CoG shift (136kg max user weight)	(NC 2201)	£1,450
ARJ1094		Seat Lifter with Powered Tilt (300mm, 30°) with CoG shift (160kg max user weight)	(NC 2201)	£1,900
ARJ2053		Seat Lifter with Powered Tilt (300mm, 45°) with CoG shift (136kg max user weight)	(NC 2201)	£2,500
SEAT TO GROUND HEIGHT*				
ARJ0354		450mm (17.75") (from floor to front of seat plate at 0° tilt)	(NC 1094, 2052, 2053)	○
ARJ0364		475mm (18.75") (from floor to front of seat plate at 0° tilt)		○
SEAT TILT OPTION				
ARJ2054		5° Fixed Anterior Tilt, Assist mount. Will reduce max tilt angle by 5°	(NC 1092, 1095)	£400
BACKREST RECLINER*				
ARJ0975		Manual Settable Recliner (0° - 30°, in 3.75° steps) Foldable		◆
ARJ0977		Powered Recliner (0° - 30°) Foldable		£500

LEG REST TYPE*				
ARJ0556		Centre Mounted: Manual Centre Mounted Legrest, angle settable. Length 280 - 385 mm	(NC 0517, 0512, 0504, 0513, 0501)	£370
ARJ0589		Centre Mounted: LNX Powered Centre Mounted Legrest, length compensating. Length 340 - 410 mm	(NC 2201, 0504, 0512, 0504, 0513, 0501)	£1,620
UNZ0592		Centre Mounted: LNX Powered Centre Mounted Legrest with Telescoping Foot Platform (lowers to the ground). Length: 380-419mm (Selectable length depends on STF height) Custom-made, requires customised setting.	(NC 2201, 0504, 0512, 0504, 0513, 0501)	POA
ARJ0526		Side Mounted: Standard hanger 80° pair. Lower leg length: 350 - 450 mm	(NC 0501)	◆
ARJ0503		Side Mounted: Vari F pair. Lower leg length: 290 - 460 mm. Angle adjustment: 70 - 0° with tools	(NC 0555, 0501)	£120
ARJ0533		Side Mounted: Vari A pair. Lower leg length: 290 - 460 mm. Angle adjustment: 70 - 0° with handwheel and calf pad	(NC 0554, 0555, 0501)	£235
ARJ0531		Side Mounted: ADE pair. Lower leg length: 290 - 460 mm. Powered angle adjustment: 80 - 0° with calf pad	(NC 2201, 0554, 0555, 0501)	£660
ARJ0501		Without Leg Rests	(NC 0504, 0513, 2365, 0559, 0516, 0517, 0559, 0516, 0517)	○
FOOTPLATE TYPE*				
ARJ2365		Angle & width adjustable footplates S (127 W x 190 mm D / 5"x7.5")	(NC 0503, 0533, 0531)	£105
ARJ0559		Angle & width adjustable footplates M (140 W x 240 mm D / 5.5"x9.5")	(NC 0503, 0533, 0531)	○
ARJ0516		Angle & width adjustable footplates L (150 W x 290 mm D / 6"x11.5")	(NC 0503, 0533, 0531)	£135
ARJ0517		Angle adjustable foot platform M (305 W x 275 mm D / 12"x10.75")	(NC 0526, 0503, 0533, 0531)	£270
ARJ0512		Fixed Footplate pair	(NC 0533, 0531)	○
ARJ0504		Angle Adjustable Footplate pair	(NC 0526)	£75
ARJ0513		Angle and depth adjustable Footplate pair		£120
ARJ0501		Without footplates for option 0501 without Leg Rests	(NC 0503, 0533, 0531)	○
AMPUTEE SUPPORT				
ARJ0554		Amputee Support RH	(NC 0556, 0589, 0592)	£195
ARJ0555		Amputee Support LH	(NC 0556, 0589, 0592)	£195
LINX CONTROLS*				
ARJ2201		DLX-REM110 (Drive Only)	(NC 1095, 1093, 1094, 0977, 0589, 0531, 2209, 2210, 2213, 2214, 2230, 2247, 2233, 2240, 2270, 2219, 2245, 1909)	◆
ARJ2203		DLX-REM211 (Drive & Actuators)	(NC 2209, 2210, 2213, 2214, 2231, 2232, 2230, 2247, 2233, 2240, 2270, 1909)	○
ARJ2204		DLX-REM216 (Drive, Actuators & Lights)	(NC 2210, 2213, 2214, 2231, 2232, 2230, 2247, 2233, 2240, 2270)	£200
ARJ2205		DLX-REM400 (Touch Screen Interface. Drive, Act. & Lights. Configurable jack sockets)	(NC 2210, 2213, 2215)	£425
ARJ2206		DLX-REM500 (Touch Screen Interface Display Unit - requires a user input/alternative joystick).	(NC 1208, 1707, 1708, 1709, 1710, 2220, 2270)	£765
REM400 TOGGLE OPTION				
ARJ2270		Assignable Toggle Switch Module for REM400 (DLX-SW400TGL)	(NC 2203, 2204)	£315

CONTROL BRACKETRY*				
ARJ1207		Standard Bracket RH		○
ARJ1208		Standard Bracket LH		○
ARJ1707		Swing Away RH Bracket		£135
ARJ1708		Swing Away LH Bracket		£135
ARJ1709		Matrx Resolve Swing in/out Bracket RH		£245
ARJ1710		Matrx Resolve Swing in/out Bracket LH		£245
ARJ2213		Centre Mount Swing Away to RH Nucleus Bracket	(NC 1405)	£495
ARJ2214		Centre Mount Swing Away to LH Nucleus Bracket	(NC 1405)	£495
ARJ2209		REM500 Swing Away RH Multi-Adjustable Bracket	(NC 2203, 2204, 2205)	○
ARJ2210		REM500 Swing Away LH Multi-Adjustable Bracket	(NC 2203, 2204, 2205)	○
ARJ2227		Bracket RH mounted for Attendant		£65
ARJ2228		Bracket LH mounted for Attendant		£65
ARZ2239		Individual mounting option per request (Special). Requires additional information/specification.		POA
DUAL CONTROLS				
ARJ2220		LiNX DLX-REM050 (Dual Control used as secondary remote)	(NC 1221, 1222, 2227, 2228, 2205, 2206, 2245, 2128, 2230, 2247, 2233, 2240)	£395
ARJ2219		LiNX DLX-ACU200 (Compact Dual Control)	(NC 1222, 2227, 2228, 2128, 2201)	£395
ARJ2245		LiNX DLX-IDC (Intuitive Dual Control)	(NC 2228, 2128, 2201)	£745
UNZ9259		mo-vis Scoot Control (Only available via Special)	(NC 2108, 0000, 1221, 1222)	£2,315
DUAL CONTROL BRACKETRY (Only applicable if Dual Control is selected)				
ARJ1207		Dual Control Standard Bracket RH		○
ARJ1208		Dual Control Standard Bracket LH		○
ARJ1221		IDC Dual Control Detachable Bracket RH	(NC 2219, 2227, 2228)	○
ARJ1222		IDC Dual Control Detachable Bracket LH	(NC 2219, 2227, 2228)	○
UNZ9260		Individual mounting option per request (Special). Requires additional information/specification.		POA
COMPACT REMOTES				
ARJ2231		LiNX DLX-CR400 (Compact Remote)	(NC 2203, 2204, 2220, 2235, 2236, 2241, 2242)	£825
ARJ2232		LiNX DLX-CR400LF (Compact Remote Low Force)	(NC 2203, 2204, 2220, 2235, 2236, 2241, 2242)	£1,075
UNZ9250		mo-vis Micro LiNX (Only available via Special)	(NC 2203, 2204, 0000, 2207, 2208, 2213, 2214, 2235, 2236, 2241, 2242, 2220)	£2,600
UNZ9251		mo-vis Micro Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204, 0000, 2207, 2208, 2213, 2214, 2235, 2236, 2241, 2242, 2220)	£2,730
UNZ9252		mo-vis Multi LiNX (Only available via Special)	(NC 2203, 2204, 0000, 2207, 2208, 2213, 2214, 2235, 2236, 2241, 2242, 2220)	£2,320
UNZ9253		mo-vis Multi Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204, 0000, 2207, 2208, 2213, 2214, 2235, 2236, 2241, 2242, 2220)	£2,520
UNZ9254		mo-vis All Round Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204, 0000, 2207, 2208, 2213, 2214, 2235, 2236, 2241, 2242, 2220)	£1,395
UNZ9255		mo-vis All Round Low Force Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204, 0000, 2207, 2208, 2213, 2214, 2235, 2236, 2241, 2242, 2220)	£2,015
UNZ9256		mo-vis Heavy Duty Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204, 0000, 2207, 2208, 2213, 2214, 2235, 2236, 2241, 2242, 2220)	£3,190
ARZ2233		ASL MEC (Proportional MEC Joystick ASL130)	(NC 2203, 2204)	£4,100
ARZ2240		ASL Heavy Duty (Proportional Joystick w mode switch ASL133)	(NC 2203, 2204)	£2,285

COMPACT REMOTE BRACKETRY (Only applicable if Compact Remote is selected)				
ARJ2207		RH Mount	(NC 2203, 2204)	○
ARJ2208		LH Mount	(NC 2203, 2204)	○
ARJ2213		Centre Mount Swing Away to RH (Nucleus)	(NC 2203, 2204)	£495
ARJ2214		Centre Mount Swing Away to LH (Nucleus)	(NC 2203, 2204)	£495
ARZ2235		Tray Mount RH (ASL603, 100mm x 150mm)	(NC 2203, 2204)	○
ARZ2236		Tray Mount LH (ASL603, 100mm x 150mm)	(NC 2203, 2204)	○
ARZ2241		Centre Mount, Swing Away to RH, Nucleus (Inc Tray ASL603, 100mm x 150mm)	(NC 2203, 2204)	£495
ARZ2242		Centre Mount, Swing Away to LH, Nucleus (Inc Tray ASL603, 100mm x 150mm)	(NC 2203, 2204)	£495
ARZ2215		Individual mounting option per request (Special). Requires additional information/specification.	(NC 2203, 2204)	POA
CHIN CONTROLS				
UNZ9250		mo-vis Micro LiNX (Only available via Special)	(NC 2201, 2203, 2204, 2220)	£2,600
UNZ9251		mo-vis Micro Omni (Only available via Special. DLX-IN500 module needed)	(NC 2201, 2203, 2204, 2220)	£2,730
UNZ9252		mo-vis Multi LiNX (Only available via Special)	(NC 2201, 2203, 2204, 2220)	£2,320
UNZ9253		mo-vis Multi Omni (Only available via Special. DLX-IN500 module needed)	(NC 2201, 2203, 2204, 2220)	£2,520
UNZ9254		mo-vis All Round Omni (Only available via Special. DLX-IN500 module needed)	(NC 2201, 2203, 2204, 2220)	£1,395
UNZ9255		mo-vis All Round Low Force Omni (Only available via Special. DLX-IN500 module needed)	(NC 2201, 2203, 2204, 2220)	£2,015
CHIN CONTROL BRACKETRY (Only applicable if Chin Control is selected)				
UNZ9261		mo-vis Chin Control Harness (One Size. Special Request Only)		£960
UNZ9262		mo-vis Manual Swing Away (Special Mounting Only. Please specify what side mounting required)		£1,380
UNZ9263		mo-vis Powered Swing Away (Special Mounting Only. Please specify what side mounting required)		£2,280
HEAD CONTROLS				
ARJ2230		ASL 104 ATOM Head Control with ultrasonic proximity sensors with standard wings	(NC 2179, 2197, 2201, 2203, 2204, 2220, 2233, 2324, 2325, 2327)	£4,220
ARJ2247		ASL 104P PROTON Head Control with ultrasonic proximity sensors with PROTON wings	(NC 2179, 2197, 2201, 2203, 2204, 2220, 2233, 2324, 2325, 2327)	£4,945
UNZ9257		mo-vis Proportional Head Control Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204, 2233, 2324, 2325, 2327, 2240, 2175, 2171, 2170, 2197, 1350, 1432, 1427, 0000, 2220)	£4,775
UNZ9258		mo-vis Switched Head Control Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204, 2233, 2324, 2325, 2327, 2240, 2175, 2171, 2170, 2197, 1350, 1432, 1427, 0000, 2220)	£3,315
FOOT CONTROL				
UNZ9264		mo-vis Foot Control LiNX (Only available via Special)	(NC 2201, 2203, 2204, 2220)	£2,880
UNZ9265		mo-vis Foot Control Omni (Only available via Special. DLX-IN500 module needed)	(NC 2201, 2203, 2204, 2220)	£2,740
FOOT CONTROL BRACKETRY (Only applicable if Foot Control is selected)				
UNZ9266		Individual mounting option per request (Special). Requires additional information/specification.		POA
LiNX 10-WAY SWITCH				
ARJ2251		10-Way Switch DLX-SW10-A for direct Powered Seating Access		£500
LiNX 10-WAY SWITCH BRACKETRY (Only applicable if 10-Way Switch is selected)				
ARJ2217		Bracket RH, mounted for user	(NC 2213)	○
ARJ2218		Bracket LH, mounted for user	(NC 2214)	○
ARJ2227		Bracket RH, mounted for attendant		£65
ARJ2228		Bracket LH, mounted for attendant		£65
LiNX ACCESSORY CABLE				
ARJ2252		Loom providing jack sockets for Drive and Mode (includes connection for 10-way switch)		£100

LINX USB CHARGER PORT			
ARJ2253		USB Charger Port DLX-USB02-A (NC 2701, 2703, 2704, 2705, 2706)	£160
LINX USB CHARGER PORT BRACKETRY (Only applicable if USB Charger is selected)			
ARJ2207		RH Mount	○
ARJ2208		LH Mount	○
FURTHER LINX COMPONENTS / ACCESSORIES			
ARJ2302		Not Necessary (NC 1686, 1687, 1690, 2254, 1711, 1712, 1713, 1714, 1715, 1716)	◆
ARJ2300		To be enclosed only	○
ARZ2301		To be connected and programmed (Special)	POA
BUDDY BUTTONS			
ARJ1685		Red Buddy Button	£90
ARJ1686		Yellow Buddy Button	£90
ARJ1687		Green Buddy Button	£90
ARJ1690		Blue Buddy Button	£90
EGG SWITCH			
ARJ2254		Egg Switch, black	£90
PIKO BUTTONS			
ARJ1711		Red Piko Button	£160
ARJ1712		Yellow Piko Button	£160
ARJ1713		Green Piko Button	£160
ARJ1714		Orange Piko Button	£160
ARJ1715		Blue Piko Button	£160
ARJ1716		Black Piko Button	£160
TOGGLES / BUTTONS			
ARZ2261		Stereo Toggle - direct control for Tilt only	£270
ARZ2262		Stereo Button - direct control for Tilt only	£255
ARZ2263		Four-Way Toggle - direct control for Tilt & Recline	£330
ARZ2264		Four-Way Button - direct control for Tilt & Recline	£355
ARZ2207		Toggle / Button RH Mount	○
ARZ2208		Toggle / Button LH Mount	○
ARZ2225		Without Mounting	○
LINX INPUT MODULE			
UNZ2273		DLX-IN200 Module (Jack port for external power button) (NC 2220)	£270
UNZ9249		DLX-IN500 Module (for DB9 Drive input. Only compatible with DLX-REM400 / DLX-REM500)	£785
BODYPOINT JOYSTICK KNOBS			
ARJ2721		3" U-Shaped, Joystick Handle, Bodypoint PC101A	£75
ARJ2722		4" U-Shaped, Joystick Handle, Bodypoint PC102A	£75
ARJ2723		Domed Rubber, Joystick Handle, Bodypoint PC107A	£75
ARJ2724		Mushroom Grooved, Joystick Handle, Bodypoint PC110A	£75
SHROUD COLOUR*			
ARJ2640		Obsidian Black Frame Shrouds	○
ARJ2641		Ruby Red Frame Shrouds	○
ARJ2642		Purple Passion Frame Shrouds	○
ARJ2643		Peacock Blue Frame Shrouds	○
ARJ2644		Arctic Blue Frame Shrouds	○

RIM COLOUR*			
ARJ2650		Obsidian Black Frame Shrouds	○
ARJ2651		Ruby Red Frame Shrouds	○
ARJ2652		Purple Passion Frame Shrouds	○
ARJ2653		Peacock Blue Frame Shrouds	○
ARJ2654		Arctic Blue Frame Shrouds	○
FORKS*			
ARJ2631		Single sided fork	£300
ARJ2632		Double sided fork	○
CENTRE WHEELS (DRIVE WHEELS)*			
ARJ0704		DW 14" (355 mm) x 3" (75 mm), Pneumatic	○
ARJ0705		DW 14" (355 mm) x 3" (75 mm), Solid inlay, dual density foam insert	○
ARJ0706		DW 14" (355 mm) x 3" (75 mm), Puncture protected, gel-sealed	○
MOTOR TYPE*			
ARJ1294		3 km/h (2mph)	○
ARJ1282		6 km/h (4mph) (4 pole)	◆
ARJ1909		12 km/h (8mph) (includes light system and indicators) G-Trac recommended	£600
BATTERY*			
ARJ1260		73.5 Ah (Gel)	◆
G TRAC			
ARJ1215		G-Trac (ensures the powerchair takes the intended direction as indicated by the user control, without the need for continual correction by the client on uneven or sloping terrain)	£400
CHARGER*			
ARJ1343		OptiCharge Charger 8A (200 - 250V)	◆
ARJ1313		OptiCharge Charger 10A (200 - 250V) Charging time will be 20% faster	£155
LIGHTS*			
ARJ1437		LED Light System with Indicators (30 LUX) (NC 1909)	£450
ARJ1254		Without light system (NC 1909)	◆
ARX1006		Lights Included with Motor Package 10 & 12kmh (NC 1282)	○
TRANSPORTATION KIT*			
ARJ1450		Preparation for DAHL-Docking station (Adequate headrest is compulsory). Includes a mounting bracket in the battery box. Docking Station for car implementation has to be ordered separately with DAHL Engineering.	£580
ARJ1494		4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19). Adequate headrest is compulsory.	◆

FIXATION BELT*			
ARJ2113		Lapbelt (adjustable both sides) adjustable mounting on seat rail extrusion (NC 2702, 2703, 2704, 2705, 2706)	○
ARJ1332		H-Belt	£75
ARJ2701		Bodypoint Evoflex Pelvic, Stabilizer, Rehab Latch size L (padded length 620mm) (NC 2121, 2122, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 2113)	£95
ARJ2702		Bodypoint Evoflex Pelvic, Stabilizer, Rehab Latch size M (padded length 460mm) (NC 2123, 2124, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 2113)	£95
ARJ2703		Bodypoint 2-Point Hip Belt, Centre pull - push button, Cinch mount size L (padded length 620mm) (NC 2121, 2122, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 2113)	£90
ARJ2704		Bodypoint 2-Point Hip Belt, Centre pull - push button, Cinch mount size M (padded length 460mm) (NC 2123, 2124, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 2113)	£90
ARJ2705		Bodypoint 2-Point Hip Belt, Dual pull Side release button, Cinch mount size L (padded length 620mm) (NC 2121, 2122, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 2113)	£90
ARJ2706		Bodypoint 2-Point Hip Belt, Dual pull Side release button, Cinch mount size M (padded length 460mm) (NC 2123, 2124, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 2113)	£90
REAR VIEW MIRROR			
ARJ1415		Rear view Mirror RH (NC 2213)	£45
ARJ1339		Rear view Mirror LH (NC 2214)	£45
BAG HOOKS			
ARJ2128		Bag Hooks (Max load up to 6kg = impact on max user weight)	£80
CANE HOLDER			
ARJ1312		Cane holder RH	£45
ARJ1363		Cane holder LH	£45
TRAY			
ARJ1404		Tray, Swing Away RH (NC 1707, 2213, 2214, 2209)	£200
ARJ1405		Tray, Swing Away LH (NC 1708, 2213, 2214, 2210)	£200

TOTAL PRICE			
-------------	--	--	--

Non Conforming Options:

ARJ0000	TDX SP2 Heavy Duty
ARJ2121	380 mm Seat Plate. Can be reduced to 330 mm in combination with inrigged armrests.
ARJ2122	430 mm Seat Plate. Can be reduced to 380 mm in combination with inrigged armrests.
ARJ2123	480 mm Seat Plate
ARJ2124	530 mm Seat Plate
ARJ2132	410 mm
ARJ2134	460 mm
ARJ2138	560 mm
ARJ2143	Matrx Flo-tech Contour Visco; (2-way stretch material, 430 mm x 430 mm, up to 125kg)
ARJ1980	Tapered Wedges (Pair)
ARJ1981	Abductor
ARJ1982	Lateral Wedge
ARJ1983	Anterior Seat Wedges
ARJ1980	Tapered Wedges (Pair)
ARJ1981	Abductor
ARJ1982	Lateral Wedge
ARJ1983	Anterior Seat Wedges
ARJ1985	Pelvic Obliquity Kit (Contains 2x positioning accessories)
ARJ2163	Backposts short (480 mm). For Tension Adjustable Upholstery/Matrx Backs, measured from seat plate to top of upholstery. Width adjustable cast bracket.
ARJ2108	Backposts WITHOUT push handles, Push handles will be cut down on the vertical section.
ARJ2171	Tension Adjustable Upholstery with Nylon overlay. Short = 480 mm / 7 straps. Long = 540 mm / 8 straps. Measured from seat plate to top of upholstery.
ARJ2179	Flex3 Tension Adjustable (540 mm or 480 mm)
ARJ2172	Flex3 adjustable Backrest Plates (without Backrest Cushion)
ARJ2173	Flex3 adjustable Backrest Plates with Laguna Backrest cushion
ARJ2174	Matrx Posture Back, height 400 mm (Mesh fabric) Corresponds to SW setting chosen
ARJ2175	Matrx Posture Back, height 500 mm (Mesh fabric) Corresponds to SW setting chosen
ARJ2176	Matrx Elite E2 Back , height 500 mm (Mesh fabric) Corresponds to SW setting chosen
ARJ2177	Matrx Elite TR Back, height 500 mm (Mesh fabric) Corresponds to SW setting chosen
ARJ2178	Matrx Elite E2 Deep Back, height 500 mm (Mesh fabric) Corresponds to SW setting chosen
ARJ2170	Without Backrest Upholstery (Back Posts/Back Frame only)
ARJ2151	Black Nylon Fabric TR18
ARJ2152	Black Multi Stretch PU Coated Wipe Down Fabric TR26
ARJ2155	Black Spacer Fabric TR36
ARJ2193	Flip-up Armrest. Height and depth adjustable (height adjustable from 230mm - 360mm). 150 kg max user weight.
ARJ2192	Following Armrest (height adjustable from 230mm - 300mm). 150 kg max user weight
ARJ2194	Following Armrest (height adjustable from 300mm - 360mm). 150 kg max user weight
ARJ2196	Long Armpad (375mm)
ARJ2197	Long Ergonomic Armpad (397mm)

ARJ2199	Soft Armpad (REA Azalea Wide Pad)
ARJ2184	Hip Support LH, pad size 140 x 100 mm
ARJ2186	Hip Support LH, pad size 230 x 100 mm
ARJ2188	Hip Support LH, pad size 330 x 100 mm
ARJ2185	Hip Support RH, pad size 140 x 100 mm
ARJ2187	Hip Support RH, pad size 230 x 100 mm
ARJ2189	Hip Support RH, pad size 330 x 100 mm
ARJ1350	Neckrest depth adjustable
ARJ1432	Headrest sideways adjustable
ARJ1427	Headrest with cheek support
ARJ2324	Elan Standard Pad 255 mm / 10" wide w/standard cover (ESP26). Complete with Elan Headrest Hardware
ARJ2325	Elan Occipital Pad 230x180 mm / 9"x7" w/standard cover (EOP23). Complete with Elan Headrest Hardware
ARJ2327	Elan 4 Point Headrest Pad 280x255 mm / 11"x10" w/standard cover. Complete with Elan Headrest Hardware
ARJ1496	Non Invacare Headrest but profile mount fitting Flex3 hard back
ARJ1445	Lateral Support with extended arm, multi-adjustable pad
ARJ1092	Powered "fixed pivot" Tilt 20° (136kg max user weight)
ARJ1095	Powered "fixed pivot" Tilt 20°, heavy duty (180kg max user weight)
ARJ1093	Powered Tilt (0° - 30°) with CoG shift (160kg max user weight)
ARJ2052	Powered Tilt (0°-50°) with CoG shift (136kg max user weight)
ARJ1094	Seat Lifter with Powered Tilt (300mm, 30°) with CoG shift (160kg max user weight)
ARJ2053	Seat Lifter with Powered Tilt (300mm, 45°) with CoG shift (136kg max user weight)
ARJ0977	Powered Recliner (0° - 30°) Foldable
ARJ0556	Centre Mounted: Manual Centre Mounted Legrest, angle settable. Length 280 - 385 mm
ARJ0589	Centre Mounted: LNX Powered Centre Mounted Legrest, length compensating. Length 340 - 410 mm
UNZ0592	Centre Mounted: LNX Powered Centre Mounted Legrest with Telescoping Foot Platform (lowers to the ground). Length: 380-419mm (Selectable length depends on STF height) Custom-made, requires customised setting.
ARJ0526	Side Mounted: Standard hanger 80° pair. Lower leg length: 350 - 450 mm
ARJ0503	Side Mounted: Vari F pair. Lower leg length: 290 - 460 mm. Angle adjustment: 70 - 0° with tools
ARJ0533	Side Mounted: Vari A pair. Lower leg length: 290 - 460 mm. Angle adjustment: 70 - 0° with handwheel and calf pad
ARJ0531	Side Mounted: ADE pair. Lower leg length: 290 - 460 mm. Powered angle adjustment: 80 - 0° with calf pad
ARJ0501	Without Leg Rests
ARJ2365	Angle & width adjustable footplates S (127 W x 190 mm D / 5"x7.5")
ARJ0559	Angle & width adjustable footplates M (140 W x 240 mm D / 5.5"x9.5")
ARJ0516	Angle & width adjustable footplates L (150 W x 290 mm D / 6"x11.5)
ARJ0517	Angle adjustable foot platform M (305 W x 275 mm D / 12"x10.75")
ARJ0512	Fixed Footplate pair
ARJ0504	Angle Adjustable Footplate pair
ARJ0513	Angle and depth adjustable Footplate pair
ARJ0501	Without footplates for option 0501 without Leg Rests
ARJ0554	Amputee Support RH
ARJ0555	Amputee Support LH
ARJ2201	DLX-REM110 (Drive Only)

ARJ2203	DLX-REM211 (Drive & Actuators)
ARJ2204	DLX-REM216 (Drive, Actuators & Lights)
ARJ2205	DLX-REM400 (Touch Screen Interface. Drive, Act. & Lights. Configurable jack sockets)
ARJ2206	DLX-REM500 (Touch Screen Interface Display Unit - requires a user input/alternative joystick).
ARJ2270	Assignable Toggle Switch Module for REM400 (DLX-SW400TGL)
ARJ1208	Standard Bracket LH
ARJ1707	Swing Away RH Bracket
ARJ1708	Swing Away LH Bracket
ARJ1709	Matrx Resolve Swing in/out Bracket RH
ARJ1710	Matrx Resolve Swing in/out Bracket LH
ARJ2213	Centre Mount Swing Away to RH Nucleus Bracket
ARJ2214	Centre Mount Swing Away to LH Nucleus Bracket
ARJ2209	REM500 Swing Away RH Multi-Adjustable Bracket
ARJ2210	REM500 Swing Away LH Multi-Adjustable Bracket
ARJ2227	Bracket RH mounted for Attendant
ARJ2228	Bracket LH mounted for Attendant
ARJ2220	LiNX DLX-REM050 (Dual Control used as secondary remote)
ARJ2219	LiNX DLX-ACU200 (Compact Dual Control)
ARJ2245	LiNX DLX-IDC (Intuitive Dual Control)
ARJ1208	Dual Control Standard Bracket LH
ARJ1221	IDC Dual Control Detachable Bracket RH
ARJ1222	IDC Dual Control Detachable Bracket LH
ARJ2231	LiNX DLX-CR400 (Compact Remote)
ARJ2232	LiNX DLX-CR400LF (Compact Remote Low Force)
ARZ2233	ASL MEC (Proportional MEC Joystick ASL130)
ARZ2240	ASL Heavy Duty (Proportional Joystick w mode switch ASL133)
ARJ2207	RH Mount
ARJ2208	LH Mount
ARJ2213	Centre Mount Swing Away to RH (Nucleus)
ARJ2214	Centre Mount Swing Away to LH (Nucleus)
ARZ2235	Tray Mount RH (ASL603, 100mm x 150mm)
ARZ2236	Tray Mount LH (ASL603, 100mm x 150mm)
ARZ2241	Centre Mount, Swing Away to RH, Nucleus (Inc Tray ASL603, 100mm x 150mm)
ARZ2242	Centre Mount, Swing Away to LH, Nucleus (Inc Tray ASL603, 100mm x 150mm)
ARZ2215	Individual mounting option per request (Special). Requires additional information/specification.
ARJ2230	ASL104 ATOM Head Control with ultrasonic proximity sensors with standard wings
ARJ2247	ASL104P PROTON Head Control with ultrasonic proximity sensors with PROTON wings
ARJ2227	Bracket RH, mounted for attendant
ARJ2228	Bracket LH, mounted for attendant
ARJ2253	USB Charger Port DLX-USB02-A
ARJ2207	RH Mount
ARJ2208	LH Mount
ARJ1686	Yellow Buddy Button
ARJ1687	Green Buddy Button

ARJ1690	Blue Buddy Button
ARJ2254	Egg Switch, black
ARJ1711	Red Piko Button
ARJ1712	Yellow Piko Button
ARJ1713	Green Piko Button
ARJ1714	Orange Piko Button
ARJ1715	Blue Piko Button
ARJ1716	Black Piko Button
ARZ2207	Toggle / Button RH Mount
ARZ2208	Toggle / Button LH Mount
ARJ1282	6 km/h (4mph) (4 pole)
ARJ1909	12 km/h (8mph) (includes light system and indicators) G-Trac recommended
ARJ2113	Lapbelt (adjustable both sides) adjustable mounting on seat rail extrusion
ARJ1332	H-Belt
ARJ2701	Bodypoint Evoflex Pelvic, Stabilizer, Rehab Latch size L (padded length 620mm)
ARJ2702	Bodypoint Evoflex Pelvic, Stabilizer, Rehab Latch size M (padded length 460mm)
ARJ2703	Bodypoint 2-Point Hip Belt, Centre pull - push button, Cinch mount size L (padded length 620mm)
ARJ2704	Bodypoint 2-Point Hip Belt, Centre pull - push button, Cinch mount size M (padded length 460mm)
ARJ2705	Bodypoint 2-Point Hip Belt, Dual pull Side release button, Cinch mount size L (padded length 620mm)
ARJ2706	Bodypoint 2-Point Hip Belt, Dual pull Side release button, Cinch mount size M (padded length 460mm)
ARJ2128	Bag Hooks (Max load up to 6kg = impact on max user weight)
ARJ1405	Tray, Swing Away LH



For a more complex customer requirement,
please email: orders@invacare.com.au

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):