



Yes, you can:

Tel: 01656 776222
 Email: ordersuk@invacare.com
 Online Spares available at: www.invacare.co.uk

Manual Retail Prescription Form

January 2026
 Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

Rea Azalea Self Propel

Standard Features

- Max User Weight 135KG/21 Stone
- Crash Tested ISO 7176-Part 19
- 24" or 22" Self Propulsion with User Manoeuvred Tyre Brakes and Solid Black Tyres
- Choice of Seat Widths 390mm, 440mm or 490mm (+50mm Width Adjustment)
- Choice of 2 Seat Heights 400mm or 450mm
- Depth Adjustable Seat 430mm-500mm
- Without Seat Cushion, Supplied with Seat Board
- Flex 3 Backrest Plate, Width and Height Adjustable
- Passad 2 Backrest Cushion with Multi-stretch Black PU coated fabric Cover and Push Handles Connected
- Assistant Operated Gas Piston Seat Tilt +1° - +24°
- Backrest Recline 30° with Tilt Lock
- Angle Adjustable Legrests with Calf Support, Height and Depth Adjustable Footplates
- Height Adjustable Ducktail Pad Armrests with Autolock
- Pearl Silver Frame Colour, Solid Black Front Castors
- Quick Release Self Propel Wheels, Solid Tyres and Aluminium Handrims
- Pelvic Belt, Adjustable both sides
- Anti Tippers



Standard Base price: £2,393

Please make your selection with a cross in the box ○ Free of charge ◆ standard feature POA special request

Customer weight = KG

Code	Description		
BASE PRODUCT*			
BND9001	Rea Azalea Self Propel		£2,393
BND8000	SPECIAL REQUEST - PLEASE INCLUDE QUOTATION REFERENCE IN SPECIAL DESCRIPTION AT THE END OF THE FORM WHERE APPLICABLE		POA
PROPULSION*			
BND0025	Self Propel 24"/22" wheels with user brakes	(NC 1076)	◆
BND0030	Self Propel 24"/22" wheels with assistant hub brakes	(NC 1074)	£140

SEAT WIDTH*				
BND0040	<input type="checkbox"/>	390 mm/15" (Small chassis)	(NC 0244, 0245, 0246, 0247)	<input type="radio"/>
BND0050	<input type="checkbox"/>	440 mm/17" (Medium chassis)	(NC 0241, 0242, 0246, 0247)	<input type="radio"/>
BND0060	<input type="checkbox"/>	490 mm/19" (Large chassis)	(NC 0241, 0242, 0243, 0244)	<input type="radio"/>
SEAT WIDTH EXTENSION*				
BND0065	<input type="checkbox"/>	+ 50 mm seat width pre-adjustment on armrests and legrests		<input type="radio"/>
BND0085	<input type="checkbox"/>	Without seat width extension, set at 0mm adjustment but can be adjusted wider if required		<input type="radio"/>
SEAT DEPTH*				
BND0100	<input type="checkbox"/>	430-500 mm / 17" - 20"		<input checked="" type="radio"/>
SEAT HEIGHT*				
BND0110	<input type="checkbox"/>	400 mm / 16"	(NC 0370, 1015, 1025, 1062)	<input type="radio"/>
BND0120	<input type="checkbox"/>	450 mm / 18"	(NC 1055)	<input type="radio"/>
SEAT ANGLE				
BND0127	<input type="checkbox"/>	Negative Seat Angle -1 To +23° (For attendant operated gas piston in negative tilt)	(NC 0360)	<input type="radio"/>
SEAT TYPE*				
BND0130	<input type="checkbox"/>	Depth Adjustable		<input checked="" type="radio"/>
SEAT CUSHION*				
BND0165	<input type="checkbox"/>	Cushion Floshape seat upholstery	(NC 1111, 0182, 0183, 0186, 0187, 0188, 0189)	£57
SWO0190	<input type="checkbox"/>	Cushion Floshape seat upholstery Shorter depth 450mm deep. Only use with seat depth in shortest position	(NC 1111, 0182, 0183, 0186, 0187, 0188, 0189)	£57
BND0167	<input type="checkbox"/>	Cushion Vicair Adjuster O2 10cm, complete with Black Comfair cover (Additional 5 day leadtime)	(NC 0182, 0183, 0186, 0187, 0188, 0189, 0191, 0201, 0202)	£547
BND0168	<input type="checkbox"/>	Cushion Matrix Libra, complete with Black cover	(NC 0191, 0201, 0202, 9993)	£410
BND0169	<input type="checkbox"/>	Cushion Flo-Tech Solution Xtra, complete with Black cover	(NC 0191, 0201, 0202, 9993)	£398
BND0179	<input type="checkbox"/>	Cushion Vicair Vector O2 100mm with 2 top covers (Additional 15 day leadtime)	(NC 0182, 0183, 0186, 0187, 0188, 0189, 0191, 0201, 0202)	£542
BND0180	<input type="checkbox"/>	Without Seat Cushion, Seat Plate only	(NC 0182, 0183, 0191, 0201, 0202)	<input type="radio"/>
SEAT DEPTH OF MATRIX LIBRA AND FLO-TECH SOLUTION XTRA CUSHIONS (OPTIMISED CUSHION DEPTH TO SUIT SEAT DEPTH SETTING)*				
BND0182	<input type="checkbox"/>	Short cushion depth: Libra delivered at 405mm deep, Solution Xtra delivered at 430mm deep		<input type="radio"/>
BND0183	<input type="checkbox"/>	Long cushion depth: Libra delivered at 455mm deep, Solution Xtra delivered at 455mm deep		<input type="radio"/>
BND9993	<input type="checkbox"/>	Only select this option when ordering Floshape, Vicair or without seat cushions		<input type="radio"/>
MATRIX LIBRA/FLOTECH SOLUTION XTRA CUSHION ACCESSORIES				
BND0186	<input type="checkbox"/>	Tapered Wedge Pair		£30
BND0187	<input type="checkbox"/>	Abductor		£23
BND0188	<input type="checkbox"/>	Lateral Side Wedge		£20
BND0189	<input type="checkbox"/>	Anterior Front Wedge		£27
SEAT UPHOLSTERY*				
BND0191	<input type="checkbox"/>	Black Plush TR35		<input type="radio"/>
BND0201	<input type="checkbox"/>	Black Multi-stretch Black PU coated fabric		<input checked="" type="radio"/>
BND0202	<input type="checkbox"/>	Black Spacer TR36		<input type="radio"/>
BND1111	<input type="checkbox"/>	No upholstery, select with options 0167, 0168, 0169, 0179 or 0180 Without seat cushion		<input type="radio"/>
BACKREST WIDTH EXTENSION*				
BND0205	<input type="checkbox"/>	+ 50 mm Backrest width pre-adjustment, set at +50mm but can be adjusted if required	(NC 0221, 0223, 0225)	<input type="radio"/>
BND0207	<input type="checkbox"/>	Without Backrest width pre-adjustment, set at 0mm but can be adjusted if required		<input type="radio"/>

BACKREST*				
BND0210		Flex 3 Backrest, Fixed Plate (Width Adjustable)	(NC 0233, 0234, 0235, 0240, 0241, 0242, 0243, 0244, 0245, 0246, 0247, 0343)	◆
BND0219		New!! Flex 3+ Backrest features tension adjustment, width adjustment, and fore/aft positioning of the vertical backpost tubes to accommodate a wider range of clinical requirements.	(NC 0233, 0235, 0234, 0241, 0242, 0243, 0244, 0245, 0246, 0247, 0277, 0343)	£225
BND0220		Flex 3 Backrest features tension adjustment and width adjustment of the vertical backpost tubes.	(NC 0233, 0235, 0234, 0241, 0242, 0243, 0244, 0245, 0246, 0247, 0277, 0343)	£102
BND0221		Matrx Elite E2 (includes hardware mounted to frame, backrest supplied in box)	(NC 0240, 0250, 0260, 0270, 0305, 0311, 0321, 0322, 0800, 0810, 0820, 0830, 0880, 0890, 0900, 0920)	£494
BND0223		Matrx Elite E2 Deep (includes hardware mounted to frame, backrest supplied in box)	(NC 0240, 0250, 0260, 0270, 0305, 0311, 0321, 0322, 0800, 0810, 0820, 0830, 0880, 0890, 0900, 0920)	£546
BND0225		Matrx Posture Back (includes hardware mounted to frame, backrest supplied in box)	(NC 0235, 0240, 0241, 0243, 0244, 0245, 0246, 0247, 0250, 0260, 0270, 0305, 0311, 0321, 0322, 0335, 0800, 0810, 0820, 0830, 0880, 0890, 0900, 0920)	£353
BND0230		Backrest Tubes Without Backrest	(NC 0233, 0234, 0235, 0240, 0241, 0242, 0243, 0244, 0245, 0246, 0247, 0250, 0260, 0270, 0305, 0311, 0321, 0322, 0345, 0740, 0800, 0810, 0820, 0830, 0880, 0890, 0900, 0920)	○
MATRX BACKREST HEIGHT				
BND0233		Matrx Backrest Height 410 mm / 16"		○
BND0235		Matrx Backrest Height 460 mm / 18"		○
BND0234		Matrx Backrest Height 510 mm / 20"		○
MATRX BACKREST WIDTH				
BND0241		Matrx Backrest Width 360mm / 14" (Chest width 265mm / 10 1/2")		○
BND0242		Matrx Backrest Width 380mm / 15" (Chest width 290mm / 11 1/2")		○
BND0243		Matrx Backrest Width 410mm / 16" (Chest width 315mm / 12 1/2")		○
BND0244		Matrx Backrest Width 430mm / 17" (Chest width 340mm / 13")		○
BND0245		Matrx Backrest Width 460mm / 18" (Chest width 370mm / 14 1/2")		○
BND0246		Matrx Backrest Width 480mm / 19" (Chest width 390mm / 15")		○
BND0247		Matrx Backrest Width 510mm / 20" (Chest width 420mm / 16 1/2")		○
BACKREST CUSHION*				
BND0240		Cover for tension adjustable upholstery	(NC 2222)	○
BND0250		Cushion, Laguna (Lateral Support)	(NC 2222)	○
BND0260		Cushion, Mistral 2 (Waist Support)	(NC 0322, 2222)	£144
BND0270		Cushion, Passad 2 (Shoulder Support)	(NC 2222)	◆
BND0276		Vicair Multifunctional O2 pressure relieving backrest cushion, supplied with Black mesh cover (Additional 5 day lead time)	(NC 0305, 0311, 0321, 0322)	£564
BND0277		Vicair Multifunctional O2 Flex pressure relieving backrest cushion, supplied with Black mesh cover (Additional 5 day lead time)	(NC 0305, 0311, 0321, 0322)	£622
BND0300		Without Backrest Cover/ Cushion	(NC 0305, 0311, 0321, 0322)	○
WIDER BACKREST COVER				
BND0305		Extra large backrest cover or cushion. One size wider backrest cushion or tension cover fitted if selected		○
BACKREST UPHOLSTERY*				
BND0311		Black Plush TR35		○
BND0321		Black Multi-stretch Black PU coated fabric		◆
BND0322		Black Spacer TR36		○
BND2222		No upholstery, select with Matrx, and Without backrest options		○

PUSH HANDLES*			
BND0335		Push Handles connected	◆
BND0343		Push Handles connected,with Rea head/neckrest attachment	£72
BND0345		Push Bar, angle adjustable	£102
SEAT TILT AND BACKREST REGULATION*			
BND0350		Assistant operated gas piston seat tilt and backrest recline (NC 5555)	◆
BND0360		Electrical operated seat tilt and backrest recline (NC 1263)	£616
TILT LOCK*			
BND1263		Tilt Lock (locking of pre-set angle)	◆
BND5555		No Tilt Lock, Select with 0360 Electric operated seat and backrest	○
CENTRAL LEGREST			
BND0370		Central Legrest, includes calfpad and footrest. Also select options 0420, 0470, 3333, 4444 (NC 0376, 0378, 0380, 0390, 0400, 0410, 0430, 0440, 0450, 0460, 0475, 0476, 0477, 0720, 0930, 0940, 0941, 0942, 0952, 0954, 0956, 1290)	£303
LEGREST HANGER ATTACHMENT*			
BND0375		Standard position for legrest mounting	◆
SWO0189		Standard position with 0-30° Legrest Abduction function (NC 0370, 0475, 0476, 0477)	£100
BND0376		High position 50mm higher for leg length compensation	○
BND0378		90° Hanger attachment to achieve a 90° legrest position on the Elevating Legrest (NC 0390, 0400, 0410, 0440, 0450, 0460)	○
LEGREST HANGER LH*			
BND0380		Elevating Legrest, angle adjustable 0 to 75° swing in/swing out with knee pad, calf support with locking clip LH (NC 0475, 0490, 0942, 0945, 6666)	◆
BND0390		Legrest Hanger fixed 80 deg, swing in/swing out LH (NC 0928, 0929, 0930)	○
BND0400		Legrest Hanger fixed 90 deg, swing in/swing out LH (NC 0520, 0928, 0929, 0930)	○
BND0410		Legrest Amputee LH (NC 0475, 0476, 0477, 0490, 0520, 0928, 0929, 0930, 0941, 0942, 0945, 0952, 0954, 0956)	£137
BND0420		Without Legrest, LH (NC 0475, 0476, 0477, 0490, 0520, 0928, 0929, 0930, 0941, 0942, 0945, 0952, 0954, 0956)	○
LEGREST HANGER RH*			
BND0430		Elevating Legrest, angle adjustable 0 to 75°swing In/swing out with knee pad and calf support with locking clip RH (NC 0475, 0550, 0942, 0945, 6666)	◆
BND0440		Legrest Hanger fixed 80 deg, swing in/swing out RH (NC 0931, 0932, 0940)	○
BND0450		Legrest Hanger fixed 90 deg, swing in/swing out RH (NC 0580, 0931, 0932, 0940)	○
BND0460		Legrest Amputee RH (NC 0475, 0476, 0477, 0550, 0580, 0931, 0932, 0940, 0941, 0942, 0945, 0952, 0954, 0956)	£137
BND0470		Without Legrest, RH (NC 0475, 0476, 0477, 0550, 0580, 0931, 0932, 0940, 0941, 0942, 0945)	○
FOOTREST LH*			
BND0475		Footrest one piece, depth and angle adjustable short tube (NC 0477, 0490, 0520, 0550, 0580, 0930, 0940, 0952, 0954, 0956)	£118
BND0476		Footrest one piece, depth and angle adjustable long tube (NC 0477, 0490, 0520, 0550, 0580, 0930, 0940, 0952, 0954, 0956)	£118
BND0477		Footbox, angle adjustable, padded with short tube (NC 0490, 0520, 0550, 0580, 0930, 0940, 0941, 0942, 0945, 0952, 0954, 0956)	£543
BND0490		Footrest lower part, depth and angle adjustable with short tube and heel strap, LH. Not compatible with Elevating legrest.	○
BND0520		Footrest lower part, depth and angle adjustable with long tube and heel strap, LH. All knee to heel ranges can be achieved when used with Elevating legrests.	○
BND3333		Without footrest, only select with option 0370 Central legrest, 0410 Amputee legrest and 0420 Without legrest	○

FOOTREST RH*			
BND0550		Footrest lower part, depth and angle adjustable with short tube and heel strap, RH. Not compatible with Elevating legrest.	○
BND0580		Footrest lower part, depth and angle adjustable with long tube and heel strap, RH. All knee to heel ranges can be achieved when used with Elevating legrests.	○
BND4444		Without footrest, only select with options 0370 Central legrest, 0460 Amputee legrest, 0470 Without legrest, 0475 and 0476 One piece footrests or 0477 Footbox	○
ARMREST LH*			
BND0600		Armrest Ducktail Pad LH, height adjustable with autolock (NC 1232)	◆
BND0607		Armrests Reinforced Wide Pad, LH, height adjustable with autolock (NC 1230, 1231, 1232, 1250, 1260)	£89
BND0610		Armrest Padded LH, height adjustable with autolock (NC 0680, 0695, 1230, 1231, 1232, 1250, 1260, 9999)	£64
BND0615		Armrest Hemiplegic LH, height adjustable with autolock (NC 0695, 1230, 1231, 1232, 1250, 1260)	£179
BND0620		Without Armrest LH Select when choosing BND1232 Half Tray LH, Narrow Armrest is included (NC 0695, 1230, 1231, 1250, 1260)	○
ARMREST RH*			
BND0630		Armrest Ducktail Pad RH, height adjustable with autolock (NC 1233)	◆
BND0637		Armrest Reinforced Wide Pad, RH, height adjustable with autolock (NC 1230, 1231, 1233, 1250, 1260)	£89
BND0640		Armrest Padded RH, height adjustable with autolock. (NC 0680, 0695, 1230, 1231, 1233, 1250, 1260, 9999)	£64
BND0645		Armrest Hemiplegic, RH, height adjustable with autolock (NC 0695, 1230, 1231, 1233, 1240, 1245, 1250, 1260)	£179
BND0650		Without Armrest RH Select when choosing BND1233 Half Tray RH, Narrow Armrest is included (NC 0695, 1230, 1231, 1240, 1245, 1250, 1260)	○
ARMREST OPTIONS			
BND0680		Side Support Cover, PU fabric cover for the sideguards (Removable with tools) (NC 0695, 9999)	£51
BND9999		Side Support Cover with foam cushion. 2cm thick cushion supplied complete with PU fabric cover. Reduces seat width by 2cm for a comfortable fit. Consists of options 0680 and 0690 (Removable with tools) (NC 0695)	£51
BND0695		Pad for sideguard, removable with hook and loop fastening	£85
SEAT OPTIONS			
BND0710		Incontinence Cover	£40
BND0720		Pommel Abduction Cushion (NC 0722)	£118
BND0722		2 cm Increased seat height plate	£64
BND0725		Tilt Scale for seat tilt	£44
HEADREST OPTIONS*			
BND0730		Headrest, Sideway Adjustable Fixed hardware included (NC 0742, 0744, 0746, 0750, 0770, 0782)	◆
BND0750		Headrest with cheek support, sideway adjustable Fixed hardware included (NC 0740, 0742, 0744, 0746, 0770, 0782)	£107
BND0731		Matrx Elan Standard Pad (254 mm / 10") (NC 0740, 0770, 0782, 0790, 8888)	£124
BND0732		Matrx Elan Standard Pad (355 mm / 14") (NC 0740, 0770, 0782, 0790, 8888)	£124
BND0733		Matrx Elan Occipital Pad (229 mm / 9") (NC 0740, 0770, 0782, 0790, 8888)	£124
BND0734		Matrx Elan Occipital Pad (305 mm / 12") (NC 0740, 0770, 0782, 0790, 8888)	£124
BND0735		Matrx Elan 4-Point Pad (229W x 177H mm / 9"W x 7"H) (NC 0740, 0770, 0782, 0790, 8888)	£232
BND0736		Matrx Elan 4-Point Pad (280W x 254H mm / 11"W x 10"H) (NC 0740, 0770, 0782, 0790, 8888)	£232
BND9991		Without headrest (NC 0742, 0744, 0746)	○

HEADREST HARDWARE*				
BND0742		Matrx Elan Hardware	(NC 0740, 0770, 0790, 1206, 1207, 1300)	£184
BND0744		Matrx Elan Extended Hardware (Extra link)	(NC 0740, 0770, 0790, 1206, 1207, 1300)	£284
BND0746		Matrx Elan LOXX Hardware	(NC 0740, 0770, 0790, 1206, 1207, 1300)	£166
BND8888		Without Elan headrest hardware (Only select with option 0730 & 0750 Headrests)		○
OTHER HEADREST OPTIONS				
BND0740		Cover for headrest		£31
BND0770		Neckrest		○
BND9992		Without neckrest		○
BND0782		Alternative headrest/neckrest mounting bracket		£23
BND0790		Handles Clinic Kit for Head-/Neckrest (Levers which replace screws for easy adjustment)		○
TRUNK SUPPORT LH				
BND0800		Trunk Support Multi Adjustable with standard arm LH	(NC 0900)	○
BND0820		Trunk Support Multi Adjustable with extended height arm LH	(NC 0900)	○
BND0880		NEW!! Trunk Support Swingaway LH, multi adjustable	(NC 0900)	£105
TRUNK SUPPORT RH				
BND0810		Trunk Support Multi Adjustable with standard arm RH	(NC 0900)	○
BND0830		Trunk Support Multi Adjustable with extended height arm RH	(NC 0900)	○
BND0890		NEW!! Trunk Support Swingaway RH, multi adjustable	(NC 0900)	£105
BACKREST OPTIONS				
BND0900		Trunk Support Attachment to mount laterals later. If trunk supports are chosen, the attachment bracket is already included		○
BND0915		Tilt Scale for backrest recline		£44
BND0920		Wedges (One Pair)		£11
ELEVATING LEGREST OPTIONS*				
BND0928		Calf pad cover, for LH elevating legrest		○
BND0929		Knee pad for LH elevating legrest		◆
BND0931		Calf pad cover, for RH elevating legrest		○
BND0932		Knee pad for RH elevating legrest		◆
BND6666		Without Elevating Legrest Options (Select with Fixed legrests, Amputee legrests and Without legrests only)		○
LEGREST OPTIONS				
BND0930		Footplate Extender LH	(NC 0941)	£34
BND0940		Footplate Extender RH	(NC 0941)	£34
BND0941		Cover Footplate, Black Plush Fabric		£54
BND0942		Calf Strap Padded		£73
BND0945		Calf Strap (for One piece footrest)		○
BODYPOINT ANKLE HUGGERS				
BND0952		Bodypoint Ankle Huggers Large (Supplied with fixing kit, unmounted)		£114
BND0954		Bodypoint Ankle Huggers Medium (Supplied with fixing kit, unmounted)		£114
BND0956		Bodypoint Ankle Huggers Small (Supplied with fixing kit, unmounted)		£114

CHASSIS COLOUR*			
BND0944	<input type="checkbox"/>	Silver Pearl	◆
BND0946	<input type="checkbox"/>	Azurite (Dark Blue)	○
BND0947	<input type="checkbox"/>	Sand (Beige)	○
BND0948	<input type="checkbox"/>	Electric Green	○
BND0949	<input type="checkbox"/>	Happy Red	○
CHASSIS SAFETY			
BND0950	<input type="checkbox"/>	Chassis extender, castors extend 50mm forward for added stability	○
CASTOR ATTACHMENT*			
BND0960	<input type="checkbox"/>	Fixed castor attachment	◆
FRONT CASTORS*			
BND0992	<input type="checkbox"/>	150x30 mm (6" x 1 1/4") Medium Hard Black tyre	○
BND1015	<input type="checkbox"/>	200x32 mm (8" x 1 1/4") Black Solid Tyre	£21
BND1025	<input type="checkbox"/>	200x45 mm (8" X 1 3/4") Black Solid tyre	◆
REAR WHEEL SIZE*			
BND1055	<input type="checkbox"/>	22" / 560 mm Self Propel rear wheel with QR attachment	○
BND1062	<input type="checkbox"/>	24" / 600 mm Self Propel rear wheel with QR attachment	○
REAR WHEEL TYPE*			
BND1074	<input type="checkbox"/>	Standard Self Propel with colourless hub, clear anodized rim and Stainless steel spokes	◆
BND1076	<input type="checkbox"/>	Hub Brake 22" or 24" Self Propel with colourless hub, clear anodized rim and Stainless steel spokes	○
REAR WHEEL TYRE*			
BND1082	<input type="checkbox"/>	Pneumatic Profile Black 1 3/8" tyres	○
BND1092	<input type="checkbox"/>	Solid Black tyres	◆
HANDRIM*			
BND1100	<input type="checkbox"/>	Aluminium handrim	◆
BND1120	<input type="checkbox"/>	Plastic coated handrim	£48
SPOKEGUARD			
BND1160	<input type="checkbox"/>	Frost spokeguard	£166
BRAKE*			
BND1170	<input type="checkbox"/>	User manoeuvred tyre brakes, includes extended handles	◆
BND1180	<input type="checkbox"/>	One Arm Brake LH	£126
BND1185	<input type="checkbox"/>	One Arm Brake RH	£126
PELVIC BELT*			
BND1200	<input type="checkbox"/>	Pelvic belt	◆
BND1201	<input type="checkbox"/>	Bodypoint 2-Point Centre-pull Push Button Medium (Mounted to chair)	£83
BND1202	<input type="checkbox"/>	Bodypoint 2-Point Centre-pull Push Button Large (Mounted to chair)	£83
BND1203	<input type="checkbox"/>	Bodypoint 4-Point Centre-pull Push Button Medium (Mounted to chair)	£98
BND1204	<input type="checkbox"/>	Bodypoint 4-Point Centre-pull Push Button Large (Mounted to chair)	£98
BODYPOINT HARNESS			
BND1206	<input type="checkbox"/>	Bodypoint Pivotfit Harness Medium (Supplied with fixing kit, unmounted)	£109
BND1207	<input type="checkbox"/>	Bodypoint Pivotfit Harness Large (Supplied with fixing kit, unmounted)	£109

TRANSPORTATION BRACKETS*			
BND1225		Transportation Brackets, mounted to the front of the wheelchair	◆
OTHER OPTIONS			
BND1230		Fixed Table Tray without elbow pads includes mounting kit (NC 1237, 1260)	£138
BND1231		Fixed Table Tray with elbow Pads includes mounting kit (NC 1237, 1260)	£138
BND1232		Half table tray, clear plastic, flip up, depth adjustable, left side (armrest with straight narrow pad included) (NC 1235, 1260)	£103
BND1233		Half table tray, clear, plastic, flip up, depth adjustable, right side (armrest with straight narrow pad included) (NC 1235, 1260)	£103
BND1240		Swing away table without elbow pads (NC 1237, 1250, 1260)	£152
BND1245		Swing away table with elbow pads (NC 1237, 1250, 1260)	£152
BND1235		Table tray padded cushion (To be ordered in combination with the fixed & swingaway table tray)	£51
BND1237		Half table tray padded cushion	£46
BND1250		Locking device for table tray	£15
BND1260		Mounting kit for table (if ordered separately)	£9
BND1290		Handles Clinic Kit for backrest, armrest, legrest and seat frame (Levers which replace screws for easy adjustment)	£55
BND1300		Harness attachment , must be selected in combination with Head/Neckrest	£81
BND1310		Attachment for positioning belt	£32

TOTAL PRICE	
--------------------	--

Non Conforming Options:

BND0182	Short cushion depth: Libra delivered at 405mm deep, Solution Xtra delivered at 430mm deep
BND0183	Long cushion depth: Libra delivered at 455mm deep, Solution Xtra delivered at 455mm deep
BND9993	Only select this option when ordering Floshape, Vicair or without seat cushions
BND0186	Tapered Wedge Pair
BND0187	Abductor
BND0188	Lateral Side Wedge
BND0189	Anterior Front Wedge
BND0191	Black Plush TR35
BND0201	Black Multi-stretch Black PU coated fabric
BND0202	Black Spacer TR36
BND1111	No upholstery, select with options 0167, 0168, 0169, 0179 or 0180 Without seat cushion
BND0221	Matrx Elite E2 (includes hardware mounted to frame, backrest supplied in box)
BND0223	Matrx Elite E2 Deep (includes hardware mounted to frame, backrest supplied in box)
BND0225	Matrx Posture Back (includes hardware mounted to frame, backrest supplied in box)
BND0233	Matrx Backrest Height 410 mm / 16"
BND0235	Matrx Backrest Height 460 mm / 18"
BND0234	Matrx Backrest Height 510 mm / 20"
BND0241	Matrx Backrest Width 360mm / 14" (Chest width 265mm / 10 1/2")
BND0242	Matrx Backrest Width 380mm / 15" (Chest width 290mm / 11 1/2")
BND0243	Matrx Backrest Width 410mm / 16" (Chest width 315mm / 12 1/2")
BND0244	Matrx Backrest Width 430mm / 17" (Chest width 340mm / 13")
BND0245	Matrx Backrest Width 460mm / 18" (Chest width 370mm / 14 1/2")
BND0246	Matrx Backrest Width 480mm / 19" (Chest width 390mm / 15")
BND0247	Matrx Backrest Width 510mm / 20" (Chest width 420mm / 16 1/2")
BND0240	Cover for tension adjustable upholstery
BND0250	Cushion, Laguna (Lateral Support)
BND0260	Cushion, Mistral 2 (Waist Support)
BND0270	Cushion, Passad 2 (Shoulder Support)
BND0277	Vicair Multifunctional O2 Flex pressure relieving backrest cushion, supplied with Black mesh cover (Additional 5 day lead time)
BND0305	Extra large backrest cover or cushion. One size wider backrest cushion or tension cover fitted if selected
BND0311	Black Plush TR35
BND0321	Black Multi-stretch Black PU coated fabric
BND0322	Black Spacer TR36
BND2222	No upholstery, select with Matrx, and Without backrest options
BND0335	Push Handles connected
BND0343	Push Handles connected,with Rea head/neckrest attachment
BND0345	Push Bar, angle adjustable
BND0360	Electrical operated seat tilt and backrest recline
BND1263	Tilt Lock (locking of pre-set angle)
BND5555	No Tilt Lock, Select with 0360 Electric operated seat and backrest
BND0370	Central Legrest, includes calfpad and footrest. Also select options 0420, 0470, 3333, 4444
SWO0189	Standard position with 0-30° Legrest Abduction function
BND0376	High position 50mm higher for leg length compensation

BND0378	90° Hanger attachment to achieve a 90° legrest position on the Elevating Legrest
BND0380	Elevating Legrest, angle adjustable 0 to 75° swing in/swing out with knee pad, calf support with locking clip LH
BND0390	Legrest Hanger fixed 80 deg, swing in/swing out LH
BND0400	Legrest Hanger fixed 90 deg, swing in/swing out LH
BND0410	Legrest Amputee LH
BND0430	Elevating Legrest, angle adjustable 0 to 75° swing In/swing out with knee pad and calf support with locking clip RH
BND0440	Legrest Hanger fixed 80 deg, swing in/swing out RH
BND0450	Legrest Hanger fixed 90 deg, swing in/swing out RH
BND0460	Legrest Amputee RH
BND0475	Footrest one piece, depth and angle adjustable short tube
BND0476	Footrest one piece, depth and angle adjustable long tube
BND0477	Footbox, angle adjustable, padded with short tube
BND0490	Footrest lower part, depth and angle adjustable with short tube and heel strap, LH. Not compatible with Elevating legrest.
BND0520	Footrest lower part, depth and angle adjustable with long tube and heel strap, LH. All knee to heel ranges can be achieved when used with Elevating legrests.
BND0550	Footrest lower part, depth and angle adjustable with short tube and heel strap, RH. Not compatible with Elevating legrest.
BND0580	Footrest lower part, depth and angle adjustable with long tube and heel strap, RH. All knee to heel ranges can be achieved when used with Elevating legrests.
BND0680	Side Support Cover, PU fabric cover for the sideguards (Removable with tools)
BND9999	Side Support Cover with foam cushion. 2cm thick cushion supplied complete with PU fabric cover. Reduces seat width by 2cm for a comfortable fit. Consists of options 0680 and 0690 (Removable with tools)
BND0695	Pad for sideguard, removable with hook and loop fastening
BND0720	Pommel Abduction Cushion
BND0722	2 cm Increased seat height plate
BND0750	Headrest with cheek support, sideway adjustable Fixed hardware included
BND0742	Matrx Elan Hardware
BND0744	Matrx Elan Extended Hardware (Extra link)
BND0746	Matrx Elan LOXX Hardware
BND8888	Without Elan headrest hardware (Only select with option 0730 & 0750 Headrests)
BND0740	Cover for headrest
BND0770	Neckrest
BND0782	Alternative headrest/neckrest mounting bracket
BND0790	Handles Clinic Kit for Head-/Neckrest (Levers which replace screws for easy adjustment)
BND0800	Trunk Support Multi Adjustable with standard arm LH
BND0820	Trunk Support Multi Adjustable with extended height arm LH
BND0880	NEW!! Trunk Support Swingaway LH, multi adjustable
BND0810	Trunk Support Multi Adjustable with standard arm RH
BND0830	Trunk Support Multi Adjustable with extended height arm RH
BND0890	NEW!! Trunk Support Swingaway RH, multi adjustable
BND0900	Trunk Support Attachment to mount laterals later. If trunk supports are chosen, the attachment bracket is already included
BND0920	Wedges (One Pair)
BND0928	Calf pad cover, for LH elevating legrest
BND0929	Knee pad for LH elevating legrest
BND0931	Calf pad cover, for RH elevating legrest
BND0932	Knee pad for RH elevating legrest
BND6666	Without Elevating Legrest Options (Select with Fixed legrests, Amputee legrests and Without legrests only)

BND0930	Footplate Extender LH
BND0940	Footplate Extender RH
BND0941	Cover Footplate, Black Plush Fabric
BND0942	Calf Strap Padded
BND0945	Calf Strap (for One piece footrest)
BND0952	Bodypoint Ankle Huggers Large (Supplied with fixing kit, unmounted)
BND0954	Bodypoint Ankle Huggers Medium (Supplied with fixing kit, unmounted)
BND0956	Bodypoint Ankle Huggers Small (Supplied with fixing kit, unmounted)
BND1015	200x32 mm (8" x 1 1/4") Black Solid Tyre
BND1025	200x45 mm (8" X 1 3/4") Black Solid tyre
BND1055	22" / 560 mm Self Propel rear wheel with QR attachment
BND1062	24" / 600 mm Self Propel rear wheel with QR attachment
BND1074	Standard Self Propel with colourless hub, clear anodized rim and Stainless steel spokes
BND1076	Hub Brake 22" or 24" Self Propel with colourless hub, clear anodized rim and Stainless steel spokes
BND1206	Bodypoint Pivotfit Harness Medium (Supplied with fixing kit, unmounted)
BND1207	Bodypoint Pivotfit Harness Large (Supplied with fixing kit, unmounted)
BND1230	Fixed Table Tray without elbow pads includes mounting kit
BND1231	Fixed Table Tray with elbow Pads includes mounting kit
BND1232	Half table tray, clear plastic, flip up, depth adjustable, left side (armrest with straight narrow pad included)
BND1233	Half table tray, clear, plastic, flip up, depth adjustable, right side (armrest with straight narrow pad included)
BND1240	Swing away table without elbow pads
BND1245	Swing away table with elbow pads
BND1235	Table tray padded cushion (To be ordered in combination with the fixed & swingaway table tray)
BND1237	Half table tray padded cushion
BND1250	Locking device for table tray
BND1260	Mounting kit for table (if ordered separately)
BND1290	Handles Clinic Kit for backrest, armrest, legrest and seat frame (Levers which replace screws for easy adjustment)
BND1300	Harness attachment , must be selected in combination with Head/Neckrest



For a more complex customer requirement,
please email: orders@invacare.com.au

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):