



**Yes, you can:**

Tel: 01656 776222  
 Email: ordersuk@invacare.com  
 Online Spares available at: www.invacare.co.uk

# Power Retail Prescription Form

July 2026  
 Customer Ref:

<b>Customer Ref*</b>	<b>Account No*</b>
<b>Contact No*</b>	<b>Order No*</b>
<b>Delivery Address*</b>	<b>Order Qty*</b>
	<b>Order date*</b>

## Esprit Action

### Standard Features

- Max User Weight 125kg (19.7 Stone)
- Crash tested ISO7176 - 19\*
- Double crossbar
- Seat widths 380mm to 605mm (15" - 24")
- Seat depths 400mm, 450mm or 500mm (16", 18" or 20")
- Seat heights 460mm or 485mm (18" or 19")
- Nylon seat upholstery
- Fixed 7° 430mm backrest
- Lumbar Support (+1 STRAP) BI Material Padded backrest
- 80° swing in / swing out hangers
- Flip up footplates including heel strap
- Removable full length armrests
- Choice of 4 standard colours (with 9 special colour options)
- Anti tippers
- Buckle belt
- 
- \* Battery MUST be removed when used as a seat in a vehicle



**Standard Base price: £3,715**



Please make your selection with a cross in the box    ○ Free of charge    ◆ standard feature    POA special request

Customer weight =                      KG

Code	Description	
<b>BASE PRODUCT*</b>		
SPE0060	ESPRIT ACTION ADULT (DOUBLE CROSSBAR)	£3,715
<b>CHAIR TYPE*</b>		
SPE0100	SPECIAL	POA
SPE0100	BASE MODEL ESPRIT ACTION	◆
<b>CHASSIS*</b>		
SPE0210	DOUBLE CROSSBAR (SW 380 MM TO 605 MM / SD FROM 400 MM TO 500 MM) USER MAX WEIGHT = 125Kg	◆

BRACKET TYPE*			
SPE0300	<input type="checkbox"/>	STANDARD SWING AWAY BRACKET	◆
SPE0310	<input type="checkbox"/>	ALBER E35 SWIVEL ARM	£380
USER REMOTE*			
SPE0400	<input type="checkbox"/>	USER REMOTE CONTROL RH	○
SPE0410	<input type="checkbox"/>	USER REMOTE CONTROL LH	○
DUAL REMOTE			
SPE0440	<input type="checkbox"/>	INTUITIVE DUAL REMOTE CONTROL RH	£900
SPE0450	<input type="checkbox"/>	INTUITIVE DUAL REMOTE CONTROL LH	£900
SPE0460	<input type="checkbox"/>	INTUITIVE SINGLE DUAL REMOTE CONTROL (ONE HANDLE) (ADJUSTABLE RH TO LH SIDE)	£640
SPEED*			
SPE0500	<input type="checkbox"/>	6 KM/H	◆
BATTERY*			
SPE0600	<input type="checkbox"/>	BATTERY 7,65AH/36V (NON SPILLABLE LEAD ACID BATTERY)	◆
POWER OPTIONS			
SPE0610	<input type="checkbox"/>	ADDITIONAL BATTERY 7,65AH/36V	£285
SPE0630	<input type="checkbox"/>	REMOTE PROTECTOR	£70
SEAT WIDTH*			
SPE1040	<input type="checkbox"/>	380 MM SEAT WIDTH (NC 1250)	○
SPE1050	<input type="checkbox"/>	405 MM SEAT WIDTH (NC 1250)	○
SPE1060	<input type="checkbox"/>	430 MM SEAT WIDTH	○
SPE1070	<input type="checkbox"/>	455 MM SEAT WIDTH	○
SPE1080	<input type="checkbox"/>	480 MM SEAT WIDTH	○
SPE1090	<input type="checkbox"/>	505 MM SEAT WIDTH	○
SPE1100	<input type="checkbox"/>	555 MM SEAT WIDTH (REINFORCED SEAT UPHOLSTERY) (NC 2000, 2010, 2030, 2200, 2210, 2220, 2230, 5020, 5040, 5080, 8200)	£110
SPE1110	<input type="checkbox"/>	605 MM SEAT WIDTH (REINFORCED SEAT UPHOLSTERY) (NC 2000, 2010, 2030, 2200, 2210, 2220, 2230, 5020, 5040, 5070, 5080, 8200)	£110
SEAT DEPTH*			
SPE1210	<input type="checkbox"/>	400 MM SEAT DEPTH (NC 5090)	○
SPE1220	<input type="checkbox"/>	400 TO 450 MM ADJUSTABLE SEAT DEPTH SET AT 430 MM	○
SPE1230	<input type="checkbox"/>	450 MM SEAT DEPTH (NC 5080, 5090)	○
SPE1240	<input type="checkbox"/>	450 TO 500 MM ADJUSTABLE SEAT DEPTH (NC 5080)	£35
SPE1250	<input type="checkbox"/>	500 MM SEAT DEPTH (NC 5080)	○
SEAT HEIGHT*			
SPE1300	<input type="checkbox"/>	460 MM SEAT HEIGHT	○
SPE1310	<input type="checkbox"/>	485 MM SEAT HEIGHT	○
SEAT UPHOLSTERY TYPE*			
SPE1410	<input type="checkbox"/>	NYLON WITH LOOP BANDS STRIPS SEAT UPHOLSTERY (DOUBLE CROSSBAR ONLY)	◆

BACKREST TYPE*			
SPE2000		FIXED 7° BACKREST 2 HEIGHTS 400/430 MM (DELIVERED 430 MM)	◆
SPE2010		FIXED 7° BACKREST 510 MM HEIGHT	○
SPE2020		FOLDING BACKREST 430 MM HEIGHT	£35
SPE2030		FOLDING BACKREST 510 MM HEIGHT	£35
SPE2040		ANGLE ADJUSTABLE AND FOLD DOWN (-15° TO +15°) FIXED 7° HEIGHT 425-450 MM / 17"-18" SET AT 425 MM / 17" (2 POSITIONS; 25 MM / 1" STEPS THROUGH UPHOLSTERY)	£90
SPE2050		RECLINING/FOLDING (GAS SPRING) 0/30° BACKREST 3 HEIGHTS 390/410/430 MM (DELIVERED 430 MM)	(NC 2210, 2220, 2230) £155
SPE2060		RECLINING (RATCHET) 0/30° BACKREST 3 HEIGHTS 390/410/430 MM (DELIVERED 430 MM)	£135
SPE2070		RECLINING (RATCHET) 0/30° BACKREST 3 HEIGHTS 470/490/510 MM (DELIVERED 510 MM)	£135
SPE2080		RECLINING (GAS SPRING) 0/30° BACKREST 3 HEIGHTS 390/410/430 MM (DELIVERED 430 MM)	£55
SPE2090		RECLINING (GAS SPRING) 0/30° BACKREST 3 HEIGHTS 470/490/510 MM (DELIVERED 510 MM)	£70
BACKREST UPHOLSTERY TYPE*			
SPE2200		MATRIX ELITE E2 BACK WITH HARDWARE 405 MM BACK HEIGHT ADDITIONAL DELIVERY TIME + 5 DAYS	(NC 2310, 2320, 5140) £494
SPE2210		MATRIX ELITE E2 BACK WITH HARDWARE 455 MM BACK HEIGHT ADDITIONAL DELIVERY TIME + 5 DAYS	(NC 2310, 2320, 5140) £494
SPE2220		MATRIX ELITE DEEP E2 BACK WITH HARDWARE 405 MM BACK HEIGHT ADDITIONAL DELIVERY TIME + 5 DAYS	(NC 2310, 2320, 5140) £546
SPE2230		MATRIX ELITE DEEP E2 BACK WITH HARDWARE 455 MM BACK HEIGHT ADDITIONAL DELIVERY TIME + 5 DAYS	(NC 2310, 2320, 5140) £546
SPE2300		LUMBAR SUPPORT (+1 STRAP) BI MATERIAL PADDED COVER	◆
SPE2310		CONTOUR ADJUSTABLE (4+1 STRAPS) BI MATERIAL COVER	○
SPE2320		CONTOUR ADJUSTABLE (FULL STRAPS) BI MATERIAL COVER	○
HANGERS RH*			
SPE3010		80° LONG, UL FROM 340 TO 470MM	○
SPE3020		ANGLE ADJUSTABLE (0° TO 75°) RH SWING IN/SWING OUT HANGER, UL FROM 340 TO 500MM (INCLUDING CALFPAD)	(NC 5150) £95
SPE3030		AMPUTEE LEGREST RH	(NC 3210, 3220) £140
HANGERS LH*			
SPE3110		80° LONG, UL FROM 340 TO 470MM	○
SPE3120		ANGLE ADJUSTABLE (0° TO 75°) LH SWING IN/SWING OUT HANGER, UL FROM 340 TO 500MM (INCLUDING CALFPAD)	(NC 5150) £95
SPE3130		AMPUTEE LEGREST LH	(NC 3300, 3310) £140
FOOTPLATES RH*			
SPE3200		FLIP UP RH FOOTPLATE INCLUDING HEEL STRAP	○
SPE3210		FLIP UP ANGLE AND DEPTH ADJUSTABLE RH FOOTPLATE, INCLUDING HEEL STRAP	£45
SPE3220		FLIP UP ONE PIECE FOOTPLATE, DEPTH & ANGLE ADJUSTABLE	(NC 3310) £80
SPZ7000		AMPUTEE LEGREST OR FOOTBOARD SELECTED	(NC 3010, 3020) ○

FOOTPLATES LH*			
SPE3300		FLIP UP LH FOOTPLATE, INCLUDING HEEL STRAP	○
SPE3310		FLIP UP ANGLE AND DEPTH ADJUSTABLE LH FOOTPLATE, INCLUDING HEEL STRAP	£45
SPZ7001		AMPUTEE LEGREST OR FOOTBOARD SELECTED	○
ARMRESTS RH*			
SPE4000		REMOVABLE RH DESK, ADJ IN HEIGHT (235-280 for FSTF 460 & 210-280 for FSTF 485MM) DEPTH (+35 MM FRONTWARD) & WIDTH (+25MM OUTWARDS)	○
SPE4010		REMOVABLE RH FULL LENGTH, ADJ IN HEIGHT (210 TO 280 MM) DEPTH (+35 MM BACKWARD) & WIDTH (+25MM OUTWARDS)	◆
SPE4020		REMOVABLE RH ARMPAD HEMI, BLACK PU TR26, MULTIPOSITION AUTOLOCK SAFETY DEVICE INCLUDED, ADJ IN HEIGHT (210 TO 345 MM) DEPTH (35 MM FRONTWARD) & WIDTH (+25MM OUTWARDS)	£35
ARMRESTS LH*			
SPE4100		REMOVABLE LH DESK, ADJ IN HEIGHT (235-280 for FSTF 460 & 210-280 for FSTF 485MM) DEPTH (+35 MM FRONTWARD) & WIDTH (+25MM OUTWARDS)	○
SPE4110		REMOVABLE LH FULL LENGTH, ADJ IN HEIGHT (210 TO 280 MM) DEPTH (35 MM BACKWARD) & WIDTH (+25MM OUTWARDS)	◆
SPE4120		REMOVABLE LH ARMPAD HEMI, BLACK PU TR26, MULTIPOSITION AUTOLOCK SAFETY DEVICE INCLUDED, ADJ IN HEIGHT (210 TO 345 MM) DEPTH (35 MM FRONTWARD) & WIDTH (+25MM OUTWARDS)	£35
SEAT CUSHION*			
SPE5010		BASIC 50 MM FOAM CUSHION - 50 MM HEIGHT, NYLON & ANTI-SLIP COVER	(NC 5092, 5093, 5091, 5092, 5093) ○
SPE5020		MATRX FLO-TECH LITE CUSHION - 78 MM HEIGHT	£120
SPE5040		MATRX CONTOUR VISCO NG CUSHION - 100 MM HEIGHT MAXIMUM USER WEIGHT 100 KG	(NC 5091, 5092, 5093, 5091, 5092, 5093) £169
SPE5070		MATRX LIBRA CUSHION - 70 MM HEIGHT, ADDITIONAL DELIVERY TIME + 5 DAYS	£410
SEAT CUSHION DEPTH*			
SPE5080		FOR 400 MM SEAT DEPTH	○
SPE5085		FOR 450 MM SEAT DEPTH	○
SPE5090		FOR 500 MM SEAT DEPTH	○
1st POSITIONING ACCESORIES FOR MATRX LIBRA			
SPE5091		Tapered Wedges (Pair)	£30
SPE5092		Lateral Wedge	£20
SPE5093		Anterior Seat Wedges	£27
2nd POSITIONING ACCESORIES FOR MATRX LIBRA			
SPE5091		Tapered Wedges (Pair)	£30
SPE5092		Lateral Wedge	£20
SPE5093		Anterior Seat Wedges	£27
SEATING OPTIONS			
SPE5100		NECKREST BLACK PU TR26	£135
SPE5110		HEADREST BLACK PU TR26	£165
SPE5120		HEADREST WITH CHEEK SUPPORT BLACK PU TR26 (INCLUDING BACKBRACE EXCEPT FOR LAGUNA)	£190
SPE5140		LATERAL POSITIONING PADS (PAIR)	£90
SPE5150		ADJUSTABLE CALF STRAP	○

<b>BACK BRACE</b>		
SPE5170	<input type="checkbox"/>	BACK BRACE FIXED ○
<b>COLOURS*</b>		
SPE6000	<input type="checkbox"/>	ANTHRACITE BLACK ○
SPE6015	<input type="checkbox"/>	DARK PURPLE ○
SPE6035	<input type="checkbox"/>	NAVY BLUE ○
SPE6040	<input type="checkbox"/>	SLATE GREY ○
SWO6925	<input type="checkbox"/>	WATER BLUE - 3 week lead time ○
SWO0806	<input type="checkbox"/>	Electric Green - 3 week lead time £175
SWO6033	<input type="checkbox"/>	Red Matt - 3 week lead time £175
SWO0607	<input type="checkbox"/>	Dark Red - 3 week lead time £175
SWO0515	<input type="checkbox"/>	Pearl Grey - 3 week lead time £175
SWO0647	<input type="checkbox"/>	Intense Yellow - 3 week lead time £175
SWO0648	<input type="checkbox"/>	Vibrant Orange - 3 week lead time £175
SWO0651	<input type="checkbox"/>	Polar Blue - 3 week lead time £175
SWO0652	<input type="checkbox"/>	Candy Pink - 3 week lead time £175
SPE6050	<input type="checkbox"/>	SPECIAL COLOUR (ON DEMAND ONLY) POA
<b>ACCENT COLOUR*</b>		
SPE6100	<input type="checkbox"/>	BLACK ○
SPE6120	<input type="checkbox"/>	LIME ○
<b>FRONT WHEELS*</b>		
SPE7000	<input type="checkbox"/>	200 X 50 MM (8" X 2") BLACK SOLID TYRE £50
SPE7010	<input type="checkbox"/>	200 X 50 MM (8" X 2") BLACK PNEUMATIC TYRE ◆
<b>REAR WHEELS*</b>		
SPE7100	<input type="checkbox"/>	355 X 70 MM (14" X 27") BLACK SOLID WITH PUNCTURE PROOF TYRE £90
SPE7110	<input type="checkbox"/>	355 X 70 MM (14" X 27") BLACK PNEUMATIC TYRE ◆
<b>ANTI TIPPER*</b>		
SPE8000	<input type="checkbox"/>	REAR ANTI-TIPPERS, FIXED, PAIR ◆
SPE8010	<input type="checkbox"/>	REAR ANTI-TIPPERS REMOVABLE WITH JACK UP, PAIR £290
<b>LAP BELT*</b>		
SPE8020	<input type="checkbox"/>	POSTURE BELT ◆
<b>TRANSPORTATION BRACKET*</b>		
SPE8030	<input type="checkbox"/>	TRANSPORTATION BRACKET ◆
<b>OTHER OPTIONS</b>		
SPE8100	<input type="checkbox"/>	CANE HOLDER £25
SPE8200	<input type="checkbox"/>	ACRYLIC TRAY FLIP UP £325

<b>TOTAL PRICE</b>	
--------------------	--

**Non Conforming Options:**

SPE1250	500 MM SEAT DEPTH
SPE2000	FIXED 7° BACKREST 2 HEIGHTS 400/430 MM (DELIVERED 430 MM)
SPE2010	FIXED 7° BACKREST 510 MM HEIGHT
SPE2030	FOLDING BACKREST 510 MM HEIGHT
SPE2200	MATRX ELITE E2 BACK WITH HARDWARE 405 MM BACK HEIGHT ADDITIONAL DELIVERY TIME + 5 DAYS
SPE2210	MATRX ELITE E2 BACK WITH HARDWARE 455 MM BACK HEIGHT ADDITIONAL DELIVERY TIME + 5 DAYS
SPE2220	MATRX ELITE DEEP E2 BACK WITH HARDWARE 405 MM BACK HEIGHT ADDITIONAL DELIVERY TIME + 5 DAYS
SPE2230	MATRX ELITE DEEP E2 BACK WITH HARDWARE 455 MM BACK HEIGHT ADDITIONAL DELIVERY TIME + 5 DAYS
SPE2310	CONTOUR ADJUSTABLE (4+1 STRAPS) BI MATERIAL COVER
SPE2320	CONTOUR ADJUSTABLE (FULL STRAPS) BI MATERIAL COVER
SPE3010	80° LONG, UL FROM 340 TO 470MM
SPE3020	ANGLE ADJUSTABLE (0° TO 75°) RH SWING IN/SWING OUT HANGER, UL FROM 340 TO 500MM (INCLUDING CALFPAD)
SPE3210	FLIP UP ANGLE AND DEPTH ADJUSTABLE RH FOOTPLATE, INCLUDING HEEL STRAP
SPE3220	FLIP UP ONE PIECE FOOTPLATE, DEPTH & ANGLE ADJUSTABLE
SPE3300	FLIP UP LH FOOTPLATE, INCLUDING HEEL STRAP
SPE3310	FLIP UP ANGLE AND DEPTH ADJUSTABLE LH FOOTPLATE, INCLUDING HEEL STRAP
SPE5020	MATRX FLO-TECH LITE CUSHION - 78 MM HEIGHT
SPE5040	MATRX CONTOUR VISCO NG CUSHION - 100 MM HEIGHT MAXIMUM USER WEIGHT 100 KG
SPE5070	MATRX LIBRA CUSHION - 70 MM HEIGHT, ADDITIONAL DELIVERY TIME + 5 DAYS
SPE5080	FOR 400 MM SEAT DEPTH
SPE5090	FOR 500 MM SEAT DEPTH
SPE5091	Tapered Wedges (Pair)
SPE5092	Lateral Wedge
SPE5093	Anterior Seat Wedges
SPE5091	Tapered Wedges (Pair)
SPE5092	Lateral Wedge
SPE5093	Anterior Seat Wedges
SPE5140	LATERAL POSITIONING PADS (PAIR)
SPE5150	ADJUSTABLE CALF STRAP
SPE8200	ACRYLIC TRAY FLIP UP



For a more complex customer requirement,  
please email: [orders@invacare.com.au](mailto:orders@invacare.com.au)

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):