



Yes, you can:

Tel: 01656 776222
 Email: ordersuk@invacare.com
 Online Spares available at: www.invacare.co.uk

Power Retail Prescription Form

July 2026
 Customer Ref:

| | |
|--------------------------|--------------------|
| Customer Ref* | Account No* |
| Contact No* | Order No* |
| Delivery Address* | Order Qty* |
| | Order date* |

Esprit Action Junior

Standard Features

- Max User Weight 80kg (12.6 Stone)
- Crash tested ISO7176 - 19*
- Single crossbar (Width Adjustable)
- Seat widths 305mm to 380mm 12" - 15"
- Seat depth 325mm (13") to 375mm (15") adjustable (set at 350mm)
- Seat height 460mm (18")
- Nylon seat upholstery
- FIXED 7° BACKREST 430mm backrest
- Standard nylon sling type backrest upholstery
- 80° swing in / swing out legrest hangers
- Flip up angle and depth adjustable footplates
- Height adjustable, removable arm armrests
- Choice of 3 standard colours (with 10 Special colour options)
-
- * Battery MUST be removed when used as a seat in a vehicle

Standard Base price: £3,715



Please make your selection with a cross in the box ○ Free of charge ◆ standard feature POA special request

Customer weight = KG

| Code | Description | |
|----------------------|--|--------|
| BASE PRODUCT* | | |
| SPD0050 | ESPRIT ACTION JUNIOR (SINGLE CROSSBAR) (NC 2300) | £3,715 |
| CHAIR TYPE* | | |
| SPDSPEC | SPECIAL | POA |
| SPD0100 | BASE MODEL ESPRIT ACTION | ◆ |
| CHASSIS* | | |
| SPD0200 | SINGLE CROSSBAR (SW 305 MM TO 380 MM / SD 325 MM TO 375 MM) USER MAX WEIGHT = 80Kg | ◆ |

| | | | |
|----------------------------------|--|--|------|
| BRACKET TYPE* | | | |
| SPD0300 | | STANDARD SWING AWAY BRACKET | ○ |
| SPD0310 | | ALBER E35 SWIVEL ARM | £380 |
| USER REMOTE* | | | |
| SPD0400 | | USER REMOTE CONTROL RH | ○ |
| SPD0410 | | USER REMOTE CONTROL LH | ○ |
| DUAL REMOTE | | | |
| SPD0440 | | INTUITIVE DUAL REMOTE CONTROL RH | £900 |
| SPD0450 | | INTUITIVE DUAL REMOTE CONTROL LH | £900 |
| SPD0460 | | INTUITIVE SINGLE DUAL REMOTE CONTROL (ONE HANDLE) (ADJUSTABLE RH TO LH SIDE) | £640 |
| SPEED* | | | |
| SPD0500 | | 6 KM/H | ◆ |
| BATTERY* | | | |
| SPD0600 | | BATTERY 7,65AH/36V (NON SPILLABLE LEAD ACID BATTERY) | ◆ |
| POWER OPTIONS | | | |
| SPD0610 | | ADDITIONAL BATTERY 7,65AH/36V | £285 |
| SPD0630 | | REMOTE PROTECTOR | £70 |
| SEAT WIDTH* | | | |
| SPD1000 | | 305 MM SEAT WIDTH | ○ |
| SPD1010 | | 330 MM SEAT WIDTH | ○ |
| SPD1020 | | 355 MM SEAT WIDTH | ○ |
| SPD1030 | | 380 MM SEAT WIDTH | ○ |
| SEAT DEPTH* | | | |
| SPD1200 | | 325 TO 375 MM ADJUSTABLE SEAT DEPTH SET AT 350 MM | ○ |
| SEAT HEIGHT* | | | |
| SPD1300 | | 460 MM SEAT HEIGHT | ○ |
| SPD1310 | | 485 MM SEAT HEIGHT | ○ |
| SEAT UPHOLSTERY TYPE* | | | |
| SPD1400 | | NYLON WITH LOOP BANDS STRIPS SEAT UPHOLSTERY, WIDTH & DEPTH ADJUSTABLE (SINGLE CROSSBAR ONLY) | ◆ |
| BACKREST TYPE* | | | |
| SPD2000 | | FIXED 7° BACKREST 2 HEIGHTS 400/430 MM (DELIVERED 430 MM) | ◆ |
| SPD2020 | | FOLDING BACKREST 430 MM HEIGHT | £35 |
| SPD2040 | | ANGLE ADJUSTABLE AND FOLD DOWN (-15° TO + 15°) FIXED 7° HEIGHT 425-450 MM / 17"-18" SET AT 425 MM / 17" (2 POSITIONS; 25 MM / 1" STEPS THROUGH UPHOLSTERY) | £90 |
| SPD2050 | | RECLINING/FOLDING (GAS SPRING) 0/30° BACKREST 3 HEIGHTS 390/410/430 MM (DELIVERED 430 MM) | £155 |
| SPD2060 | | RECLINING (RATCHET) 0/30° BACKREST 3 HEIGHTS 390/410/430 MM (DELIVERED 430 MM) | £135 |
| SPD2080 | | RECLINING (GAS SPRING) 0/30° BACKREST 3 HEIGHTS 390/410/430 MM (DELIVERED 430 MM) | £135 |
| BACKREST UPHOLSTERY TYPE* | | | |
| SPD2300 | | LUMBAR SUPPORT (+1 STRAP) BI MATERIAL PADDED COVER | ◆ |
| SPD2320 | | CONTOUR ADJUSTABLE (FULL STRAPS) BI MATERIAL COVER | ○ |

| HANGERS RH* | | |
|------------------------|--------------------------|---|
| SPD3000 | <input type="checkbox"/> | 80° SHORT, UL FROM 200 TO 330MM ◆ |
| HANGERS LH* | | |
| SPD3100 | <input type="checkbox"/> | 80° SHORT, UL FROM 200 TO 330MM ◆ |
| FOOTPLATES RH* | | |
| SPD3200 | <input type="checkbox"/> | FLIP UP RH FOOTPLATE INCLUDING HEEL STRAP ○ |
| SPD3210 | <input type="checkbox"/> | FLIP UP ANGLE AND DEPTH ADJUSTABLE RH FOOTPLATE, INCLUDING HEEL STRAP £45 |
| SPD3220 | <input type="checkbox"/> | FLIP UP ONE PIECE FOOTPLATE, DEPTH & ANGLE ADJUSTABLE (NC 3310) £80 |
| FOOTPLATES LH* | | |
| SPD3300 | <input type="checkbox"/> | FLIP UP LH FOOTPLATE, INCLUDING HEEL STRAP ○ |
| SPD3310 | <input type="checkbox"/> | FLIP UP ANGLE AND DEPTH ADJUSTABLE LH FOOTPLATE, INCLUDING HEEL STRAP £45 |
| SPD9925 | <input type="checkbox"/> | No LH footplate required due to FLIP UP ONE PIECE BEING SELECTED (NC 3210) ○ |
| ARMRESTS RH* | | |
| SPD4000 | <input type="checkbox"/> | REMOVABLE RH DESK, ADJ IN HEIGHT (235-280 for FSTF 460 & 210-280 for FSTF 485MM) DEPTH (+35 MM FRONTWARD) & WIDTH (+25MM OUTWARDS) ○ |
| SPD4010 | <input type="checkbox"/> | REMOVABLE RH FULL LENGTH, ADJ IN HEIGHT (210 TO 280 MM) DEPTH (+35 MM BACKWARD) & WIDTH (+25MM OUTWARDS) ◆ |
| ARMRESTS LH* | | |
| SPD4100 | <input type="checkbox"/> | REMOVABLE LH DESK, ADJ IN HEIGHT (235-280 for FSTF 460 & 210-280 for FSTF 485MM) DEPTH (+35 MM FRONTWARD) & WIDTH (+25MM OUTWARDS) ○ |
| SPD4110 | <input type="checkbox"/> | REMOVABLE LH FULL LENGTH, ADJ IN HEIGHT (210 TO 280 MM) DEPTH (35 MM BACKWARD) & WIDTH (+25MM OUTWARDS) ◆ |
| SEAT CUSHION* | | |
| SPD5000 | <input type="checkbox"/> | ESSENTIAL CUSHION 375 MM DEPTH ○ |
| SEATING OPTIONS | | |
| SPD5110 | <input type="checkbox"/> | HEADREST BLACK PU TR26 £165 |
| SPD5120 | <input type="checkbox"/> | HEADREST WITH CHEEK SUPPORT BLACK PU TR26 (INCLUDING BACKBRACE EXCEPT FOR LAGUNA) £190 |
| SPD5140 | <input type="checkbox"/> | LATERAL POSITIONING PADS (PAIR) £90 |
| SPD5150 | <input type="checkbox"/> | ADJUSTABLE CALF STRAP ○ |
| BACK BRACE | | |
| SPD5160 | <input type="checkbox"/> | BACK BRACE WIDTH ADJUSTABLE ○ |

| COLOURS* | | | |
|--------------------------------|--------------------------|--|------|
| SPD6000 | <input type="checkbox"/> | ANTHRACITE BLACK | ○ |
| SPD6015 | <input type="checkbox"/> | DARK PURPLE | ○ |
| SPD6020 | <input type="checkbox"/> | WATER BLUE | ○ |
| SWO6529 | <input type="checkbox"/> | POP BLUE - 3 week lead time | £175 |
| SWO6259 | <input type="checkbox"/> | SLATE GREY - 3 week lead time | £175 |
| SWO0806 | <input type="checkbox"/> | Electric Green - 3 week lead time | £175 |
| SWO6033 | <input type="checkbox"/> | Red Matt - 3 week lead time | £175 |
| SWO0607 | <input type="checkbox"/> | Dark Red - 3 week lead time | £175 |
| SWO0515 | <input type="checkbox"/> | Pearl Grey - 3 week lead time | £175 |
| SWO0647 | <input type="checkbox"/> | Intense Yellow - 3 week lead time | £175 |
| SWO0648 | <input type="checkbox"/> | Vibrant Orange - 3 week lead time | £175 |
| SWO0651 | <input type="checkbox"/> | Polar Blue - 3 week lead time | £175 |
| SWO0652 | <input type="checkbox"/> | Candy Pink - 3 week lead time | £175 |
| SPD6050 | <input type="checkbox"/> | SPECIAL COLOUR (ON DEMAND ONLY) | POA |
| ACCENT COLOUR* | | | |
| SPD6100 | <input type="checkbox"/> | BLACK | ○ |
| SPD6120 | <input type="checkbox"/> | LIME | ○ |
| FRONT WHEELS* | | | |
| SPD7000 | <input type="checkbox"/> | 200 X 50 MM (8" X 2") BLACK SOLID TYRE | £50 |
| SPD7010 | <input type="checkbox"/> | 200 X 50 MM (8" X 2") BLACK PNEUMATIC TYRE | ◆ |
| REAR WHEELS* | | | |
| SPD7100 | <input type="checkbox"/> | 355 X 70 MM (14" X 27") BLACK SOLID WITH PUNCTURE PROOF TYRE | £90 |
| SPD7110 | <input type="checkbox"/> | 355 X 70 MM (14" X 27") BLACK PNEUMATIC TYRE | ◆ |
| ANTI TIPPER* | | | |
| SPD8000 | <input type="checkbox"/> | REAR ANTI-TIPPERS, FIXED, PAIR | ◆ |
| SPD8010 | <input type="checkbox"/> | REAR ANTI-TIPPERS REMOVABLE WITH CRUTCHES (E25), PAIR | £290 |
| LAP BELT* | | | |
| SPD8020 | <input type="checkbox"/> | POSTURE BELT | ◆ |
| TRANSPORTATION BRACKET* | | | |
| SPD8030 | <input type="checkbox"/> | TRANSPORTATION BRACKET | ◆ |
| OTHER OPTIONS | | | |
| SPD8100 | <input type="checkbox"/> | CANE HOLDER | £25 |
| SPD8200 | <input type="checkbox"/> | ACRYLIC TRAY FLIP UP | £325 |
| TOTAL PRICE | | | |

Non Conforming Options:

| | |
|---------|---|
| SPD2300 | LUMBAR SUPPORT (+1 STRAP) BI MATERIAL PADDED COVER |
| SPD3210 | FLIP UP ANGLE AND DEPTH ADJUSTABLE RH FOOTPLATE, INCLUDING HEEL STRAP |
| SPD3310 | FLIP UP ANGLE AND DEPTH ADJUSTABLE LH FOOTPLATE, INCLUDING HEEL STRAP |



For a more complex customer requirement,
please email: orders@invacare.com.au

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):