

# Invacare Aquatec Ocean Shower Chair Commode

## Product Data Sheet

---



### Product Range Available

- Ocean Transit  
Product Code: 1517044
- Ocean Self Propel  
Product Code: 1517045
- Ocean XL Transit  
Product code: 1517058
- Ocean XL Self Propel  
Product Code: 1517059

All products are delivered with soft seat, insert for soft seat, pan and lid as standard

Cross compatibility on a wide range of accessories available on all products

### Standard Features

- Modular shower chair commode
- Tension adjustable back
- Stable, powder coated stainless steel frame – 100% corrosion free
- 4 Locking castors
- Angle adjustable, flip up and textured footplates with heel loops
- Height adjustable, swing away and detachable footrests
- Ergonomic swing away armrests, can be adjusted to two different height

### Technical Specification

#### Ocean Transit & Self Propel and Ocean XL

- Armrests height 210 or 250mm
- Height, footrest to seat 255 – 495mm (no tools required)
- Adjustable seat to ground height 475 – 600mm (no tools required)
- **Ocean XL Only:** Increased width between armrests 540mm
- 3 Year Warranty

- 100% corrosion free frame
- Maximum user weight 150kg (23 stone) (Ocean, Ocean VIP, Ocean VIP Dual, and Ocean VIP XL)
- **Ocean XL and Ocean XL Self Propel** Maximum user weight 180kg (28 stone)

#### **Optional accessories available including:**

- Sot seat and Visco soft seat
- Variable multi adjustable soft seat (can be retrofitted – no tools required)
- Solid plastic backrest – clips onto frame
- Hip and chest belts
- Armrest lock – lock in upright and down position
- Lateral supports
- Amputee and calf supports – no tools required
- 24” self-propel wheels

All accessories can be found on the product brochure

#### **Supporting Documentation**

- Aquatec Ocean range brochure PN: 1504395
- Ocean A5 quick guide PN: 1575042
- Ocean Poster PN: 1575974
- Ocean Order guide\*
- Ocean packs are also available email [UK@invacare.com](mailto:UK@invacare.com) to order your pack/packs
- Invacare’s Business Development Managers can provide support with client assessments and provide training for OT’s and carers.

*\*Contact marketing department*

---

#### **Manufacturers Contact Details**

Invacare Ltd, Unit 4 Pencoed Technology Park, Bridgend, CF35 5AQ

**Tel:** 01656 776222 **E-mail:** [UK@invacare.com](mailto:UK@invacare.com)

**Web:** [www.invacare.co.uk](http://www.invacare.co.uk)