

Invacare® | **Essential Comfort**



Cost effective comfort mattress

The *Invacare Essential Comfort* is a basic entry level pressure redistributing mattress designed for the community healthcare environment. Manufactured from a highly resilient foam, this mattress provides a good level of comfort at a cost effective price.





Features and Options

Comfort and stability

The *Invacare Essential Comfort* mattress is constructed from a highly resilient foam, providing a good level of comfort and a stable support surface for the patient.

Cover is resistant to fluid ingress

The outer cover is a multi stretch, water resistant, polyurethane fabric providing excellent recovery and durability which help to reduce shear and friction forces. The cover has a white inner liner allowing easy inspection for strikethrough or other cover damage.



Turn and rotate

Mattress can be turned or rotated for increased product longevity.



Available in the following sizes:

L2000 x W880 x D120mm
Part no: 1546004

L2000 x W855 x D120mm
Part no: 1538897

L1970 x W880 x D120mm
Part no: 1538896




Zip around three sides


Allows easy inspection of the foam core and easy removal for laundering.

Technical data

Invacare Essential Comfort



Mattress sizes :
2000 x 855 mm
2000 x 880mm
and 1970 x 880 mm



Height:
120 mm



Weight:
9 kg



Max. User Weight:
110 kg / 17 stone



Cleaning:
80°

Mattress is manufactured to:

- BS EN 597-1 Cig, BS EN 597-2 Match and BS 7177: Crib Source 5

Colour



Baltic Blue



1 year warranty

Invacare Limited
Pencoed Technology Park - Pencoed -
Bridgend CF35 5AQ - United Kingdom
Tel: +44 1 656 776 222
Fax: +44 1656 776 220
E-mail: uk@invacare.com
Sales Order E-mail:
ordersuk@invacare.com
www.invacare.co.uk

Invacare Ireland Ltd
Unit 5 - Seatown Business Campus
Seatown Road - SWORDS - County
Dublin - Ireland
Tel. +353 1 8107084
Fax +353 1 8107085
Email: Ireland@invacare.com
www.invacare.ie

© 2014 Invacare International Sàrl
All rights reserved.

All information quoted is believed to be correct at time of print.
Invacare reserves the right to alter product specifications without prior consultation.

Essential Comfort - UK - 07/2014



Yes, you can.®