

# *ProFlex Mattress Overlay*



Good value pressure  
redistributing overlay



The **ProFlex Overlay** has been designed to be used in conjunction with an existing mattress to provide a cost effective Pressure Ulcer Management solution. Designed for those at 'At Risk/Increased Vulnerability'\* of developing pressure ulcers, the **Proflex overlay** is constructed from a contoured foam to maximise pressure reduction. The cover is a 2 way stretch waterproof and vapour permeable fabric.

### Effective pressure re-distribution

The contoured foam affords pressure re-distribution and cradles bony prominences. Made from high quality foam, the **ProFlex Overlay** offers a good life span and the possibility of the overlay "bottoming out" is reduced. **ProFlex Overlay** sits on top of the existing mattress to provide a comfortable surface for the patient.

### Even weight distribution

The **ProFlex Overlay** has cut foam castellations to help to reduce shear and friction forces.

## Features and Options



### Cover helps to prevent a build up of shear and friction forces

Cover is multi-stretch, water resistant and vapour-permeable.



### Covers have an anti-slip base

Ensures safe use of the **Pro Flex Overlay** and helps to prevent slippage.

## Technical data

For more information about this product, including the product's user manual, please visit [www.invacare.co.uk](http://www.invacare.co.uk).

Size range	Height	Total product weight	Max. user weight	Max. cleaning temperature
				
Size: 1875 x 876 mm	85 mm	5.5 kg	108 kg 17 Stone	80°

Manufactured to European Standards: BS 7177 Crib 5

### Colour

Baltic Blue



### Warranty

1 Year

\* Essential nursing care is pivotal in pressure ulcer prevention. This product will positively contribute to the outcome of a pressure ulcer prevention care plan. Education, clinical judgement and action based planning based on vulnerability are fundamental factors in the prevention of pressure ulcers. A range of assessment scales can be used as a formal method of assessing risk from pressure ulcer development, and should be used in conjunction with an informal assessment (informed nursing judgement). Informal assessment is considered to be of greater importance and clinical value.

Invacare Limited  
Pencoed Technology Park - Pencoed  
Bridgend CF35 5AQ - United Kingdom  
Tel: +44 1656 776 222  
Fax: +44 1656 776 220  
E-mail: [uk@invacare.com](mailto:uk@invacare.com)  
Sales Order E-mail:  
[ordersuk@invacare.com](mailto:ordersuk@invacare.com)  
[www.invacare.co.uk](http://www.invacare.co.uk)

Invacare Ireland Ltd  
Unit 5 - Seatown Business Campus  
Seatown Road - SWORDS -  
County Dublin - Ireland  
Tel. +353 1 8107084  
Fax +353 1 8107085  
Email: [Ireland@invacare.com](mailto:Ireland@invacare.com)  
[www.invacare.ie](http://www.invacare.ie)

